

Distel DS laboratory disinfectant.

PRODUCT APPLICATION FOR THE DECONTAMINATION OF PIPETTES

The following are recommended decontamination procedures for all models. Please note that these instructions are intended as a **GUIDE ONLY** and all pipettes requiring decontamination must be decontaminated according to a) your institution's own guidelines and b) the type of contaminant present.

GENERAL DAILY CLEANING/DECONTAMINATION

- Use Distel DS at the recommended working solution
- Spray or wipe handle and shaft
(insuring the tip ejector button and volumeter are also cleaned)
- Leave for a minimum contact time of 10 minutes



1 L DISTEL DS STARTER PACK	1 L DISTEL DS CONCENTRATE
<ul style="list-style-type: none">● 1 x 1 L Bottle Distel DS Concentrate● 5 x 1 L Empty Reusable Dosing Bottles● 6 x 5.0 ml Dosing Heads● 1 x Working Solution Bottle Label	<ul style="list-style-type: none">● 1 x 1 L Bottle Distel DS Concentrate● 1 x 5.0 ml Dosing Head● 1 x Working Solution Bottle Label
ARTICLE NO.	ARTICLE NO.
CCH050101	CCH050201

COMPLETE DECONTAMINATION

Pipettes requiring full decontamination can be dismantled and cleaned:

- Use Distel DS at recommended working solution
- Use a container that is deep enough to immerse the parts completely.
(NOTE: Electronic pipettes must not be submerged in any liquids)
- Dismantle the pipette (as specified by the manufacturer) and **place only parts that can be immersed**, eg tip holder, piston, tip ejector etc completely into the solution for 30 minutes
- Remove parts and rinse with de-mineralised water
- Lie on a lint-free cloth and place in a dry oven at 60°C for at least one hour or until dry. Alternatively, place on the bench on a lint-free cloth and air dry for approximately 4 hours
- Reassemble the pipette following the manufacturer's instructions
- If complete decontamination is required prior to servicing, please ensure all parts are completely dry before reassembly and presenting for servicing

NOTE: All pipettes for servicing must be accompanied by a completed Decontamination Declaration stating that the pipette is safe for us to work on, irrespective of whether the pipette required decontamination or not.

DECONTAMINATION DECLARATIONS

starlab.click/decontamination