Passionate for pipetting.

valid to 26 April 2024

SPECIAL OFFERS

ErgoOne® Pipettes

Pipette Starter Packs

Microdilution Tubes

Reagent Reservoirs

PCR Plates & Seals

Plate Seals

Plate Centrifuge

Minicentrifuge

Improve your pipetting experience with Starlab ErgoOne® pipettes and complement with a selection of high quality products for your daily pipetting tasks.





Starlab ErgoOne® manual pipettes.

STYLE AND COMFORT: YOU CAN HAVE BOTH!

Reliability is key in the laboratory and your pipette must be trusted to perform well, day in, day out. ErgoOne® manual pipettes provide accuracy and precision you can trust, alongside ergonomic comfort and a stylish look you will love.

- Ergonomic and lightweight
- Volume lock
- Universal tip fit



starlab.click/tipone



Single-Channel Pipettes

S7100-0125	S7100-0510	S7100-0221	\$7100-0220
0,1-2,5 µl	0,5 – 10 µl	2–20 µl Micro	2 – 20 µl Standard
130 €*	130 €*	130 €*	130 €*
S7100-1100	S7100-2200	S7110-1000	S7150-5000
10 – 100 µl	20 – 200 µl	$100 - 1.000 \mu l$	500 – 5.000 µl
130 €*	130 €*	130 €*	130 €*

8-Channel Pipettes

S7108-0510	S7108-1100	S7108-3300
0,5 – 10 µl	10 – 100 µl	30 – 300 µl
380 €*	380 €*	380 €*

12-Channel Pipettes

S7112-0510	S7112-1100	S7112-3300
0,5 – 10 µl	10 – 100 µl	$30 - 300 \mu l$
415 €*	465 €*	465 €*



^{*}Offer prices, additional counts not applicable. Valid until 26.04.2024.

ErgoOne® starter packs.

CREATE THE PERFECT PIPETTING PARTNERSHIP

Combine the reliability of Starlab's ErgoOne® pipettes with the trusted quality and performance of Starlab's TipOne® pipette tips for the perfect pipetting partnership.



CONTENTS		STARTERPACK	STARTERPACK42.5	STARTERPACK410
0,1 – 2,5 μl ErgoOne® Single-Ch Pipette S7100-0125			•	
0,5 – 10 μl ErgoOne® Single-Ch Pipette S7100-0510				•
2,0 – 20 μl ErgoOne® Single-Ch Pipette, Stan S7100-0220	dard Cone	•	•	•
20 – 200 μl ErgoOne® Single-Ch Pipette S7100-2200		•	•	•
100 – 1000 μl ErgoOne® Single-Ch Pipette S7110-1000		•	•	•
TipOne® Pipette Tips (three racks of sterile tips per pack)		•	•	•
Self-Adhesive Pipette Holder (one holder per pipette supplied)		•	•	•
Laboratory Timer (one timer per pack)		•	•	•
	OFFER*	445 €	595 €	595 €

*Offer prices, additional discounts not applicable. Valid until 26.04.2024.





Microdilution tubes.

1.2 ml TUBES IN REUSABLE RACKS

Choose between individual tubes and 8- or 12-strips. These autoclavable, SBS footprint racks are ideal for automation and long term storage. You can reuse the racks too! The tubes and caps are available separately for refilling the racks



MICRODILU [.]	TION TUBE	S	
Description	Pack Size	Article No.	OFFER*
1.2 ml Individua	al Tubes, Racke	ed	
Individual Tubes in Rack	10 racks	E1750-0400	42,75 €*
Individual Tubes in Rack, Sterile	10 racks	E1750-0410	54,45 €*
Individual Reaction Tubes only	1000 tubes	E1710-0000	21,60 €*
1.2 ml 8-Strip	Tubes, Racked		
8-Strip Tubes in Rack	10 racks	E1760-8400	42,75 €*
8-Strip Tubes in Rack, Sterile	10 racks	E1760-8410	49,05 €
8-Strip Reaction Tubes only	125 strips	E1720-8000	24,08 €
1.2 ml 12-Strip	Tubes, Racked	i	
12-Strip Tubes in Rack	10 racks	E1770-9400	49,05 €*
12-Strip Tubes in Rack, Sterile	10 racks	E1770-9410	54,45 €
12-Strip Reaction Tubes only	80 strips	E1730-9000	24,08 €*
Strip Caps fo	r 1.2 ml Tubes		
8-Strip Caps	125 strips	E1702-8400	15,98 €*
8-Strip Caps (Sterile)	125 strips	E1702-8410	20,48 €*
12-Strip Caps	80 strips	E1703-9400	19,80 €*
12-Strip Caps (Sterile)	80 strips	E1703-9410	23,63 €*
	Individual Tubes in Rack Individual Tubes in Rack, Sterile Individual Tubes in Rack, Sterile Individual Reaction Tubes only 1.2 ml 8-Strip 8-Strip Tubes in Rack 8-Strip Tubes in Rack, Sterile 8-Strip Reaction Tubes only 1.2 ml 12-Strip 12-Strip Tubes in Rack 12-Strip Tubes in Rack 12-Strip Tubes in Rack 8-Strip Tubes in Rack 12-Strip Tubes in Rack	DescriptionPack Size1.2 ml Individual Tubes, RacketIndividual Tubes in Rack10 racksIndividual Tubes in Rack, Sterile10 racksIndividual Reaction Tubes only1000 tubes1.2 ml 8-Strip Tubes, Racked8-Strip Tubes in Rack10 racks8-Strip Tubes in Rack, Sterile10 racks8-Strip Reaction Tubes only125 strips12-Strip Tubes in Rack10 racks12-Strip Tubes in Rack, Sterile10 racks12-Strip Tubes in Rack, Sterile10 racks12-Strip Tubes in Rack, Sterile10 racks12-Strip Caps for 1.2 ml Tubes8-Strip Caps125 strips8-Strip Caps (Sterile)125 strips12-Strip Caps80 strips	1.2 ml Individual Tubes, Racked Individual Tubes in Rack 10 racks E1750-0400 Individual Tubes in Rack, Sterile 10 racks E1750-0410 Individual Reaction Tubes only 1000 tubes E1710-0000 1.2 ml 8-Strip Tubes, Racked 8-Strip Tubes in Rack 10 racks E1760-8400 8-Strip Reaction Tubes only 125 strips E1720-8000 1.2 ml 12-Strip Tubes, Racked 12-Strip Tubes in Rack 10 racks E1770-9400 12-Strip Tubes in Rack, Sterile 10 racks E1770-9410 12-Strip Reaction Tubes only 80 strips E1730-9000 Strip Caps 8-Strip Caps 125 strips E1702-8400 8-Strip Caps (Sterile) 125 strips E1702-8410 12-Strip Caps 80 strips E1703-9400

₩ starlab.click/tubes

Reagent reservoirs.

*Offer prices, additional discounts not applicable. Valid until 26.04.2024.

REUSABLE POLYPROPYLENE RESERVOIR Description Article No. Pack Size **OFFER*** 50/60 ml StarTub, PP Reservoir. Non-sterile. V-shaped bottom. 10 S4026-5806 12,00 €* Reusable, Autoclavable,

SINGLE-USE PVC RESERVOIR				
Harry Harry	Description	Pack Size	Article No.	OFFER*
	55 ml StarTub, PVC Reservoir. Non-sterile. V-shaped bottom. Made from clear PVC.	10	E2310-1000	18,00 €*

SINGLE-USE POLYSTYRENE RESERVOIRS				
Description	Pack Size	Article No.	OFFER*	
25 ml StarTub, PS Reservoir. Sterile. Individually wrapped. V-shaped bottom. Graduations.	50	E2310-1025	42,25 €*	
55 ml StarTub PS reservoir. Sterile. Individually wrapped. V-shaped bottom. Graduations.	50	E2310-1010	38,50 €*	
100 ml StarTub PS Reservoir. Sterile. Individually wrapped. Trough-shaped bottom. Graduations. Corner pour spout.	50	E2310-1100	43,25 €*	

See our full range here!

☐ starlab.click/reservoirs
☐



PCR plates & plate seals.

CERTIFIED PCR-INHIBITOR FREE

- Ultra-thin, uniform wells for optimal heat transfer and high reaction efficiency
- Raised well edges help prevent sample evaporation when sealed
- Plates certified RNase, DNase, DNA & PCR Inhibitor free. Tested Pyrogen free.

Universal cycler compatibility



96-Well Non-Skirted Universal PCR Plate compatible with most 0.2 ml cyclers (excl. FAST)

Art. No. E1403-1200 OFFER: 18,75 €* / 10 plates

Great for ABI® FAST® systems!



96-Well Semi-Skirted PCR Plate

Low Profile (for FAST® Systems)

Art. No. E1403-7700

OFFER: 21,75 €* / 10 plates



Great for ABI® cyclers!



96-Well Semi-Skirted PCR Plate

Standard height

Art. No. I1402-9700

OFFER: 19,20 €* / 10 plates

Plate for Roche Lightcycler®



96-Well Semi-Skirted White PCR Plate for LightCycler® 96 or LightCycler® 480

Art. No. I1402-9909

OFFER: 20,40 €* / 10 plates

BUY ONLINE

starlab.click/pcrconsumables

*Offer prices, additional discounts not applicable. Valid until 26.04.2024.

Does our plate fit your cycler?



Starlab PCR plates have a precise fit for most thermal cyclers. Check your cycler on our compatibility chart.

starlab.click/pcr-compatibility

Great for ABI® cyclers!



96-Well Semi-Skirted PCR Plate

Straight edges

Art. No. I1402-9800

OFFER: 20,40 €* / 10 plates

Great for BioRad® cyclers!



96-Well Skirted PCR Plate

Low Profile

Art. No. E1403-5200 OFFER: 19,50 €* / 10 plates

See our PCR Plate Sealing Guide

starlab.click/pcrsealsguide



Aluminium StarSeal (PCR)

100 self-adhesive seals

Art. No. E2796-9792

OFFER: 42,90 €*



Polyolefin StarSeal (PCR)

100 self-adhesive seals

Art. No. E2796-9793 OFFER: 42,90 €*



X-Clear Advanced Polyolefin StarSeal (qPCR)

100 self-adhesive seals

Art. No. E2796-9795

OFFER: 130,00 €*



Polyolefin Sealing Film (qPCR)

100 self-adhesive seals

Art. No. E2796-9895

OFFER: 200.20 €*



Plate seals.

IDEAL FOR STORAGE!

Clear, self-adhesive polyester plate seals ideal for long-term storage and incubation.

■ RNase, DNase and Pyrogen free

■ Temperature use: -40 °C to 40 °C

Seal thickness: 50 μm

Sterile option

Polyester Seal (50 μm) Art. No. E2796-0101

OFFER: 47 €* / 100 seals

Polyester Seal (50 μm) Sterile Art. No. E2796-0102

OFFER: 64 €* / 100 seals



starlab.click/platesealing

Plate centrifuge.

YOU SPIN ME ROUND!

Sealed plates are loaded vertically and surface tension keeps samples in place so there is no danger of spilling samples. The Plate Centrifuge accepts all plates up to 25 mm high and is less than quarter the size of most plate centrifuges.

FEATURES

- 2500 rpm / 500 x g
- ▶ Fits all plates up to 25 mm high
- Rapid start to stop

USE ME FOR:

- Safe and quick spin down of droplets and condensation
- Spinning before and after thermal cycling to increase PCR yield

Plate Centrifuge.

Art. No. N2631-0008

OFFER: 766 €*

*Offer prices, additional discounts not applicable. Valid until 26.04.2024.





starlab.click/platecentrifuge



The minicentrifuge.

QUICK START just close the lid to start spinning!



AHEAD OF THE CURVE!

Ideal for all your quick spin-downs, the Starlab Minicentrifuge is not all about style, it also has substance! This is a compact and solid unit that delivers a quiet, smooth operation and a maximum speed of 6,500 rpm.


- Quick start on lid closure
- Magnetic lid lock
- One-touch STOP for fast braking and automatic lid release





Tube Rotor holds up to eight microcentrifuge tubes or eight conical or skirted screw cap microtubes





PCR Strip Rotor holds two, 0.2 ml 8-strips or $32 \times 0.2 \text{ ml single tubes}$.



starlab.click/minicentrifuge

QUICK & EASY TOOL FREE ROTOR EXCHANGE

> ONE-TOUCH STOP for fast braking and lid release

MAGNETIC LID LOCK aids easy, one-handed operation



WATCH VIDEO

Minicentrifuge.

Art. No. S8030-5000

OFFER: 190 €*

*Offer prices, additional discounts not applicable. Valid until 26.04.2024.



Passionate for science.

Passionate for sustainability in science.

TipOne® REFILLS

As the next step in our journey to make science more sustainable, the **clear plastic packaging** on all TipOne® refills is now being made from **100% recycled plastic.**

Nothing else has changed!

- the same trusted tips
- the same reliable quality and usability
- all plastic packaging can still be recycled in the same way

Appearance wise, there is little difference betwen the old and new packaging, so look out for the labels on the lines already converted to the recycled plastic.

#EcolutionMovement



100% RECYCLED PLASTIC

Starlab GmbH

Neuer Höltigbaum 38 22143 Hamburg Germany

T: +49 (0)40 675 99 39 0 F: +49 (0)40 675 99 39 20 info@starlab.de www.starlab.de

Starlab FRANCE SARL

30 Rue Jean Rostand 91400 Orsay France

T: +33 (0)1 60 13 71 70 F: +33 (0)1 69 41 48 65 info@starlab.fr www.starlab.fr

Starlab S.r.l.

Via Roberto Bracco 6 20159 Milano Italy

T: +39 (0)2 7020 1040 info@starlab.it www.starlab.it

