

**Your PCR can't wait,
so you grab a deal now.
You are passionate
for science.**

SPECIAL OFFERS

Consumables & Equipment for PCR

Valid until 25.04.2025



WATCH VIDEO

Passionate for Science.
Starlab inspires and
supports the passion
driving your research.



Passionate for science.

 **starlab**

PCR plates.

CERTIFIED PCR-INHIBITOR FREE

- Ultra-thin, uniform wells for optimal heat transfer and high reaction efficiency
- Raised well edges help prevent sample evaporation when sealed
- Plates certified RNase, DNase, DNA & PCR Inhibitor free. Tested Pyrogen free

Does our plate fit your cycler?



Starlab PCR plates have a precise fit for most thermal cyclers. Check your cycler on our compatibility chart.

starlab.click/pcr-compatibility

PCR Plates

FROM*: 2,00 € per plate

starlab.click/pcrplates

Universal cycler compatibility 3 options!



96-Well Non-Skirted PCR Plate
compatible with most 0,2 ml cyclers (excl. FAST)

- ▶ **Art. No. E1403-1200**
10 plates. List: 38,50€ **OFFER: 20 €***
- ▶ **Art. No. E1403-1209 (Opaque white plate)**
10 plates. List: 48,00€ **OFFER: 26 €***
- ▶ **Art. No. E1403-0100 (Cuttable plate)**
10 plates. List: 41,00€ **OFFER: 21 €***

Reduce contamination



96-Well Non-Skirted PCR Plate
with elevated wells

- ▶ **Art. No. E1403-2200**
10 plates. List: 47,00€ **OFFER: 25 €***

Plate for TAQMan®



96-Well Semi-Skirted PCR Plate
Raised Rim (for TAQMan®)

- ▶ **Art. No. E1403-8200**
10 plates. List: 38,50€ **OFFER: 20 €***

Great for ABI® cyclers!



96-Well Semi-Skirted PCR Plate
Standard height

- ▶ **Art. No. I1402-9700**
10 plates. List: 49,00€ **OFFER: 20 €***

BARCODED skirted and semi-skirted plates also available. Contact Starlab for details.

Great for ABI® FAST® systems!



96-Well Semi-Skirted PCR Plate
Low Profile (for FAST® Systems)

- ▶ **Art. No. E1403-7700**
10 plates. List: 44,50€ **OFFER: 23 €***

Plate for Roche Lightcycler®



96-Well Semi-Skirted PCR Plate
for LightCycler® 96 or LightCycler® 480

- ▶ **Art. No. I1402-9909 (Opaque white plate)**
10 plates. List: 52,50€ **OFFER: 22 €***

Universal cycler compatibility



96-Well Semi-Skirted PCR Plate
Straight edges

- ▶ **Art. No. I1402-9800**
10 plates. List: 52,50€ **OFFER: 22 €***

Great for BioRad® cyclers!



96-Well Skirted PCR Plate
Low Profile

- ▶ **Art. No. E1403-5200**
10 plates. List: 40,00€ **OFFER: 20,50 €***



***OFFER REF: Q25000746809**
must be stated on your PO (excl. eShop).
Offer valid to 25.04.2025.

Plate sealing.

♥ SEAL FOR PCR



Aluminium StarSeal, 36 µm
100 self-adhesive seals

Art. No. E2796-9792

List: 67,50€ **OFFER: 39,50 €***

♥ SEALING FOR qPCR



X-Clear Advanced Polyolefin StarSeal, 50 µm
100 self-adhesive seals

Art. No. E2796-9795

List: 205,00€ **OFFER: 120 €***



X-Clear Polyolefin Sealing Film, 96 µm (incl. plastic seal applicator!)
100 self-adhesive seals

Art. No. E2796-9895

List: 315,50€ **OFFER: 192 €***



8-Strip PCR Caps, Flat (Xtra-Clear)
125 Strips

Art. No. I1400-0900

List: 25,00€ **OFFER: 12,50 €***


starlab.click/pcrseals



See our
PCR Plate
Sealing
Guide

starlab.click/pcrsealsguide

PCR strips.



OPTIMAL HEAT TRANSFER, OPTIMUM CERTIFICATION

- Ultra-thin, uniform walls ensure optimal heat transfer.
- Supplied in resealable bags to minimise contamination risk.
- Certified RNase, DNase, DNA & PCR Inhibitor free. Tested Pyrogen free.

8-Strip Tube & Caps, Combo Packs



Description	Pack Size	Article No.	List	OFFER*
0,2 ml 8-Strip PCR Tubes and 8-Strip PCR Caps, Domed	125 pairs	I1402-6700	130,50€	69 €
0,2 ml, 8-Strip PCR Tubes and 8-Strip PCR Caps, Flat (X-Clear)	125 pairs	I1402-5700	131,00€	72 €

Non-Flex 8-Strip PCR Tubes with Individually Attached Flat Caps, Xtra-Clear



Description	Pack Size	Article No.	List	OFFER*
0,2 ml 8-Strip Non-Flex PCR Tubes, Individually Attached Flat Caps (X-Clear)	120 strips	I1402-3700	117,50€	60 €

The rigid, triple-link between the tubes provides increased stability for easier & faster handling

Rotor-Gene® Q strips.



Discs and tube strips for use with the
Qiagen Rotor-Gene® Q Real-Time Rotary Analyzer

Description	Pack Size	Article No.	List	OFFER*
0,1 ml 4-Strip Rotor-Gene® Style Tubes & Caps	250 pairs	I1402-0400	181,00€	88 €

starlab.click/pcrstrips




***OFFER REF: Q25000746809**
must be stated on your PO (excl. eShop).
Offer valid to 25.04.2025.

5.0 ml Tubes

**WHEN 2,0 ml IS TOO SMALL AND 15 ml TOO BIG,
WE HAVE THE PERFECT ANSWER!**

- Ultra clear tubes for excellent sample visibility
- Same diameter as 15 ml tubes, so will fit your existing 15 ml bench racks
- Certified RNase, DNase, DNA & Pryrogen free
- Use down to -196 °C (gas phase only)



5.0 ml Centrifuge Tube with Flip Cap
2 x 100

Art. No. E1450-5000

List: 70,00€ **OFFER: 42 €***

5.0 ml Centrifuge Tube with Flip Cap, (Sterile)
10 x 20

Art. No. E1450-5010

List: 81,00€ **OFFER: 49 €***

5.0 ml Centrifuge Tube with Screw Cap, (Sterile)
10 x 20

Art. No. E1450-5050





List: 91,50€- **OFFER: 55 €***

starlab.click/5mltubes


PCR sample racks.


StarChill CHILLER RACKS

The centre of the rack changes colour when frozen and then gradually turns back to the original colour (~7 °C) as it warms up, to provide you with a visual indication of the rack's temperature.

Description	Pack Size	Colour	Article No.	List	OFFER*
 StarChill PCR For use with one 0,2ml PCR Plate or 0,2ml PCR Tubes or Strips	2	Pink/Purple	E2396-2004	169,00€	66 €
		Yellow/Green	E2396-2006	161,00€	63 €
 StarChill MCT For use with 0,5, 1,5 and 2,0 ml Microcentrifuge or Screw Cap Tubes	1	Pink/Purple	E2396-5634	145,50€	57 €
		Yellow/Green	E2396-5636	145,50€	57 €
 StarChill SBS For use on AUTOMATED SYSTEMS with one 0,2 ml PCR Plate or 0,2ml PCR Tubes or Strips	2	Pink/Purple	E2396-5654	186,00€	73 €
		Yellow/Green	E2396-5656	186,00€	73 €
 StarChill qPCR For 64 x 0,1 ml Rotor-Gene® tubes, 16 x 0,2 ml tubes and 4 x 0,5 ml tubes	2	White/Blue	E2396-3008	186,00€	73 €

Description	Pack Size	Colour	Article No.	List	OFFER*
 StarPCR Work Racks For use with one 0,2ml PCR Plate or 0.2ml PCR tubes or strips	5	Mixed (Blue, Green, Red, Violet, Yellow)	E2396-5240	68,00€	27 €

Description	Pack Size	Colour	Article No.	List	OFFER*
 Double-Sided PCR Work Rack For use with 168 x 0,2 ml tubes and 12 x 1,5 ml tubes on one side, 40 x 0,5 ml tubes & 12 x 1,5 ml tubes on the other	5	Mixed (Blue, Green, Orange, Pink, Yellow)	E2396-0599	179,00€	70 €

Description	Pack Size	Colour	Article No.	List	OFFER*
 PCR Workstation Rack Rack with lid for use with one 0,2 ml PCR Plate or 0,2 ml PCR tubes or strips	5	Mixed Neon (Blue, Green, Orange, Pink, Yellow)	E2396-9649	70,00€	27,50 €

***OFFER REF: Q25000746809**
must be stated on your PO (excl. eShop).
Offer valid to 25.04.2025.

starlab.click/pcrracks

Reagent reservoirs.

SEE OUR FULL RANGE

starlab.click/reservoirs

Reusable White Polypropylene Reagent Reservoir



Description	Pack Size	Article No.	List price	OFFER*
50/60 ml StarTub, PP Reservoir. Non-sterile. V-shaped bottom. Reusable. Autoclavable.	10	S4026-5806	20,50€	14 €

Single-Use Clear PVC Reagent Reservoir



Description	Pack Size	Article No.	List price	OFFER*
55 ml StarTub, PVC Reservoir. Non-sterile. V-shaped bottom. Made from clear PVC.	100	E2310-1000	37,00€	18,50 €

Single-Use White Polystyrene Reagent Reservoirs

Description	Pack Size	Article No.	List price	OFFER*
25 ml StarTub, PS Reservoir. Sterile. Individually wrapped. V-shaped bottom. Graduations.	50	E2310-1025	79,00€	39,50 €
55 ml StarTub PS reservoir. Sterile. Individually wrapped. V-shaped bottom. Graduations.	50	E2310-1010	86,50€	43,25 €
100 ml StarTub PS Reservoir. Sterile. Individually wrapped. Trough-shaped bottom. Graduations. Corner pour spout.	50	E2310-1100	88,50€	44,25 €
100 ml StarTub, PS Reservoir. Non-sterile. Bulk packed. Trough-shaped bottom. Graduations. Corner pour spout.	100	E2310-2000	90,90€	66 €

NEW!



***OFFER REF: Q25000746809**
must be stated on your PO (excl. eShop).
Offer valid to 25.04.2025.

StarGuard® Comfort.



WATCH VIDEO



StarGuard® Comfort have snug fitting fingers and are 245 mm long to provide full cover over the wrist.

***ONLY: 5,32 €**
per 100 gloves (XL: 5,78 €)



NITRILE GLOVES WITH AN EXCEPTIONAL FIT

- Lightweight, ergonomic glove
- Provides dexterity and tactile sensitivity
- Helps reduce hand fatigue caused by gloves
- Combines comfort and chemical protection

StarGuard® Comfort Nitrile Gloves

Description	Pack Size	Article No.	List	OFFER*
StarGuard® Comfort, X-Small	10 x 250	SG-C-XS	224€	133 €
StarGuard® Comfort, Small	10 x 250	SG-C-S	224€	133 €
StarGuard® Comfort, Medium	10 x 250	SG-C-M	224€	133 €
StarGuard® Comfort, Large	10 x 250	SG-C-L	224€	133 €
StarGuard® Comfort, X-Large	10 x 230	SG-C-XL	224€	133 €

starlab.click/sg-comfort

ErgoOne® Pipettes.

ACCURACY AND PRECISION YOU CAN TRUST!



PRODUCT FEATURES

- Ergonomic and lightweight
- Minimal operating forces
- Volume lock
- Universal tip fit (excl. 5,0 ml)
- Colour coded

ErgoOne®

from: 133 €*

 starlab.click/ergoone



WATCH VIDEO

Single-Channel Manual Pipettes

0,1 – 2,5 µl (red) S7100-0125 €254,00 133 €*	0,5 – 10 µl (red) S7100-0510 €254,00 133 €*	2 – 20 µl Micro (red) S7100-0221 €254,00 133 €*	2 – 20 µl (yellow) S7100-0220 €254,00 133 €*
10 – 100 µl (yellow) S7100-1100 €254,00 133 €*	20 – 200 µl (yellow) S7100-2200 €254,00 133 €*	100 – 1.000 µl (blue) S7110-1000 €254,00 133 €*	500 – 5.000 µl (violet) S7150-5000 €254,00 133 €*

8-Channel Manual Pipettes

0,5 – 10 µl (red) S7108-0510 €732,50 388 €*	10 – 100 µl (yellow) S7108-1100 €732,50 388 €*	30 – 300 µl (green) S7108-3300 €732,50 388 €*
--	---	--

12-Channel Manual Pipettes

0,5 – 10 µl (red) S7112-0510 €875,00 401 €*	10 – 100 µl (yellow) S7112-1100 €873,00 446 €*	30 – 300 µl (green) S7112-3300 €873,00 446 €*
--	---	--



***OFFER REF: Q25000746809**
must be stated on your PO (excl. eShop).
Offer valid to 25.04.2025.

GuardOne® Workstation.

FROM: 2.652 €*



Contents not included.



WATCH VIDEO

AVOID CROSS-CONTAMINATION DURING AMPLIFICATION OF DNA/RNA

Protect your valuable samples and reagents¹⁾ to save time and costs by avoiding the need to repeat contaminated work.

Starlab's GuardOne® Workstation provides

- An ideal clean and controlled environment that can help you achieve accurate PCR results
- Convenience! All functions are controlled via the glove compatible, easy to use touchscreen
- Quick, tool-free HEPA filter exchange and UV lamp exchange

Perfect for your PCR / qPCR applications or any other application that requires a controlled, clean environment.



DEAD AIR 32"

S8040-1000

List: 4.398,50-€*

OFFER: 2.652 €*



DEAD AIR 48"

S8040-1001

List: 5.579,50-€*

OFFER: 3.420 €*



LAMINAR FLOW 32"

S8040-0000

List: 5.226,50-€*

OFFER: 3.150 €*



LAMINAR FLOW 48"

S8040-0001

List: 6.966,00-€*

OFFER: 4.200 €*



www.starlab.click/guardone



- One-touch QuickStart / QuickStop settings
- Efficient UV_{245 nm} decontamination with timer function
- Efficient, class H14 HEPA filter with pre-filter provides ISO class 5 working environment²⁾
- Real-time airflow measurement²⁾
- Life-span monitoring of consumables
- Virtually maintenance-free
- Easy to clean
- Easy, tool-free 15-minute assembly

¹⁾GuardOne® is not designed to protect the user.

²⁾Laminar flow models only.



***OFFER REF: Q25000746809**

must be stated on your PO (excl. eShop).
Offer valid to 25.04.2025.


Starlab Minicentrifuge.

IDEAL FOR YOUR QUICK SPIN DOWNS!

The Starlab Minicentrifuge is a compact and solid unit that delivers a quiet, smooth operation and a maximum speed of 6.500 rpm.

- Quick start on lid closure
- Magnetic lid lock
- One-touch STOP for fast braking and automatic lid release

SUPPLIED WITH TWO ROTORS



TUBE ROTOR holds up to eight microcentrifuge tubes or eight conical or skirted screw cap microtubes



PCR STRIP ROTOR holds two, 0,2 ml 8-strips or 32 x 0,2 ml single tubes



MAGNETIC LID LOCK aids easy, one-handed operation

QUICK START
just close the lid
to start spinning!

**QUICK & EASY
TOOL FREE
ROTOR EXCHANGE**

**ONE-TOUCH
STOP** for fast
braking and
lid release



WATCH VIDEO

Minicentrifuge.

Art. No. S8030-5000

List: 278,50-€

OFFER: 194 €*

 starlab.click/minicentrifuge



***OFFER REF: Q25000746809**
must be stated on your PO (excl. eShop).
Offer valid to 25.04.2025.

Starlab Vortex Mixer.

THE TOUCH OF A VORTEXER. THE POWER OF A MIXER

Ideal for the gentle mixing of samples right up to vigorous vortex mixing and the resuspending of cells or chemical pellets.

- Precise vortexing speeds
- Speed range: 300 to 3.500 rpm
- Unique, glove compatible touch screen
- One-touch pre-set options for speed / time
- Touch or continuous operation
- End of run audible alarm
- Power-saving sleep mode

TOUCH MODE

Minimal pressure required to maintain vibration. Mix single or multiple tubes up to 3.500 rpm.

GLOVE COMPATIBLE TOUCHSCREEN for precise speed and time settings



SPEED RANGE
300 – 3.500 RPM



WATCH VIDEO

Vortex Mixer.

Art. No. S8010-0000

List: 409,00€

OFFER: 285 €*

 starlab.click/vortexmixer



Adapter for Tubes

For mixing up to 15 x 1,5/2,0 ml tubes up to 1.500 rpm

Art. No. S8010-0011

List: 28,50€

OFFER: 16,50 €*



Adapter for Microplate

For mixing one microplate up to 1.500 rpm

Art. No. S8010-0012

List: 17,50€

OFFER: 10,50 €*

Tubes and plate are not included with the adapters.



***OFFER REF: Q25000746809**

must be stated on your PO (excl. eShop). Offer valid to 25.04.2025.

Thermomixer-Mixer HC.

The tiny footprint of the Mixer HC fits on any bench to provide a convenient, all-in-one solution that simplifies PCR workflows:



WATCH VIDEO

PRE-PCR STEPS

- Isolation of genomic DNA and RNA
- Plasmid DNA purification
- Reverse Transcription of mRNA
- Mixing of PCR setups

POST-PCR STEPS

- Purification of DNA from agarose gels
- Cloning of DNA fragments (Restriction digestion, A-tailing, Phosphorylation, Ligation, Transformation into host cells, Cultivation of host cells)

FUNCTIONAL!
INTERVAL MIX: alternate between mixing & resting
PROGRAM MODE: create your own protocols
SHORT MIX: instant mixing at the push of a button



DIGITAL DISPLAY
 for precise adjustment.
 Clear menu navigation and
 simple parameter setting

Thermomixer-Mixer HC

Art. No. **S8012-0000** (excludes Thermoblock)

List: 2.678,50 €

OFFER: 1.470 €*

starlab.click/mixer-hc

HEATING OR COOLING
 Precise temperature control
 up to 99°C or down to
 13°C below room temperature

MIXING
 Up to 1.500 rpm

TIMING
 1 min to 99:59 hours

USE WITH ALL POPULAR TUBES AND PLATES



24 x 0,5 ml
Microcentrifuge Tubes¹⁾
S8012-0010
OFFER: 386 €*



24 x 1,5 ml
Microcentrifuge Tubes¹⁾
S8012-0011
OFFER: 386 €*



24 x 2,0 ml Microcentrifuge
Tubes, Round Bottoms¹⁾
S8012-0012
OFFER: 386 €*



24 x 2,0 ml Microcentrifuge
Tubes Ø 12 mm, V-Bottoms
S8012-0013
OFFER: 478 €*



8 x 5,0 ml
Preparation Tubes
S8012-0014
OFFER: 386 €*



8 x 15 ml
Conical Centrifuge Tubes
S8012-0015
OFFER: 386 €*



4 x 50 ml
Conical Centrifuge Tubes
S8012-0016
OFFER: 386 €*



24 x 1,5/2,0 ml
Cryo Tubes
S8012-0017
OFFER: 386 €*



Micro/Deepwell Plate,
with Lid
S8012-0018
OFFER: 478 €*



Adapter Plate 96 x 0,2 ml
Tubes or a PCR plate
(Use with S8012-0018)
S8012-0019
OFFER: 131 €*

¹⁾ Includes ISO-Rack. Tubes not included.
 ISO Racks also available to purchase separately.
 See website for details.



***OFFER REF: Q25000746809**
 must be stated on your PO (excl. eShop).
 Offer valid to 25.04.2025.

Clearance deals!

ONLY AVAILABLE WHILST STOCK LASTS!

Microflex Supreno® Blue Nitrile Gloves



Glove Size	Pack Size	Article No.	List	OFFER*
Medium	10 x 100	SU-INT-M	212,50€	89 €
Large	10 x 100	SU-INT-L	212,50€	89 €
X-Large	10 x 100	SU-INT-XL	212,50€	89 €

PPE Complex Design Category III.

Ansell TouchNTuff® Green Nitrile Gloves



Glove Size	Pack Size	Article No.	List	OFFER*
Small	10 x 100	92-600-S	310,50€	167 €
Medium	10 x 100	92-600-M	310,50€	167 €
Large	10 x 100	92-600-L	310,50€	167 €
X-Large	10 x 100	92-600-XL	310,50€	167 €

PPE Complex Design Category III.

StarGuard® Sensitive Blue Nitrile Gloves



Glove Size	Pack Size	Article No.	List	OFFER*
X-Large	10 x 200	SG-N-XL	202,00€	115 €

PPE Complex Design Category III.
Low Dermatitis Potential: ideal for users with sensitive skin.
Contains no contaminants to interfere with experiments.



250 µl Sterile Tip for Biomek FX

Art. No. E1076-0410
10 x 96-Tip Racks
List: 109,50€
OFFER: 43 €*



50 µl Conductive Sterile Filter Tip for Tecan/Qiagen

Art. No. E1172-1418
10 x 96-Tip Hanging Wafers
List: 432,50€
OFFER: 169 €*

96-Well Semi-Skirted PCR Plates

Great for ABI® cyclers!



Colour	Pack Size	Article No.	List	OFFER*
Blue	10	I1402-9701	58,50€	14,75 €
Green	10	I1402-9702	67,50€	16,90 €
Red	10	I1402-9704	67,50€	16,90 €
Violet	10	I1402-9705	67,50€	16,90 €
Yellow	10	I1402-9706	67,50€	16,90 €

384-Well Skirted qPCR Plate, Opaque White

384-Well ABI® / Universal Plate



Colour	Pack Size	Article No.	List	OFFER*
White	10	E1042-3849	69,50€	21,50 €

0.1 ml 8-Strip qPCR Tubes, Ind-Attached Flat Caps (X-Clear)



Colour	Pack Size	Article No.	List	OFFER*
Clear	120 strips	A1402-3800	122€	42,00 €

0.1 ml 8-Strip PCR Tubes, White AND 8-Strip Flat Caps (X-Clear)



Colour	Pack Size	Article No.	List	OFFER*
White (tubes) Clear (caps)	125 pairs	I1409-4700	138€	44,00 €

0.1 ml 8-Strip PCR Tubes



Colour	Pack Size	Article No.	List	OFFER*
Clear	125 strips	I1402-3600	120€	60,00 €
White	125 strips	I1402-3609	126€	63,00 €

All PCR products certified RNase, DNase, DNA & PCR Inhibitor free. Tested Pyrogen free.



***OFFER REF: Q25000746809**

must be stated on your PO (excl. eShop).
Offer valid to 25.04.2025 / while stock lasts.

TipOne® is coming home!

We have exciting news in answer to our customers' most asked question:

Q: What are you doing with the plastic from my TipOne® recycling?

A: We are closing the loop!

Starting on 2nd January, we began the process of moving TipOne® manufacturing from the USA back to Germany. The move to Germany will bring numerous economic and sustainability benefits.

The manufacturing move also signals the next milestone on the sustainable journey for TipOne®. It means we will now start to **use the plastic returned via the TipOne® recycling service to make the base of our TipOne® racks!** We are proud that we can say that we have started to close the loop. We could not have achieved this milestone without you, our customers, and your continued support and use of the TipOne® recycling collection service.

For more information on the incredible sustainability journey of TipOne® so far:

starlab.click/tipone

Are you a TipOne® tip user but not using the Recycling Service? Sign up here:

starlab.click/recycling



T: +49 (0)40 675 99 39 0
info@starlab.de
www.starlab.de



T: +33 (0)1 60 13 71 70
info@starlab.fr
www.starlab.fr



T: +39 (0)2 7020 1040
info@starlab.it
www.starlab.it

Passionate for science.

