

Cell culture challenges, but
you continue to grow.
**You are passionate
for science.**

SPECIAL OFFERS

Cell Culture Consumables & Equipment

Valid until 27.08.2025



WATCH VIDEO

Passionate for Science.
Starlab inspires and
supports the passion
driving your research.

Passionate for science.



CytoOne® cell culture ware.

YOUR CELLS DESERVE THE BEST

Designed for optimal performance and manufactured to the highest quality standards. Outfit your cells for true performance with CytoOne® cell culture products.

CytoOne® products are Gamma sterilised (SAL 10⁻⁵).

CytoOne® Dishes.

CytoOne® dishes feature stacking rings with vented design for thermal equalisation. Vented lids improve gas exchange.

FROM 41,04 €*



Size	Pack Size	Art. No.	List	OFFER*
TC Treated Dishes				
35 x 10 mm	30 x 10	CC7682-3340	75,50€	43,04 €
60 x 15 mm	30 x 10	CC7682-3359	72,00€	41,04 €
100 x 20 mm	30 x 10	CC7682-3394	158,50€	90,35 €
150 x 20 mm	12 x 5	CC7682-3614	93,50€	53,30 €

starlab.click/dishes

Convenient, small pack sizes and easy to open packaging!

CERTIFIED



CytoOne® Flasks.

All flasks are 100% pressure tested for leak-free assurance. Features stacking rims with vented design for thermal equalisation.



FROM 52,44 €*

Flasks with Vented Filter Caps

Size	Pack Size	Art. No.	List	OFFER*
TC Treated Flasks, Vented Filter Caps				
T-25	30 x 10	CC7682-4825	397,00€	226,29 €
T-75	20 x 5	CC7682-4875	235,00€	133,95 €
T-150	8 x 5	CC7682-4815	152,50€	86,93 €
T-225	5 x 5	CC7682-4822	147,50€	84,08 €



CytoOne® vented filter caps have a 0.2 µm hydrophobic membrane for gas exchange and to exclude bacteria and fungal contaminants

Flasks with 2-Position Plug Seal Caps

Size	Pack Size	Art. No.	List	OFFER*
TC Treated Flasks, 2-Position Plug Seal Caps				
T-150	8 x 5	CC7682-4415	124,50€	70,97 €
T-225	5 x 5	CC7682-4225	112,50€	64,13 €
Non-Treated Flasks, 2-Position Plug Seal Caps				
T-75	20 x 5	CC7672-4175	151,00€	86,07 €
T-225	5 x 5	CC7672-4225	92,00€	52,44 €



CytoOne® two-position plug seal caps allow you to choose secure manual venting or an airtight seal

starlab.click/flasks



CERTIFIED



CytoOne® Plates.

FROM 40,19 €*

Wells	Pack Size	Art. No.	List	OFFER*
TC Treated Plates				
6-Well	50	CC7682-7506	70.50 €	40,19 €
12-Well	50	CC7682-7512	94.00 €	53,58 €
24-Well	50	CC7682-7524	82.00 €	46,74 €
48-Well	50	CC7682-7548	101.50 €	57,86 €
96-Well	50	CC7682-7596	98.00 €	55,86 €

starlab.click/multiwellplates




The complete 360° chimney well design on all CytoOne® plates prevents cross-contamination and ensures uniform temperature transfer across all wells, resulting in a more uniform cell layer. Watch our video to see how CytoOne® plates are designed to prevent the edge effect so you can use all 96-wells of a plate.



WATCH VIDEO

PLATE SEALING



Porous, 50 µm self-adhesive film that allows gas exchange for cellular and bacteria cultivation. Sterile. RNase, DNase & Pyrogen free. Temp use: 4° to 40° C.

Breathable Film

Pack size: 100

Art. No. E2796-3005

List: 157,50 €

91,35 €*

Breathable Film (Sterile)

Pack size: 50

Art. No. E2796-3015

List: 125,00 €

72,50 €*



Clear, 50 µm self-adhesive polyester seal ideal for long-term storage & incubation. Temp use: -4.0 °C to 120 °C.

Polyester Seal

Pack size: 100

Art. No. E2796-0101

List: 79,00 €

45,82 €*

Polyester Seal (Sterile)

Pack size: 100

Art. No. E279-0102

List: 107,50 €

62,35 €*

starlab.click/plateseals



Cell Scrapers.

All cell scrapers feature flexible ABS handles, and are sterile and individually-wrapped.



220 mm straight handle, 11 mm fixed blade

Ideal for stem cells due to reduced shearing from gentle silicone blade. Use to gently lift cells from surface rather than scrape. Use for plates and dishes.

Art. No. CC7600-0220

Pack Size: 100

List: 253,00 € **146,74 €***



250 mm angled handle, 10 mm pivoting blade

Ideal for stem cells due to reduced shearing from gentle silicone blade. Flexibility allows you to reach all corners of your small flask.

Art. No. CC7600-0250

Pack Size: 100

List: 241,00 € **139,78 €***



300 mm straight handle, 20 mm pivoting blade

Large pivoting polyethylene blade for easy access and use in large flasks.

Art. No. CC7600-0300

Pack Size: 150

List: 204,50 € **118,61 €***

starlab.click/cellscrapers



***OFFER REF: Q25000794669**

must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

Aspirator.

ONLY: 1.448,45 €*



See how
to add
a second
hand
controller.

- The powerful pump can operate two hand controllers without any loss of performance
- Chemically-resistant, high-performance pump with a long service life
- Auto ON-OFF pump switch with low noise emission for a quieter working environment
- Ergonomic, fully autoclavable hand controller
- Intuitive, reproducible adjustment of suction power via touch panel



Aspirator

Incl. everything you need to get started! 1-year warranty.

N2400-9000

2.455,00€

1.448,45 €



Hand Controller for Aspirator

Incl. 1x Adapter for Pipettes, 1x Adapter for Pipette Tips, 1x Wall Mount, and 2.5 m Tubing.

N2400-9001

229,50€

133,11 €



Connection Kit for Second Hand Controller

N2400-9006

57,00€

33,06 €



starlab.click/aspirator

Aspiration pipettes.

- Accuracy: $\pm 2\%$
- Colour coded, paper-plastic peel-back wraps
- Gamma sterilised to SAL 10^{-6} (ISO 11137)
- Certified RNase, DNase, Pyrogen and latex free
- Non-cytotoxic and non-haemolytic USP Class VI Polystyrene

2.0 ml

5.0 ml

10 ml

Lot number and expiry date printed on every aspiration wrap

Volume	Pack Size	Art. No.	List	OFFER*
2.0 ml	500	E4861-0002	119,00€	59,50 €
5.0 ml	200	E4861-0005	55,50€	27,75 €
10 ml	200	E4861-0010	71,00€	35,50 €



starlab.click/aspirationpipettes



Aspiration pipettes are supplied in convenient dispenser-style boxes for easy access when stacked



***OFFER REF: Q25000794669**

must be stated on your PO (excl. eShop). Offer valid to 27.08.2025.

ErgoOne® FAST pipette controller.

Use with all
glass / plastic
pipettes from
1.0 to 100 ml.



FAST!
Fills a 25 ml
pipette in
less than
4 seconds.

Pressure-responsive speed
controls make the ErgoOne®
FAST ideal for cell culture
work such as resuspension
of bacteria or cell pellets,
nucleic acid precipitates
and aspiration of a cell layer.



WATCH VIDEO

ErgoOne® FAST

Art. No. S7166-0010

List: 451 €

OFFER*: 261,58 €

starlab.click/ergoonefast

See our range of serological pipettes starlab.click/serologicals


Completely charged,
the battery can
perform 2000 x 25 ml
dosing cycles. The
battery is also very
easy to replace!



StarGuard® Protect.



ONLY:
7,51 €
per 100 gloves



AQL
0.65



WATCH VIDEO

TRUSTED, RELIABLE PROTECTION

- Combines strength and durability for everyday protection
- Comfort, tactile sensitivity and consistent grip
- Excellent chemical splash protection
- Longer, 250 mm glove to cover the wrist

Description	Pack Size	Article No.	List	OFFER*
X-Small	10 x 100	SG-P-XS	129,50-€	75,11 €
Small	10 x 100	SG-P-S	129,50-€	75,11 €
Medium	10 x 100	SG-P-M	129,50-€	75,11 €
Large	10 x 100	SG-P-L	129,50-€	75,11 €
X-Large	10 x 100	SG-P-XL	129,50-€	75,11 €

starlab.click/sg-protect



***OFFER REF: Q25000794669**

must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

ErgoOne® Starter Packs

Reliability is key! Your pipette must be trusted to provide accurate and precise pipetting, day in, day out. ErgoOne® manual pipettes provide you with both the reliability you need and an ergonomic design that prevents fatigue over long periods to help keep your pipetting consistent. ErgoOne® pipettes have tighter specifications than ISO 8655, and feature a volume lock to ensure the correct volumes in your reaction.



CONTENTS	STARTERPACK	STARTERPACK42.5	STARTERPACK410
0.1 – 2.5 µl ErgoOne® Single-Ch Pipette S7100-0125		●	
0.5 – 10 µl ErgoOne® Single-Ch Pipette S7100-0510			●
2.0 – 20 µl ErgoOne® Single-Ch Pipette, Standard Cone S7100-0220	●	●	●
20 – 200 µl ErgoOne® Single-Ch Pipette S7100-2200	●	●	●
100 – 1000 µl ErgoOne® Single-Ch Pipette S7110-1000	●	●	●
TipOne® Pipette Tips (three racks of sterile tips per pack)	●	●	●
Self-Adhesive Pipette Holder (one holder per pipette supplied)	●	●	●
Laboratory Timer (one timer per pack)	●	●	●
List:	List: 445 €	List: 595 €	List: 595 €
OFFER:	378,25 €* 	505,75 €* 	505,75 €*

ErgoOne® pipettes paired with Starlab's trusted TipOne® pipette tips create the perfect pipetting partnership!

 starlab.click/starterpacks



***OFFER REF: Q25000794669**
must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

Cryogenic Vials.

FROM: 78,78 €*

- 100% polypropylene tubes
- For storage down to -196 °C; **should only be used in the gas phase of liquid nitrogen**
- Internal threaded vials have a silicone seal in the cap for the highest possible impermeability
- White writing area
- Highly visible black graduations at 100 µl increments
- Sterile (SAL 10⁻⁶)
- Certified RNase, DNase, DNA and Endotoxin free



Supplied in convenient, resealable, self-standing pouches.




EXTERNAL Threaded Tube Self-Sealing Cap, Sterile



INTERNAL Threaded Tube, Silicone Seal in Cap Sterile

CRYOGENIC VIALS (STERILE)

Images in ordering table show internal threaded tubes.



Volume:	1.0 ml	1.8 ml	1.8 ml	3.6 ml	3.6 ml	4.5 ml	4.5 ml
Base:	Skirted	Conical	Skirted	Conical	Skirted	Conical	Skirted
	E3090-6212 10 x 50 271,50 € 105,89 €* 	E3090-6321 10 x 50 292,50 € 114,08 €* 	E3090-6222 10 x 50 253,50 € 98,87 €* 	n/a	n/a	n/a	E3090-6242 6 x 50 202,00 € 78,78 €*
	E3110-6112 10 x 50 252,50 € 123,73 €* 	E3110-6121 10 x 50 283,00 € 138,67 €* 	E3110-6122 10 x 50 237,00 € 116,13 €* 	E3110-6131 8 x 50 207,50 € 101,68 €* 	E3110-6132 8 x 50 234,50 € 114,91 €* 	E3110-6141 6 x 50 216,50 € 106,09 €* 	E3110-6142 6 x 50 216,50 € 106,09 €*

starlab.click/cryovials

Cryogenic Storage Boxes.



81-place & 100-place
Polycarbonate boxes with fixed dividers. Transparent lid printed with reference grid and has a one-way fit so your inventory is safe.
130 x 130 x 53 mm.

I2381-2305	48,50 €
81-place box (x5)	22,80 €*
I2510-5005	48,50 €
100-place box (x5)	22,80 €*



81-place tall
Polycarbonate boxes with fixed dividers. Matrix printed on dividers. Lift-off transparent lid.
133 x 133 x 95 mm.

E2382-6002	149,50 €
Green (x5)	49,50 €*
E2382-6003	149,50 €
Orange (x5)	49,50 €*
E2382-6006	149,50 €
Yellow (x5)	49,50 €*

Available only whilst stock lasts!



25-place
Polycarbonate boxes with fixed dividers. Matrix printed on dividers. Lift-off transparent lid.
133 x 133 x 54 mm.

E2381-5002	76,50 €
Green (x5)	29,85 €*
E2381-5003	76,50 €
Orange (x5)	29,85 €*
E2381-5004	76,50 €
Red (x5)	29,85 €*

Available only whilst stock lasts!



***OFFER REF: Q25000794669**
must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

Erlenmeyer Flasks.

Autoclavable
and reusable
a minimum of
10 times!

 starlab.click/erlenmeyers

DESIGNED TO BE REUSED

Ideal for the cultivation of bacteria, yeast and fungi, as well as plant and animal cells in suspension.

PRODUCT FEATURES:

- Optically clear polycarbonate flasks
- Moulded graduations
- Leakproof and shatterproof
- Dual-function, vented polypropylene cap provides sterile air exchange or a tight seal
- Suitable for both aerobic and anaerobic cultures
- 0.22 µm PTFE cap membrane
- Individually-wrapped
- Sterility level (SAL 10⁻⁶)
- Certified RNase, DNase, DNA & Pyrogen free
- Autoclavable
- Reusable a minimum of 10 times

	Volume	Pack Size	Art. No.	List	OFFER*
FLAT BASE 	125 ml	24	E5000-1012	244,00 €	151,28 €
	250 ml	12	E5000-1025	161,00 €	99,82 €
	500 ml	12	E5000-1050	211,00 €	130,82 €
	1000 ml	6	E5000-1100	147,50 €	91,45 €
	2000 ml	6	E5000-1200	186,50 €	115,63 €
BAFFLED BASE 	3000 ml	3	E5000-1300	223,50 €	138,57 €
	125 ml	24	E5001-1012	244,00 €	151,28 €
	250 ml	12	E5001-1025	161,00 €	99,82 €
	500 ml	12	E5001-1050	211,00 €	130,82 €
	1000 ml	6	E5001-1100	147,50 €	91,45 €
	2000 ml	6	E5001-1200	186,50 €	115,63 €
	3000 ml	3	E5001-1300	223,50 €	138,57 €

All products are sterile (SAL 10⁻⁶) and individually wrapped.

DUAL-FUNCTION CAPS

Polypropylene caps with 0.22 µm PTFE membrane, ideal for both aerobic and anaerobic cultures. Remove the top section to allow sterile exchange, and then replace when ready. Avoid contamination risk and the hassle of switching caps and stock only one flask type for all cultures.



TOTALLY TUBULAR

Four, 29-place cube racks that hold a variety of tube sizes, in a support tray. Overall WxDxH: 93 x 93 x 48 mm



12344-1000
72,00 €
33,84 €*

COMBI MINI RACKS

80-place polypropylene racks for use with a variety of tube sizes. Ideal for use in the hood. WxDxH: 171 x 91 x 54 mm



E2345-1000	12,50 €
Natural x1	5,38 €*
E2345-1041	13,00 €
Neon Blue x1	4,55 €*
E2345-1042	13,00 €
Neon Green x1	4,55 €*
E2345-1043	13,00 €
Neon Orange x1	4,55 €*
E2345-1044	13,00 €
Neon Pink x1	4,55 €*
E2345-1046	13,00 €
Neon Yellow x1	4,55 €*
E2345-1099	64,00 €
Mixed Neon x5	28,80 €*

NOTE: Individual neon colour racks only available whilst stock lasts!



***OFFER REF: Q25000794669**

must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

Bottle top filtration.

Certified
RNase, DNase,
Pyrogen,
non-cytotoxic
& heavy-metal
free.



PES Filters are ideal for pre-filtering before PFAS analysis

BOTTLE TOP FILTRATION SYSTEM

	150 ml	250 ml	500 ml	1000 ml
--	--------	--------	--------	---------

Full Assembly Unit

0.22 µm

PES membrane

Pack size: 12

CC8221-3226

List: 102,50 €

59,45 €*

CC8221-4226

List: 114,50 €

66,41 €*

CC8221-5226

List: 144,50 €

83,81 €*

CC8221-6226

List: 231,00 €

133,98 €*

0.45 µm

PES membrane

Pack size: 12

CC8221-3456

List: 102,50 €

59,45 €*

CC8221-4456

List: 114,50 €

66,41 €*

CC8221-5456

List: 144,50 €

83,81 €*

CC8221-6456

List: 231,00 €

133,98 €*

Funnel Tops only

0.22 µm

PES membrane

Pack size: 24

CC8222-3226

List: 98,50 €

68,95 €*

CC8222-4226

List: 109,50 €

77,75 €*

CC8222-5226

List: 138,50 €

98,34 €*

CC8222-6226

List: 221,00 €

161,33 €*

0.45 µm

PES membrane

Pack size: 24

CC8222-3456

List: 98,50 €

68,95 €*

CC8222-4456

List: 109,50 €

77,75 €*

CC8222-5456

List: 138,50 €

98,34 €*

CC8222-6456

List: 221,00 €

161,33 €*

Media Bottles only

Bottle & Cap

Pack size: 24

CC8221-3000

List: 81,50 €

47,27 €*

CC8221-4000

List: 92,00 €

56,12 €*

CC8221-5000

List: 102,00 €

62,22 €*

CC8221-6000

List: 160,00 €

97,60 €*

All products are sterile (SAL 10⁻⁶) and individually wrapped.

starlab.click/filtration



***OFFER REF: Q25000794669**

must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

15 ml & 50 ml centrifuge tubes.

CERTIFIED



DESIGNED FOR THE CENTRIFUGATION OF CELLS

USP Class VI polymers and certified non-toxicity ensure the viability and integrity of your cells. The sterility assurance level (SAL) of 10^{-6} fulfills the highest demands on sterility, meeting EN 556-1 standard (Sterilisation of Medical Devices).

- White writing surface and clear, black print graduations
- Longer length screw caps with sealing ring prevents leakage



15 ML CENTRIFUGE TUBES
FROM 56,16 €*



50 ML CENTRIFUGE TUBES
FROM 73,13 €*

Description	Pack Size	Art. No.	List	OFFER*
15 ml Centrifuge Tubes				
15 ml Centrifuge Tube, Conical, Loose (Sterile)	20 x 25	E1415-0100	144,00 €	56,16 €
15 ml Centrifuge Tube, Conical, Racked (Sterile)	20 x 25	E1415-0700	165,00 €	64,35 €
50 ml Centrifuge Tubes				
50 ml Centrifuge Tube, Conical, Loose (Sterile)	20 x 25	E1450-0100	187,50 €	73,13 €
50 ml Centrifuge Tube, Conical, Racked (Sterile)	20 x 25	E1450-0700	211,00 €	82,29 €



starlab.click/centrifugetubes

Centrifuge Tube Storage.



15-place rack for 15 ml tubes

Chemically resistant epoxy-coated wire rack, blue.
127 x 76 x 176 mm

E2396-1515
(1 Rack)

List: 104,50 €
OFFER: 40,76 €*



8-place rack for 50 ml tubes

Chemically resistant epoxy-coated wire rack, blue.
147 x 80 x 89 mm

E2396-1550
(1 Rack)

List: 100,00 €
OFFER: 39,00 €*



60-place rack for 15 ml tubes

Acetal-polymer rack. White/Various colours.
293 x 115 x 65 mm

E2396-9711 Blue x1 37,50 €
17,25 €*
E2396-9712 Green x1 37,50 €
17,25 €*
E2396-9713 Orange x1 37,50 €
17,25 €*
E2396-9715 Purple x1 37,50 €
17,25 €*
E2396-9716 Yellow x1 37,50 €
17,25 €*



18-place rack for 50 ml tubes

Acetal-polymer rack. White/Various colours.
293 x 115 x 65 mm

E2396-9721 Blue x1 37,00 €
17,00 €*
E2396-9722 Green x1 37,00 €
17,00 €*
E2396-9723 Orange x1 37,00 €
17,00 €*
E2396-9725 Purple x1 37,00 €
17,00 €*
E2396-9726 Yellow x1 37,00 €
17,00 €*



***OFFER REF: Q25000794669**
must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

Robotic Tips.

Biomek® i-Series • Biomek® FX • Biomek® NX •
Biomek® 3000 • Biomek® 4000 • Multimek®



RACKED TIPS (10 x 96)

20 µl Tip, Racked	E1076-4400	97,50 €	40,95 €* *
20 µl Tip, Racked (Sterile)	E1076-4410	109,50 €	45,99 €* *
20 µl Filter Tip, Racked (Sterile)	E1076-4418	257,50 €	108,15 €* *
50 µl Tip, Racked	E1076-2400	97,50 €	40,95 €* *
50 µl Tip, Racked (Sterile)	E1076-2410	109,50 €	45,99 €* *
50 µl Filter Tip, Racked (Sterile)	E1076-2418	255,50 €	107,31 €* *
250 µl Tip, Racked	E1076-0400	97,50 €	40,95 €* *
250 µl Tip, Racked (Sterile)	E1076-0410	109,50 €	45,99 €* *
175 µl Filter Tip, Racked (Sterile)	E1076-0418	257,50 €	108,15 €* *

REFILLABLE RACKS (10 x 96)

20 µl Tip, Refillable Rack	E1076-4800	102,00 €	42,84 €* *
250 µl Tip, Refillable Rack	E1076-0800	102,00 €	42,84 €* *


REFILLS (10 x 96)

20 µl Tip, Refill	E1076-4700	97,50 €	40,95 €* *
250 µl Tip, Refill	E1076-0700	97,50 €	40,95 €* *

Non-sterile tips certified RNase, DNase and DNA free. Tested Pyrogen free. Sterile tips certified RNase, DNase, DNA & Pyrogen free. Non-refillable racks also compatible with Gilson® Cyberlab. Available whilst stock lasts

starlab.click/biomek

Packard MultiProbe II, MultiProbe IIHT/EX.
Tecan (with LiHa) Freedom EVO. Genesis.
MiniPrep.



Supplied in individual
blister packs (24 x 96)

TipOne® ROBOTIC TIPS (24 x 96)

50 µl Tip	R1181-5700	324,00 €	187,92 €* *
50 µl Filter Tip	R1128-5810	345,50 €	200,39 €* *
200 µl Tip	R1181-0700	262,50 €	152,25 €* *
200 µl Filter Tip	R1128-0810	286,00 €	165,88 €* *
1000 µl Tip	R1181-2700	289,00 €	167,62 €* *
1000 µl Filter Tip	R1128-2810	311,50 €	180,67 €* *

TipOne® CONDUCTIVE ROBOTIC TIPS (24 x 96)

50 µl Conductive Tip	R1181-5709	412,00 €	238,96 €* *
50 µl Conductive Filter Tip	R1128-5819	433,00 €	251,14 €* *
200 µl Conductive Tip	R1181-0709	259,00 €	150,22 €* *
200 µl Conductive Filter Tip	R1128-0819	321,00 €	186,18 €* *
1000 µl Conductive Tip	R1181-2709	284,50 €	165,01 €* *
1000 µl Conductive Filter Tip	R1128-2819	377,00 €	218,66 €* *

All tips are non-sterile and certified RNase, DNase, DNA & Pyrogen free. Also compatible with Qiagen® QIAgility® (CAS-1200) with single-channel head and Aurora Biomed® Systems with 8-channel head. Available only whilst stocks last.



Seated in white hanging racks for use with robots
with Varispan 4 or 8 head probes

starlab.click/tipone-tecan



***OFFER REF: Q25000794669**
must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

TipOne® 'closed-loop products' getting closer!

We are thrilled to announce that TipOne® racks made with recycled plastic from the TipOne® closed-loop recycling service will soon be available. The base of all TipOne® racks made in Germany¹⁾ will contain **75% TipOne® recyclate** plus 25 % recyclate from another source.

Providing a closed-loop recycling service is part of our commitment to sustainability and reducing environmental impact. With TipOne®, you not only receive the high-quality pipette tips you have come to expect, you know you are using a truly sustainable and plastic-saving pipette tip system.

How will I know when my TipOne® rack contains TipOne® recycled plastic?

Look out for the new design shelf boxes for racked tips made in Germany. Just like the blue shelf box for products made in USA²⁾, the latest shelf boxes are also made with 100 % recycled cardboard, but now also boast a design that uses a lot less ink, and have labels made from recycled paper. More small but important steps on the TipOne® sustainability journey.

starlab.click/tipone

¹⁾ TipOne® manufacturing is moving from USA to Germany. German production started January 2025; there will be a lengthy transition phase.

²⁾ During the transition phase you will receive TipOne® products from either manufacturing site.



TipOne® racks made in Germany that contain TipOne® recyclate will be supplied in a new-design shelf box!




TipOne® shelf boxes of racks made in USA.



T: +49 (0)40 675 99 39 0
info@starlab.de
www.starlab.de



T: +33 (0)1 60 13 71 70
info@starlab.fr
www.starlab.fr



T: +39 (0)2 7020 1040
info@starlab.it
www.starlab.it

Passionate for science.

