

Stay prepared because  
breakthroughs don't wait.  
**You are passionate  
for science.**

## SPECIAL OFFERS

Pipettes & Lab Equipment

[starlab.click/specialoffers](https://starlab.click/specialoffers)

Valid until 17.12.2025



Passionate for science.



# Consumables for your lab.

## Parafilm® sealing film.

Seals labware tightly to prevent spills, evaporation and contamination.  
Good for securing beakers, flasks, petri dishes and more.  
It comes in three sizes and works from -45 °C to 50 °C.



Description	Pack Size	Art. No.	OFFER
Parafilm 50 mm Wide x 75 m Long	1	<b>I3080-5075</b>	<b>51.50 €* </b>
Parafilm 100 mm Wide x 38 m Long	1	<b>I3080-1038</b>	<b>62.75 €* </b>
Parafilm 100 mm Wide x 75 m Long	1	<b>I3080-1075</b>	<b>105.25 €* </b>

## Spreaders.

Designed for any petri dish without damaging media or contact with dish wall or while spreading. Particularly convenient for extensive test series.



Description	Pack Size	Art. No.	OFFER
L-Shaped Spreader, PS (Sterile), Blue Individually Wrapped	500	<b>E1412-1001</b>	<b>76.80 €* </b>
L-Shaped Spreader, PS (Sterile), Blue Supplied in peel-packs of 5 spreaders	100 x 5	<b>E1412-1005</b>	<b>72 €* </b>
L-Shaped Spreader, PP (Sterile), Blue Supplied in peel-packs of 25 spreaders	20 x 25	<b>E1410-1025</b>	<b>68.40 €* </b>

## Bijou Containers.

Ultra-clear sterile bijou containers with secure screw caps, often used to store small volumes of liquid samples, such as chemical reagents, solvents, or biological specimens..



Description	Pack Size	Art. No.	OFFER
7.0 ml Polystyrene Bijou Container, without Label (Sterile)	700	<b>E1412-0710</b>	<b>79.46 €* </b>
7.0 ml Polystyrene Bijou Container, with Label (Sterile)	700	<b>E1412-0711</b>	<b>86.29 €* </b>
Dimensions (HxØ) incl. cap: 50 x 23 mm.			



# Pasteur pipettes.

A range of polyethylene (PE) pipettes, manufactured to ISO standards.

- Sterile and non-sterile options
- Low affinity surface
- Attached bulb avoids cross-contamination
- One-step biohazard disposal
- Shatter-proof low-density PE
- Can be frozen in liquid nitrogen



Sterile pasteur pipettes supplied in peel-back packs of 20 pipettes, with Lot. No. & expiry printed on pack.




E1414-1311 supplied individual peel-back packs with Lot. No & expiry printed on pack.


Description	Graduations	Length	Bulb Draw	Drop Volume	Stem Diameter	Pack Size	Art No.	OFFER
<b>Standard Pasteur Pipettes</b>								
1.0ml Graduated Pasteur Pipette	0.25 ml	160 mm	3.5 ml	15.2 µl	5.0 mm	500	<b>E1414-0100</b>	<b>10.88 €* </b>
1.0ml Graduated Pasteur Pipette (Sterile)						25 x 20	<b>E1414-0110</b>	<b>25.64 €* </b>
3.0ml Graduated Pasteur Pipette	0.5 ml	160 mm	3.5 ml	15.4 µl	7.5 mm	500	<b>E1414-0300</b>	<b>10.88 €* </b>
3.0ml Graduated Pasteur Pipette (Sterile)						25 x 20	<b>E1414-0310</b>	<b>25.64 €* </b>
5.0ml Graduated Pasteur Pipette	1.0 ml	210 mm	6.2 ml	Data not available	8.0 mm	400	<b>E1414-0500</b>	<b>68.75 €* </b>
<b>Fine Tip Pasteur Pipettes</b>								
1.0 ml Micro Fine Tip Pasteur Pipette	n/a	104 mm	1.0 ml	15.2 µl	5.9 mm	400	<b>E1414-1100</b>	<b>37.84 €* </b>
3.0 ml Fine Tip Pasteur Pipette						500	<b>E1414-1300</b>	<b>36.52 €* </b>
3.0 ml Fine Tip Pasteur Pipette (Sterile)	n/a	150 mm	3.5 ml	15.4 µl	4.4 mm	25 x 20	<b>E1414-1310</b>	<b>42.07 €* </b>
3.0 ml Fine Tip Pasteur Pipette (Sterile), Individually Wrapped						500	<b>E1414-1311</b>	<b>76.72 €* </b>

**NOTE:** Volumetric tolerance of all pipettes is ±10 %.

## FINE TIP



## STANDARD PIPETTES



**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop).  
Offer valid to 17.12.2025 / while stock lasts.


# ErgoOne® manual pipettes.

## ACCURACY AND PRECISION YOU CAN TRUST!

- Ergonomic and lightweight
- Minimal operating forces
- Volume lock
- Universal tip fit (excl. 5.0 ml)
- 3 year warranty

**ErgoOne®**

**from: 133.35 €\***



WATCH VIDEO

### Single-Channel Manual Pipettes

0.1 – 2.5 µl (red) <b>S7100-0125</b> <del>254€</del> <b>133.35 €* </b>	0.5 – 10 µl (red) <b>S7100-0510</b> <del>254€</del> <b>133.35 €* </b>	2 – 20 µl Micro (red) <b>S7100-0221</b> <del>254€</del> <b>133.35 €* </b>	2 – 20 µl (yellow) <b>S7100-0220</b> <del>254€</del> <b>133.35 €* </b>
10 – 100 µl (yellow) <b>S7100-1100</b> <del>254€</del> <b>133.35 €* </b>	20 – 200 µl (yellow) <b>S7100-2200</b> <del>254€</del> <b>133.35 €* </b>	100 – 1,000 µl (blue) <b>S7110-1000</b> <del>254€</del> <b>133.35 €* </b>	500 – 5,000 µl (violet) <b>S7150-5000</b> <del>254€</del> <b>133.35 €* </b>

### 8-Channel Manual Pipettes

0.5 – 10 µl (red) <b>S7108-0510</b> <del>732.50€</del> <b>388.23 €* </b>	10 – 100 µl (yellow) <b>S7108-1100</b> <del>732.50€</del> <b>388.23 €* </b>	30 – 300 µl (green) <b>S7108-3300</b> <del>732.50€</del> <b>388.23 €* </b>
---	--	---

### 12-Channel Manual Pipettes

0.5 – 10 µl (red) <b>S7112-0510</b> <del>873€</del> <b>400.35 €* </b>	10 – 100 µl (yellow) <b>S7112-1100</b> <del>873€</del> <b>445.23 €* </b>	30 – 300 µl (green) <b>S7112-3300</b> <del>873€</del> <b>445.23 €* </b>
--	---	--



**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop).  
Offer valid to 17.12.2025 / while stock lasts.

# ErgoOne® E electronic pipettes.

## INCREASE EFFICIENCY IN YOUR LAB

- Easy to use intuitive menu structure
- Lightweight, ergonomic design
- Adjustable finger rest for individual comfort
- Significantly reduced tip loading and ejection forces
- Universal tip fit
- 3 year warranty

**ErgoOne® E**  
**from: 398.57 €\***

### ErgoOne® E OPERATING MODES

Pipetting | Reverse Pipetting | Mixing  
Dispense | Electrophoresis<sup>1)</sup>



**NEW!**  
50 – 1250 µl  
multi-channel  
pipettes

### Single-Channel ErgoOne® E Pipettes

G9001-0010	G9001-0020	G9001-0200	G9001-1000	G9001-5000
0.5 – 10 µl	2.0 – 20 µl	10 – 200 µl	50 – 1,000 µl	250 – 5,000 µl
535€	535€	535€	535€	535€
<b>398.57 €* </b>	<b>398.57 €* </b>	<b>398.57 €* </b>	<b>398.57 €* </b>	<b>398.57 €* </b>

### 8-Channel ErgoOne® E Pipettes

G9008-0010	G9008-0020	G9008-0100	G9008-0300	G9008-1250
0.5 – 10 µl	1.0 – 20 µl	5.0 – 100 µl	15 – 300 µl	50 – 1250 µl
981€	981€	981€	981€	981€
<b>698.47 €* </b>	<b>698.47 €* </b>	<b>698.47 €* </b>	<b>698.47 €* </b>	<b>698.47 €* </b>

### 12-Channel ErgoOne® E Pipettes

G9012-0010	G9012-0020	G9012-0100	G9012-0300	G9012-1250
0.5 – 10 µl	1.0 – 20 µl	5.0 – 100 µl	15 – 300 µl	50 – 1250 µl
1171€	1171€	1171€	1171€	1171€
<b>839.61 €* </b>	<b>839.61 €* </b>	<b>839.61 €* </b>	<b>839.61 €* </b>	<b>839.61 €* </b>

<sup>1)</sup> GEL mode is not available on ErgoOne® E models: 1000 µl, 1250 µl, 5000 µl.

**MULTI-CHANNEL COMBO PACKS!** 10 µl to 300 µl multi-channel pipettes supplied with an individual stand (G9999-1004) and one rack of TipOne® pipette tips. 1250 µl models supplied with a pipette shelf holder (G9999-1009) and two racks of TipOne® pipette tips.

# Wide orifice pipette tips.

## PROTECT DELICATE SAMPLES

Protect your delicate samples from shear force and reduce flow resistance with Starlab's wide orifice tips. These high quality pipette tips made from 100% virgin polypropylene are ideal for use in a range of applications including next generation sequencing, cell culture, and plant sciences.


- Eliminate cell fragmentation caused by mechanical shearing
- Maintain the integrity of high molecular weight DNA like gDNA or rDNA
- Ideal for pipetting viscous solutions such as glycerol, Triton™ X100, blood, milk, syrup etc

## EASY-OPEN RACKS

No fiddly latch to catch gloves, simply press and lift the lid




200 µl



Orifice  
inside  
diameter  
1.5 mm

200 µl XL



Orifice  
inside  
diameter  
1.0 mm

1000 µl XL



Orifice  
inside  
diameter  
1.63 mm

Our XL tip range.

Description	Pack Type	Sterile	Pack Size	Art. No.	OFFER
200µl Graduated Tip / Wide Orifice	96-Tip Rack	No	10 x 96	E1011-8400	40.98 €*
	96-Tip Rack	Yes	10 x 96	E1011-8410	45.38 €*
200µl Graduated Filter Tip / Wide Orifice	96-Tip Rack	Yes	10 x 96	E1011-8420	67.92 €*
200µl XL Graduated Tip / Wide Orifice	96-Tip Rack	No	2 x 96	E1011-7400	22.20 €*
	96-Tip Rack	Yes	2 x 96	E1011-7410	24.90 €*
200µl XL Graduated Filter Tip / Wide Orifice	96-Tip Rack	Yes	2 x 96	E1011-7420	31.50 €*
	Bag	No	1000	E1011-9000	23.38 €*
1000µl XL Graduated Tip / Wide Orifice	96-Tip Rack	No	10 x 96	E1011-9400	42.35 €*
	96-Tip Rack	Yes	10 x 96	E1011-9410	46.48 €*
1000µl XL Graduated Filter Tip / Wide Orifice	96-Tip Rack	Yes	10 x 96	E1011-9420	69.03 €*

**Tip length.** 200µl: 50 mm long. 200µl XL: 91.1 mm long. 1000µl XL: 86.4 mm long. NOTE: Filter tips are not autoclavable.



**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop).  
Offer valid to 17.12.2025 / while stock lasts.



# Repeat dispensing.

## UP TO 48 DISPENSES FROM THE SAME SAMPLE!

Use the RepeatOne® S to dispense precise volumes of liquid in a series, without the need to aspirate in between each step.

### RepeatOne S®

Art. No. G4026-1500

**OFFER: 254.10 €\***

Dispenser tips not included.



## TipOne® Repeat dispenser tips.

Tip volume	Pack Size	Article No.	OFFER*	Tip Volume	Pack Size	Article No.	OFFER*
Non-Sterile				Sterile, Individually Wrapped			
0.1 ml	100	S4751-0010	55.55 €* <sup>1</sup>	0.1 ml	100	S4761-0010	92.95 €* <sup>1</sup>
0.5 ml	100	S4751-0050	40.15 €* <sup>1</sup>	0.5 ml	100	S4761-0050	82.22 €* <sup>1</sup>
1.0 ml	100	S4751-0100	40.15 €* <sup>1</sup>	1.0 ml	100	S4761-0100	82.22 €* <sup>1</sup>
2.5 ml	100	S4751-0250	40.15 €* <sup>1</sup>	2.5 ml	100	S4761-0250	82.22 €* <sup>1</sup>
5.0 ml	100	S4751-0500	40.15 €* <sup>1</sup>	5.0 ml	100	S4761-0500	82.22 €* <sup>1</sup>
10 ml	100	S4751-1000	40.15 €* <sup>1</sup>	10 ml	100	S4761-1000	82.22 €* <sup>1</sup>
25 ml <sup>1)</sup>	100	S4751-2500	69.30 €* <sup>1</sup>	25 ml <sup>1)</sup>	100	S4761-2500	118.80 €* <sup>1</sup>
50 ml <sup>1)</sup>	100	S4751-5000	82.22 €* <sup>1</sup>	50 ml <sup>1)</sup>	100	S4761-5000	122.10 €* <sup>1</sup>
				25 ml Adapter	7	S4765-2500	39 €* <sup>1</sup>
				50 ml Adapter	7	S4765-5000	36.90 €* <sup>1</sup>

<sup>1)</sup> 25 ml and 50 ml dispenser tips are supplied with four adapters.  
Note: TipOne® Repeat Dispenser Tips are not autoclavable.






# Aspirator.

- Chemically resistant, long-life high-performance pump
- Auto ON/OFF switch with low-noise operation
- Ergonomic, autoclavable hand controller
- The powerful pump can operate two hand controllers without any loss of performance

**OFFER: 1448.45 €\*<sup>1</sup>**

See how to add a second hand controller.



Description	Art. No.	List	OFFER*
 <b>Aspirator</b> Incl. everything you need to get started! 1-year warranty.	N2400-9000	£2230.40	<b>1448.45 €*<sup>1</sup></b>
 <b>Hand Controller for Aspirator</b> Incl. 1x Adapter for Pipettes, 1x Adapter for Pipette Tips, 1x Wall Mount, and 2.5 m Tubing.	N2400-9001	£202.95	<b>133.11 €*<sup>1</sup></b>
 <b>Connection Kit for Second Hand Controller</b>	N2400-9006	£50.21	<b>33.06 €*<sup>1</sup></b>




**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop).  
Offer valid to 17.12.2025 / while stock lasts.

# Essential lab equipment.

## Orbital Shaker.



Speed range: 100 - 500 rpm.  
LCD displays for speed and time.  
Continuous or timed mixing.  
Supplied without platform.


**Orbital Shaker**

**Art. No: N2400-8030**  
(excludes platform)

List: 1347.50€

**OFFER: 761.34 €\***

## 3D / Rocking Shakers.



See saw rocking or 3D mixing.  
Speed range: 10 - 70 rpm.  
LCD displays for speed and time.  
Continuous or timed mixing.  
Supplied with platform with bungee cords.

**Rocking Shaker**

**Art. No: N2400-8010**

List: 1532.50€

**OFFER: 865.86 €\***

**3D Rocking Shaker**

**Art. No: N2400-8020**

List: 1532.50€


**OFFER: 861.26 €\***



Platform with clamping rolls  
(for use with Orbital, 3D and Rocking Shakers)

**Art. No: N2400-8011**

**OFFER: 211.70 €\***



Platform with rubber mat and bungee cords  
(for use with Orbital, 3D and Rocking Shakers)

**Art. No: N2400-8021**

**OFFER: 211.70 €\***



## Rotator Mixer.

**Art. No. N2400-5000**

List: 1297.50€

**OFFER: 583.88 €\***

Controlled sample inverting.  
Speed range: 1 – 100 rpm.  
Programmable.  
Three different motions.  
Variable speed and timer.  
Low profile and small footprint.



### Vertical Rotation

360° vertical motion for time period of 0–250 seconds, or non-stop.



### Reciprocal Motion

Platform direction changes in turns (turning angle 1–90°) for 0–250 seconds, or non-stop.



### Vibration

Same as reciprocal motion, plus additional vibration (turning angle 1–5° for 1–5 seconds.



### Pause


Platform direction changes in turns, incl. a temporary motion stop of 1–5 seconds (in conjunction with reciprocal motion only).



**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop).  
Offer valid to 17.12.2025 / while stock lasts.





Use for mixing and shaking any skirted plate with a SBS footprint, and for resuspending pellets. Speed range: 0 – 1500 rpm. Stepless speed control. Range of accessories available.

### Microplate Mixer.

Art. No. N2400-8040

List: 533.50€

**OFFER: 304.10 €\***



An economical magnetic stirrer. Speed range: 100 – 1500 rpm. Supplied with a 40 mm stir bar.

### Magnetic Stirrer.

Art. No. N2400-2010

List: 237€

**OFFER: 133.91 €\***



Fast temperature ramping saves time! Speed range: 200 – 1500 rpm. Temp. range: ambient to 280°C (1°C increments). Supplied with a 40 mm stir bar.

### Magnetic Stirrer with Heating.

Art. No. N2400-3010

List: 459€

**OFFER: 259.33 €\***

# The Starlab Minicentrifuge.

WAY AHEAD OF THE CURVE



WATCH VIDEO

Load, close the lid, and spin! One-touch STOP brakes quickly and opens the lid for easy, one-handed use. Compact, quiet, and smooth— with a maximum speed of 6,500 rpm.

### Supplied with two rotors:

Tube rotor for 8 x microcentrifuge tubes.

PCR strip rotor for 32 single PCR tubes or 4 PCR strips.

### Minicentrifuge

Art. No. S8030-5000

List: 278.50€

**OFFER: 194.11 €\***



**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop). Offer valid to 17.12.2025 / while stock lasts.

# Starlab Vortex Mixer.

## THE TOUCH OF A VORTEXER.

## THE POWER OF A MIXER

Ideal for the gentle mixing of samples right up to vigorous vortex mixing and the resuspending of cells or chemical pellets.

- Precise vortexing speeds
- Speed range: 300 to 3,500 rpm
- Unique, glove compatible touch screen
- One-touch pre-set options for speed / time
- Touch or continuous operation
- End of run audible alarm
- Power-saving sleep mode

## Vortex Mixer.

Art. No. S8010-0000

List: 409€

**OFFER: 285.07 €\***



### TOUCH MODE

Minimal pressure required to maintain vibration. Mix single or multiple tubes up to 3,500 rpm.

GLOVE COMPATIBLE TOUCHSCREEN for precise speed and time settings

SPEED RANGE  
300 – 3,500 RPM



WATCH VIDEO

## ACCESSORIES



### Adapter for Tubes

For mixing up to 15 x 1.5/2.0 ml tubes up to 1,500 rpm

Art. No. S8010-0011

List: 28.50€

**OFFER: 16.53 €\***



### Adapter for Microplate

For mixing one microplate up to 1,500 rpm

Art. No. S8010-0012

List: 17.50€

**OFFER: 10.50 €\***

Tubes and plate are not included with the adapters.

**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop).  
Offer valid to 17.12.2025 / while stock lasts.

# The Plate Centrifuge.

**NEW!**

## GET SPINNING. IMPROVE PCR YIELDS.

Starlab's innovative Plate Centrifuge is tailored for the rapid and effortless centrifugation of samples in PCR plates. Use before and after thermal cycling to enhance your yield.

- Holds two plates
- Compatible with non-skirted, semi-skirted and skirted 96- and 384-well plates
- 2,500 rpm
- Two operating modes
- Starts / stops in seconds
- Convenient small footprint
- Optional adapter available for half plates or tubes
- 2 year guarantee

## Plate Centrifuge

Art. No. N2631-0100

List: ~~1136€~~

**OFFER: 795.20 €\***

### SHORT mode

Press and hold to run. When maximum speed is reached, the timer will start counting up in seconds. Release button to stop

### ACCEPTS TWO PLATES

any non-skirted, semi-skirted or skirted plate, up to 25 mm high

Optional adapter available for half plates or 0.2 ml single or strip tubes

### INSTANT ACCESS

for quick loading and unloading of plates

OPEN LID

### TIMER SETTING

Set the timer using the arrows. Timer can be set to seconds or minutes

Press START button and the timer will count down. When the set time has elapsed, a beep will sound once and the unit will automatically stop spinning

SMALL FOOTPRINT



**WATCH VIDEO**



**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop).  
Offer valid to 17.12.2025 / while stock lasts.

# GuardOne® Workstation.



Contents not included.

## AVOID CROSS-CONTAMINATION DURING AMPLIFICATION OF DNA/RNA

Protect your valuable samples and reagents<sup>1)</sup> to save time and costs by avoiding the need to repeat contaminated work.

Starlab's GuardOne® Workstation provides

- An ideal clean and controlled environment that can help you achieve accurate PCR results
- Convenience! All functions are controlled via the glove compatible, easy to use touchscreen
- Quick, tool-free HEPA filter exchange and UV lamp exchange

Perfect for your PCR / qPCR applications or any other application that requires a controlled, clean environment.

- One-touch QuickStart / QuickStop settings
- Efficient UV<sub>245 nm</sub> decontamination with timer function
- Efficient, class H14 HEPA filter with pre-filter provides ISO class 5 working environment<sup>2)</sup>
- Real-time airflow measurement<sup>2)</sup>
- Life-span monitoring of consumables
- Virtually maintenance-free
- Easy to clean
- Easy, tool-free 15-minute assembly

<sup>1)</sup>GuardOne® is not designed to protect the user.

<sup>2)</sup>Laminar flow models only.



**DEAD AIR 32"**

**S8040-1000**

List: 4398.50€

**OFFER: 2652.30 €\***



**DEAD AIR 48"**

**S8040-1001**

List: 5579.50€

**OFFER: 3425.81 €\***



**LAMINAR FLOW 32"**

**S8040-0000**

List: 5226.50€

**OFFER: 3151.58 €\***



**LAMINAR FLOW 48"**

**S8040-0001**

List: 6966€

**OFFER: 4200.50 €\***



**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop).

Offer valid to 17.12.2025 / while stock lasts.

# Freezer racking.

## SUSTAINABILITY JUST GOT COOLER

- Freezer racking helps optimise capacity by keeping freezers tidy and maximising storage space for efficient energy use
- High quality racking prevents the need for regular repair or exchange
- Made from 85% recycled stainless steel
- Manufactured in a factory powered by 100% sustainable energy
- Made in Germany. Short journey from factory to customer reduces carbon footprint

## Cardboard storage boxes.

- Moisture resistant with removable grid insert
- Ideal for use in the freezer racking on this page



Description	Art. No	OFFER*
<b>49-place</b> cardboard storage box with removable 7x7 dividers for tubes up to 16 mm diameter, eg. 1.2 / 2.0 ml cryogenic vials	<b>E2023-4949</b>	<b>3.65 €*</b>
<b>81-place</b> cardboard storage box with removable 9x9 dividers for tubes up to 13 mm diameter, eg. 1.5 / 2.0 ml microcentrifuge tubes, screw cap microtubes, 1.2 / 2.0 ml cryogenic vials with internal thread	<b>E2023-8181</b>	<b>3.65 €*</b>
<b>100-place</b> cardboard storage box with removable 10x10 dividers for tubes up to 12 mm diameter, eg. screw cap microtubes	<b>E2023-3100</b>	<b>3.65 €*</b>

Dimensions (WxDxH): 133 x 133 x 53 mm (with lid). Max. tube diam. stated includes tube's cap. Tubes not included.

[www.starlab.click/cardboardstorage](https://www.starlab.click/cardboardstorage)

## RACKS WITH SIDE ACCESS



Boxes not included.

## RACKS WITH PULL OUT SHELVES



Boxes not included.

## RACKS WITH PULL OUT SHELVES



Boxes not included.

## Freezer racking for UPRIGHT FREEZERS for use with all **standard storage boxes** up to 136 x 136 x 53 mm.

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No	OFFER*
9-Place	3 deep x 3 high	139 x 426 x 169 mm	<b>E6000-0033</b>	<b>105.30 €*</b>
12-Place	3 deep x 4 high	139 x 426 x 224 mm	<b>E6000-0034</b>	<b>120.58 €*</b>
15-Place	3 deep x 5 high	139 x 426 x 297 mm	<b>E6000-0035</b>	<b>134.55 €*</b>
12-Place	4 deep x 3 high	139 x 565 x 169 mm	<b>E6000-0043</b>	<b>120.58 €*</b>
16-Place	4 deep x 4 high	139 x 565 x 224 mm	<b>E6000-0044</b>	<b>120.58 €*</b>
20-Place	4 deep x 5 high	139 x 565 x 279 mm	<b>E6000-0045</b>	<b>137.80 €*</b>
24-Place	4 deep x 6 high	139 x 565 x 334 mm	<b>E6000-0046</b>	<b>157.95 €*</b>

Racking suitable for use with storage boxes up to maximum (WxDxH) 136 x 136 x 53 mm.

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No	OFFER*
9-Place	3 deep x 3 high	140 x 450 x 173 mm	<b>E6001-0033</b>	<b>145.93 €*</b>
12-Place	3 deep x 4 high	140 x 450 x 231 mm	<b>E6001-0034</b>	<b>151.12 €*</b>
15-Place	3 deep x 5 high	140 x 450 x 288 mm	<b>E6001-0035</b>	<b>156.65 €*</b>
12-Place	4 deep x 3 high	140 x 569 x 174 mm	<b>E6001-0043</b>	<b>156.65 €*</b>
16-Place	4 deep x 4 high	140 x 569 x 231 mm	<b>E6001-0044</b>	<b>161.53 €*</b>
20-Place	4 deep x 5 high	140 x 569 x 281 mm	<b>E6001-0045</b>	<b>167.05 €*</b>
24-Place	4 deep x 6 high	140 x 569 x 337 mm	<b>E6001-0046</b>	<b>179.40 €*</b>
28-Place	4 deep x 7 high	140 x 569 x 403 mm	<b>E6001-0047</b>	<b>191.10 €*</b>

Racking suitable for use with storage boxes up to maximum (WxDxH) 136 x 136 x 53 mm.

## Racking for UPRIGHT FREEZERS for use with Starlab's **StarStore 100** storage boxes only.

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No	OFFER*
9-Place	3 deep x 3 high	147 x 500 x 179 mm	<b>E6001-0133</b>	<b>163.80 €*</b>
12-Place	3 deep x 4 high	147 x 500 x 239 mm	<b>E6001-0134</b>	<b>169.97 €*</b>
16-Place	4 deep x 4 high	147 x 655 x 239 mm	<b>E6001-0144</b>	<b>175.50 €*</b>

Racking suitable for use with StarStore 100 storage boxes (WxDxH) 140 x 155 x 55 mm.



**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop).  
Offer valid to 17.12.2025 / while stock lasts.



# Clearance deals!

ONLY AVAILABLE WHILST STOCK LASTS!

## Microflex Supreno® Blue Nitrile Gloves



Glove Size	Pack Size	Article No.	OFFER*
Large	10 x 100	SU-INT-L	75.15 €*
X-Large	10 x 100	SU-INT-XL	75.15 €*

PPE Complex Design Category III.

## Ansell TouchNTuff® Green Nitrile Gloves



Glove Size	Pack Size	Article No.	OFFER*
Medium	10 x 100	92-600-M	131.42 €*
Large	10 x 100	92-600-L	131.42 €*

PPE Complex Design Category III.

## StarGuard® Sensitive Blue Nitrile Gloves



Glove Size	Pack Size	Article No.	OFFER*
X-Large	10 x 200	SG-N-XL	83.80 €*

PPE Complex Design Category III.  
Low Dermatitis Potential: ideal for users with sensitive skin.  
Contains no contaminants to interfere with experiments.

## CytoOne® Flasks with 2-Position Plug Seal Caps - Non Treated



Description	Pack Size	Article No.	OFFER*
T-75	20 x 5	CC7672-4175	62.72 €*
T-225	5 x 5	CC7672-4225	38.23 €*

Sterile, and certified RNase, DNase, DNA and Pyrogen free.

## Breathable Self-Adhesive Film (Sterile)



Colour	Pack Size	Article No.	OFFER*
White	50	E2796-3015	51.08 €*

Sterile, and certified RNase, DNase and Pyrogen free.



**\*OFFER REF: Q25000853348**

must be stated on your PO (excl. eShop).  
Offer valid to 17.12.2025 / while stock lasts.

Please note that delivery may  
take a little longer than usual  
for some clearance products.

# Reducing environmental impact with the TipOne® closed-loop.



## The circle is complete!

We are thrilled to announce that the first TipOne® racks made using post-consumer recycled material from TipOne® products are now on sale.

This marks the next milestone for TipOne® and our closed-loop recycling service: turning lab waste back into high-quality products.



TipOne® racks<sup>1)</sup> produced in Germany, are made entirely from recycled plastic and include 75 % TipOne® recylcate in the base.

[starlab.click/tiponefiltertips](https://starlab.click/tiponefiltertips)

Together, with over 2000 customers using the service, the closed-loop saved 38.72 tonnes of CO<sub>2</sub> equivalent during the year 2024. We can all be proud of the positive impact our efforts are making in reducing both waste and CO<sub>2</sub> emissions.

For more information on the TipOne® closed-loop recycling service and to join, go to [www.starlab.click/recycling](https://www.starlab.click/recycling)



**WIN!**  
**A TipOne®**  
**Lego® set.**  
See back page.



The TipOne® closed-loop recycling service.

## RECYCLING SERVICE IMPACT STATISTICS



**~117**  
**tonnes PP<sup>2)</sup>**

By end of Jan-25, ~ 117 tonnes of raw material had been collected since the service began. This means our need for virgin material has been reduced by 117 tonnes, resulting in a direct saving of natural resources



**38.72**  
**tonnes CO<sub>2</sub>**

The amount of CO<sub>2</sub> equivalent saved by the TipOne® recycling program in the year 2024.<sup>3)</sup>



**~1575**  
**trees**

38.72 tonnes of CO<sub>2</sub> equivalent is the annual carbon sequestration of ~1575 trees.<sup>4)</sup>



**~40%**  
**Hyde Park**

1575 trees equates to almost 40 % of the number of trees (4000) currently growing in Hyde Park, London.<sup>5)</sup>



Starlab will repeatedly reuse the recylcate in its own TipOne® products, so the actual savings will increase with every cycle through the closed loop.

<sup>1)</sup> Rack base and lid.

<sup>2)</sup> Polypropylene.

<sup>3)</sup> Calculated by REMONDIS using calculation model developed by the Fraunhofer Institute for Environmental, Safety & Energy Technology (UMSICHT).

<sup>4)</sup> Calculated using average data per tree from [www.fortomorrow.eu/de/blog/co2-baum](https://www.fortomorrow.eu/de/blog/co2-baum)

<sup>5)</sup> [www.visitlondon.com/de/sehen-und-erleben/ort/619718-hyde-park](https://www.visitlondon.com/de/sehen-und-erleben/ort/619718-hyde-park)

# WIN a limited edition TipOne® Lego® set.


**TipOne® racks made using post-consumer recycled material from TipOne® products are now on sale!**

To celebrate this closed-loop milestone we are giving away 10 limited edition TipOne® Rack Lego® sets.

To get your hands on one of these fantastic prizes, follow the link and answer three questions to enter the prize draw.

Answers can be found on our website at:  
[starlab.click/recycling](https://starlab.click/recycling)

Prize draw closes 17-Dec-2025.



**ENTER HERE**

[starlab.click/win-lego](https://starlab.click/win-lego)



T: +49 (0)40 675 99 39 0  
[info@starlab.de](mailto:info@starlab.de)  
[www.starlab.de](https://www.starlab.de)



T: +44 (0)1908 283800  
[info@starlab.co.uk](mailto:info@starlab.co.uk)  
[www.starlab.co.uk](https://www.starlab.co.uk)



T: +33 (0)1 60 13 71 70  
[info@starlab.fr](mailto:info@starlab.fr)  
[www.starlab.fr](https://www.starlab.fr)



T: +39 (0)2 7020 1040  
[info@starlab.it](mailto:info@starlab.it)  
[www.starlab.it](https://www.starlab.it)



T: +41 56 6645 980  
[info@starlab.ch](mailto:info@starlab.ch)  
[www.starlab.ch](https://www.starlab.ch)

Passionate for science.

