

Starlab Special Offers



Screw Cap Tubes with caps

Art. No.	Description	Pack Size	List Price	Your Price
E1405-0301	0.5 ml Plain Skirted Tube with EasyGrip Screw Cap White/Blue (Sterile)	5 x 50	€85.17	€43.44
E1405-0302	0.5 ml Plain Skirted Tube with EasyGrip Screw Cap White/Green (Sterile)	5 x 50	€85.17	€43.44
E1405-0304	0.5 ml Plain Skirted Tube with EasyGrip Screw Cap White/Red (Sterile)	5 x 50	€85.17	€43.44
E1405-0306	0.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Yellow (Sterile)	5 x 50	€85.17	€43.44
E1405-2141	0.5 ml Plain Skirted Tube, Natural Standard Screw Cap (Sterile)	5 x 50	€78.54	€40.06
E1415-0301	1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Blue (Sterile)	5 x 50	€85.17	€43.44
E1415-0302	1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Green (Sterile)	5 x 50	€85.17	€43.44
E1415-0304	1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Red (Sterile)	5 x 50	€85.17	€43.44
E1415-0306	1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Yellow (Sterile)	5 x 50	€85.17	€43.44
E1415-0311	1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Blue (Sterile)	5 x 50	€85.17	€43.44
E1415-0312	1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Green (Sterile)	5 x 50	€85.17	€43.44
E1415-0314	1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Red (Sterile)	5 x 50	€85.17	€43.44
E1415-0316	1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Yellow (Sterile)	5 x 50	€85.17	€43.44
E1415-2231	1.5 ml Graduated Conical Tube, Natural Standard Screw Cap (Sterile)	5 x 50	€78.54	€40.06
E1415-2241	1.5 ml Plain Skirted Tube, Natural Standard Screw Cap (Sterile)	5 x 50	€78.54	€40.06
E1415-2243	1.5 ml Plain, Skirted Tube, Natural EasyGrip Cap (Sterile)	5 x 50	€85.17	€43.44
E1420-0301	2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Blue (Sterile)	5 x 50	€116.28	€59.30
E1420-0302	2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Green (Sterile)	5 x 50	€116.28	€59.30
E1420-0304	2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Red (Sterile)	5 x 50	€116.28	€59.30
E1420-0306	2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Yellow (Sterile)	5 x 50	€89.25	€45.52
E1420-0311	2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Blue (Sterile)	5 x 50	€85.17	€43.44
E1420-0312	2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Green (Sterile)	5 x 50	€85.17	€43.44
E1420-0314	2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Red (Sterile)	5 x 50	€85.17	€43.44
E1420-0316	2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Yellow (Sterile)	5 x 50	€85.17	€43.44
E1420-2331	2.0 ml Graduated Conical Tube, Natural Standard Screw Cap (Sterile)	5 x 50	€77.01	€39.28
E1420-2341	2.0 ml Graduated Skirted Tube, Natural Standard Screw Cap (Sterile)	5 x 50	€70.38	€35.89
E1420-2343	2.0 ml Graduated Skirted Tube, Natural EasyGrip Cap (Sterile)	5 x 50	€77.01	€39.28

Please use Q26000928528 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Starlab Special Offers



Screw Cap Tubes only

Art. No.	Description	Pack Size	List Price	Your Price
E1405-2130	0.5 ml Screw Cap Tube, Plain, Conical, Natural	500	€44.27	€22.63
E1405-2140	0.5 ml Screw Cap Tube, Plain, Skirted, Natural	500	€40.80	€20.81
E1405-2147	0.5 ml Screw Cap Tube, Plain, Skirted, Amber	500	€45.90	€23.41
E1415-2230	1.5 ml Screw Cap Tube, Graduated, Conical, Natural	500	€35.19	€17.95
E1415-2240	1.5 ml Screw Cap Tube, Plain, Skirted, Natural	500	€34.68	€17.69
E1420-2330	2.0 ml Screw Cap Tube, Graduated, Conical, Natural	500	€35.19	€17.95
E1420-2340	2.0 ml Screw Cap Tube, Graduated, Skirted, Natural	500	€33.66	€17.17



Ultra High Recovery Microcentrifuge Tubes

Art. No.	Description	Pack Size	List Price	Your Price
E1405-2600	0.5 ml Ultra High Recovery Microcentrifuge Tube	500	€39.27	€20.03
E1415-2600	1.5 ml Ultra High Recovery Microcentrifuge Tube	250	€22.95	€11.70
I1420-2600	2.0 ml Ultra High Recovery Microcentrifuge Tube	250	€26.52	€14.59



Centrifuge Tubes

Art. No.	Description	Pack Size	List Price	Your Price
E1415-0100	15 ml Centrifuge Tube, Conical, Loose (Sterile)	500	€146.08	€64.63
E1415-0700	15 ml Centrifuge Tube, Conical, Racked (Sterile)	500	€169.30	€79.10
E1450-0100	50 ml Centrifuge Tube, Conical, Loose (Sterile)	500	€101.25	€84.15
E1450-0400	50 ml Centrifuge Tube, Skirted, Loose (Sterile)	500	€212.67	€99.95
E1450-0700	50 ml Centrifuge Tube, Conical, Racked (Sterile)	500	€215.22	€101.15

Please use Q26000928528 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Starlab Special Offers



Screw Caps only


Art. No.	Description	Pack Size	List Price	Your Price
E1480-0100	Standard Caps for Screw Cap Microtubes, natural	500	€63.24	€32.25
E1480-0101	Standard Cap for Screw Cap Microtubes, Blue	500	€68.34	€34.85
E1480-0102	Standard Cap for Screw Cap Microtubes, Green	500	€68.34	€34.85
E1480-0103	Standard Cap for Screw Cap Microtubes, Orange	500	€68.34	€34.85
E1480-0104	Standard Caps for Screw Cap Microtubes, red	500	€68.34	€34.85
E1480-0105	Standard Cap for Screw Cap Microtubes, Violet	500	€68.34	€34.85
E1480-0106	Standard Cap for Screw Cap Microtubes, Yellow	500	€68.34	€34.85
E1480-0107	Standard Cap for Screw Cap Microtubes, Amber	500	€68.34	€34.85
E1480-0108	Standard Cap for Screw Cap Microtubes, Opaque Black	500	€81.60	€41.62
E1480-0109	Standard Cap for Screw Cap Microtubes, Opaque White	500	€68.34	€34.85
E1480-0300	EasyGrip Cap, Natural / White	500	€68.34	€34.85
E1480-0301	Easy Grip Cap, White / Blue	500	€71.40	€36.41
E1480-0302	Easy Grip Cap, White / Green	500	€71.40	€36.41
E1480-0303	Easy Grip Cap, White / Orange	500	€71.40	€36.41
E1480-0304	Easy Grip Cap, White / Red	500	€71.40	€36.41
E1480-0305	Easy Grip Cap, White / Violet	500	€71.40	€36.41
E1480-0306	Easy Grip Cap, White / Yellow	500	€71.40	€36.41
E1480-0309	Easy Grip Cap, White / White	500	€71.40	€36.41
E1480-0399	Easy Grip Caps, Mixed	500	€66.30	€33.81

Cryogenic Vials - Cap Inserts

Art. No.	Description	Pack Size	List Price	Your Price
E1494-6001	Coloured Cap Inserts for Cryogenic Vials, Blue	500	€39.78	€19.89
E1494-6002	Coloured Cap Inserts for Cryogenic Vials, Green	500	€39.78	€19.89
E1494-6003	Coloured Cap Inserts for Cryogenic Vials, Orange	500	€39.78	€19.89
E1494-6004	Coloured Cap Inserts for Cryogenic Vials, Red	500	€39.78	€19.89
E1494-6005	Coloured Cap Inserts for Cryogenic Vials, Violet	500	€39.78	€19.89
E1494-6006	Coloured Cap Inserts for Cryogenic Vials, Yellow	500	€39.78	€19.89
E1494-6007	Coloured Cap Inserts for Cryogenic Vials, Brown	500	€39.78	€19.89
E1494-6009	Coloured Cap Inserts for Cryogenic Vials, White	500	€39.78	€19.89
E1494-6099	Coloured Cap Inserts for Cryogenic Vials, Mixed	500	€39.78	€19.89
E1494-6104	Coloured Cap Inserts for Cryogenic Vials, Pink	500	€39.78	€19.89
E1494-6108	Coloured Cap Inserts for Cryogenic Vials, Grey	500	€39.78	€19.89


Please use Q26000928528 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Starlab Special Offers



Cardboard Storage Boxes

Art. No.	Description	Pack Size	List Price	Your Price
E2023-1615	Cardboard Storage Box for 15 ml Centrifuge Tubes, 36 places	1	€12.24	€5.50
E2023-1650	Cardboard Storage Box for 50 ml Centrifuge Tubes, 16 places	1	€12.00	€5.50
E2023-3100	Cardboard Storage Box, 133x133 Footprint, White, 100-Place Grid	1	€6.10	€3.67
E2023-8181	Cardboard Storage Box, 133x133 Footprint, White, 81-Place Grid	1	€5.10	€3.62



Cryogenic Storage Boxes

Art. No.	Description	Pack Size	List Price	Your Price
E2350-4001	50 Place Lock-Rack, Blue	4	€100.47	€51.24
E2350-4004	50 Place Lock-Rack, Red	4	€117.30	€59.82
E2350-4005	50 Place Lock-Rack, Violet	4	€117.30	€59.82
I2510-5005	100-place Polycarbonate Storage Box	5	€49.47	€26.71



Tube Racks

Art. No.	Description	Pack Size	List Price	Your Price
E2380-5000	StarRack 80 polypropylene fraction collection rack, Natural	5	€27.54	€13.77
E2380-5048	StarRack 80 polypropylene fraction collection rack, Neon Mixed	5	€34.17	€17.09
E2396-0010	StarRack 96 double-sided polypropylene rack, Natural	Each	€22.50	€10.13
E2396-1515	Racks for 15ml centrifuge tubes (15-place)	1	€106.59	€45.83
E2396-1550	Racks for 50 ml centrifuge tubes (8-place)	1	€102.00	€43.86
I2396-5048	StarRack 96 double-sided polypropylene rack, Neon Mixed	5	€93.84	€51.61



Cryogenic Vials

Art. No.	Description	Pack Size	List Price	Your Price
E3090-6212	1.0 ml Cryovial Vials with External Thread, Self-Sealing Cap, Skirted (Sterile)	10 x 50	€276.93	€121.85
E3090-6222	1.8 ml Cryovial Vials with External Thread, Self Sealing Seal Cap, Skirted (Sterile)	10 x 50	€258.57	€111.19
E3090-6242	4.5 ml Cryovial Vials with External Thread, Self-Sealing Cap, Skirted (Sterile)	6 x 50	€206.04	€90.66
E3090-6321	1.8 ml Cryovial Vials with External Thread, Self-Sealing Cap Conical (Sterile)	10 x 50	€298.35	€131.27
E3110-6112	1.0 ml Cryovial with Internal Thread, Silicone Seal Cap, Skirted (Sterile)	10 x 50	€257.55	€144.23
E3110-6121	1.8 ml Cryovial with Internal Thread, Silicone Seal Cap Conical (Sterile)	10 x 50	€288.66	€161.65
E3110-6122	1.8 ml Cryovial with Internal Thread, Silicone Seal Cap, Skirted (Sterile)	10 x 50	€241.74	€132.96
E3110-6131	3.6 ml Cryovial with Internal Thread, Silicone Seal Cap, Conical (Sterile)	8 x 50	€211.65	€118.52
E3110-6132	3.6 ml Cryovial with Internal Thread, Silicone Seal Cap, Skirted (Sterile)	8 x 50	€239.19	€133.95
E3110-6141	4.5 ml Cryovial with Internal Thread, Silicone Seal Cap, Conical (Sterile)	6 x 50	€220.83	€123.66
E3110-6142	4.5 ml Cryovial with Internal Thread, Silicone Seal Cap, Skirted (Sterile)	6 x 50	€220.83	€123.66



Please use Q26000928528 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Starlab Special Offers




Polypropylene Storage Boxes

Art. No.	Description	Pack Size	List Price	Your Price
E9025-1650	Centrifuge Tube Storage / Transport Box for 50ml Tubes, Blue	set of 2	€47.94	€27.81
E9025-3615	Centrifuge Tube Storage / Transport Box for 15ml Tubes, Green	set of 2	€48.96	€28.40
I2310-5800	StarStore 100-Place Storage Box, Natural	5	€52.02	€28.61
I2310-5848	StarStore 100-Place Storage Box, Neon Mixed*	5	€48.45	€26.65
I2381-5040	StarStore 81, Natural	5	€46.92	€25.81
I2381-5041	StarStore 81, Neon Blue	5	€48.45	€26.65
I2381-5042	StarStore 81, Neon Green	5	€48.45	€26.65
I2381-5043	StarStore 81, Neon Orange	5	€48.45	€26.65
I2381-5044	StarStore 81, Neon Pink	5	€48.45	€26.65
I2381-5046	StarStore 81, Neon Yellow	5	€48.45	€26.65
I2381-5048	StarStore 81, Neon Mixed	5	€53.04	€29.17



Startags

Art. No.	Description	Pack Size	List Price	Your Price
E9185-0001	StarTags Laboratory Labels (9.5 mm Diameter spot), blue	2000	€35.19	€22.52
E9187-1009	StarTags Laboratory Labels, Large (32.5 mm x 12.5 mm)	1000	€35.19	€22.52
E9188-1009	StarTags Laboratory Labels, Medium (24 mm x 12.5 mm)	1000	€35.19	€22.52
E9188-8699	StarTag Laboratory Labels, Small (20 mm x 6 mm)	1000	€35.19	€22.52



Cryogenic Storage

Art. No.	Description	Pack Size	List Price	Your Price
I2381-2305	81-place Polycarbonate Storage Box	5	€49.47	€27.21



TubeOne®

Art. No.	Description	Pack Size	List Price	Your Price
S1605-0000	0.5 ml TubeOne® Microcentrifuge Tubes, Natural (non-sterile)	5 x 200	€43.35	€26.88
S1615-5500	1.5 ml TubeOne® Microcentrifuge Tubes, Natural (non-sterile)	500	€14.79	€9.17
S1615-5510	1.5 ml TubeOne® Microcentrifuge Tubes, Natural (Sterile)	5 x 100	€37.23	€23.45
S1615-5550	1.5 ml TubeOne® Microcentrifuge Tubes, Natural (non-sterile)	5 x 100	€17.85	€11.07
S1620-2700	2.0ml TubeOne® Microcentrifuge Tubes, Natural	500	€62.02	€32.25

Please use Q26000928528 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Starlab Special Offers



Manual Pipettes

Art. No.	Description	Pack Size	List Price	Your Price
S7100-0125	ErgoOne® Single Channel Pipette, 0.1 - 2.5µl	1	€259.08	€145.00
S7100-0220	ErgoOne® Single Channel Pipette (Standard Cone), 2 - 20µl	1	€259.08	€145.00
S7100-0221	ErgoOne® Single Channel Pipette (Micro Cone), 2 - 20µl	1	€259.08	€145.00
S7100-0510	ErgoOne® Single Channel Pipette, 0.5 - 10µl	1	€259.08	€145.00
S7100-1100	ErgoOne® Single Channel Pipette, 10 - 100µl	1	€259.08	€145.00
S7100-2200	ErgoOne® Single Channel Pipette, 20 - 200µl	1	€259.08	€145.00
S7108-0510	ErgoOne® 8-Channel Pipette, 0.5 - 10µl	1	€747.15	€445.00
S7108-1100	ErgoOne® 8-Channel Pipette, 10 - 100µl	1	€747.15	€445.00
S7108-3300	ErgoOne® 8-Channel Pipette, 30 - 300µl	1	€747.15	€445.00
S7110-1000	ErgoOne® Single Channel Pipette, 100 - 1000µl	1	€259.08	€145.00
S7112-0510	ErgoOne® 12-Channel Pipette, 0.5 - 10µl	1	€890.79	€535.00
S7112-1100	ErgoOne® 12-Channel Pipette, 10 - 100µl	1	€890.46	€535.00
S7112-3300	ErgoOne® 12-Channel Pipette, 30 - 300µl	1	€890.46	€535.00
S7150-5000	ErgoOne® Single Channel Pipette, 0.5 - 5ml	1	€259.08	€145.00

Equipment

Art. No.	Description	Pack Size	List Price	Your Price
S8010-0000	Starlab Vortex Mixer	1	€421.27	€320.00
S8030-5000	Starlab Minicentrifuge	1	€286.86	€218.00

Please use Q26000928528 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Passionate for science.

