

**Built for your lab,
made to reduce impact.**
**You are passionate
for science.**

SPECIAL OFFER HIGHLIGHTS

Valid until 29.05.2026



**SCAN TO
BUY ONLINE**

**IN THIS EDITION:
PIPETTE TIPS
MANUAL PIPETTES
ERLENMEYER FLASKS
FREEZER RACKING**

 **starlab**

Repeat dispensing.



Made in Germany. "Local for local" production reduces carbon footprint.

Reusable saving on plastic waste.

RepeatOne S®

Art. No. G4026-1500

OFFER: 248 €

Dispenser tips not included.

YOUR FIRST CHOICE FOR CHALLENGING LIQUIDS

With up to 48 dispenses from the same sample, you can use the **RepeatOne® S** to dispense precise volumes of liquid in a series, without the need to aspirate in between each step. 1-year warranty.

TipOne® Repeat dispenser tips.

| Tip volume | Pack Size | Article No. | OFFER* | Tip Volume | Pack Size | Article No. | OFFER* |
|---------------------|-----------|-------------|---------|--------------------------------------|-----------|-------------|----------|
| Non-Sterile | | | | Sterile, Individually Wrapped | | | |
| 0.1 ml | 100 | S4751-0010 | 55.63 € | 0.1 ml | 100 | S4761-0010 | 93.09 € |
| 0.5 ml | 100 | S4751-0050 | 40.21 € | 0.5 ml | 100 | S4761-0050 | 82.34 € |
| 1.0 ml | 100 | S4751-0100 | 40.21 € | 1.0 ml | 100 | S4761-0100 | 82.34 € |
| 2.5 ml | 100 | S4751-0250 | 40.21 € | 2.5 ml | 100 | S4761-0250 | 82.34 € |
| 5.0 ml | 100 | S4751-0500 | 40.21 € | 5.0 ml | 100 | S4761-0500 | 82.34 € |
| 10 ml | 100 | S4751-1000 | 40.21 € | 10 ml | 100 | S4761-1000 | 82.34 € |
| 25 ml ¹⁾ | 100 | S4751-2500 | 69.40 € | 25 ml ¹⁾ | 100 | S4761-2500 | 118.97 € |
| 50 ml ¹⁾ | 100 | S4751-5000 | 82.34 € | 50 ml ¹⁾ | 100 | S4761-5000 | 122.28 € |

¹⁾ 25 ml and 50 ml dispenser tips are supplied with four adapters.
Note: TipOne® Repeat Dispenser Tips are not autoclavable.

ErgoOne® manual pipettes.

MADE FOR DURABILITY!

ErgoOne® manual pipettes provide accuracy and precision you can trust, alongside ergonomic comfort and a stylish look you will love. ErgoOne pipettes have a greased system, which means they rarely need their sealing system changed. Less parts = less waste! 3-year warranty.



Made in Germany. "Local for local" production reduces carbon footprint.

Single-Channel Pipettes

| | | | |
|--|---|---|--|
| S7100-0125 0.1–2.5µl 145 € | S7100-0510 0.5 – 10µl 145 € | S7100-0221 2–20µl Micro 145 € | S7100-0220 2–20µl Standard 145 € |
| S7100-1100 10–100µl 145 € | S7100-2200 20–200µl 145 € | S7110-1000 100–1,000µl 145 € | S7150-5000 500–5,000µl 145 € |

8-Channel Pipettes

| | | |
|---|---|---|
| S7108-0510 0.5–10µl 445 € | S7108-1100 10–100µl 445 € | S7108-3300 30–300µl 445 € |
|---|---|---|



ErgoOne® carousel
Holds up to 8 pipettes. For use with ErgoOne®, Biohit®, Finn®, Gilson® and Rainin®, single and multi-channel pipettes. UV resistant.

Art. No: S7200-0000
94.61 €



ErgoOne® Fast pipette controller.
The ErgoOne® FAST is one of the lightest pipette controllers on the market. Designed for your comfort, it balances perfectly in the hand.

Art. No: S7166-0010
276.01 €

RACKS WITH PULL OUT SHELVES

For use with all **standard storage boxes** up to 136 x 136 x 53 mm.

| Capacity | Format | Ext. Dimensions (WxDxH) | Art. No | OFFER* |
|----------|-----------------|-------------------------|------------|----------|
| 9-Place | 3 deep x 3 high | 140 x 450 x 173 mm | E6001-0033 | 144.26 € |
| 12-Place | 3 deep x 4 high | 140 x 450 x 231 mm | E6001-0034 | 149.40 € |
| 15-Place | 3 deep x 5 high | 140 x 450 x 288 mm | E6001-0035 | 154.87 € |
| 12-Place | 4 deep x 3 high | 140 x 569 x 174 mm | E6001-0043 | 154.87 € |
| 16-Place | 4 deep x 4 high | 140 x 569 x 231 mm | E6001-0044 | 159.69 € |
| 20-Place | 4 deep x 5 high | 140 x 569 x 281 mm | E6001-0045 | 165.15 € |
| 24-Place | 4 deep x 6 high | 140 x 569 x 337 mm | E6001-0046 | 177.36 € |
| 28-Place | 4 deep x 7 high | 140 x 569 x 403 mm | E6001-0047 | 188.92 € |

For use with Starlab's StarStore 100 storage boxes only.

| Capacity | Format | Ext. Dimensions (WxDxH) | Art. No | OFFER* |
|----------|-----------------|-------------------------|------------|----------|
| 9-Place | 3 deep x 3 high | 147 x 500 x 179 mm | E6001-0133 | 161.94 € |
| 12-Place | 3 deep x 4 high | 147 x 500 x 239 mm | E6001-0134 | 168.04 € |
| 16-Place | 4 deep x 4 high | 147 x 655 x 239 mm | E6001-0144 | 173.50 € |

Racking suitable for use with StarStore 100 storage boxes (WxDxH) 140 x 155 x 55 mm.

Freezer racking

SUSTAINABILITY JUST GOT COOLER

- Optimise space in your freezer, reduce energy use and cut your running costs
- Manufactured using 100% sustainable energy
- Manufactured in Germany — shorter supply chain links reduce carbon footprint
- Racks specifically designed for Eppendorf, Thermo, PHCbi, Haier and Nordic Labs freezers

MADE FROM
85% RECYCLED MATERIAL

RACKS WITH SIDE ACCESS



Boxes not included.

For use with **standard storage boxes** up to 136 x 136 x 53 mm.

| Capacity | Format | Ext. Dimensions (WxDxH) | Art. No | OFFER* |
|----------|-----------------|-------------------------|------------|----------|
| 9-Place | 3 deep x 3 high | 139 x 426 x 169 mm | E6000-0033 | 104.10 € |
| 12-Place | 3 deep x 4 high | 139 x 426 x 224 mm | E6000-0034 | 119.20 € |
| 15-Place | 3 deep x 5 high | 139 x 426 x 297 mm | E6000-0035 | 133.02 € |
| 12-Place | 4 deep x 3 high | 139 x 565 x 169 mm | E6000-0043 | 119.20 € |
| 16-Place | 4 deep x 4 high | 139 x 565 x 224 mm | E6000-0044 | 119.20 € |
| 20-Place | 4 deep x 5 high | 139 x 565 x 279 mm | E6000-0045 | 136.23 € |
| 24-Place | 4 deep x 6 high | 139 x 565 x 334 mm | E6000-0046 | 156.15 € |



Boxes not included.



Boxes not included.

RACKING FOR CHEST FREEZERS

For use with all **standard storage boxes** up to 136 x 136 x 53 mm.

| Capacity | Ext. Dimensions (WxDxH) | Art. No | OFFER* |
|------------------------------|-------------------------|------------|----------|
| 9-High SS Rack, Locking Rod | 141 x 148 x 501 mm | E6010-0009 | 89.96 € |
| 10-High SS Rack, Locking Rod | 141 x 148 x 556 mm | E6010-0010 | 89.96 € |
| 11-High SS Rack, Locking Rod | 141 x 148 x 611 mm | E6010-0011 | 95.10 € |
| 12-High SS Rack, Locking Rod | 141 x 148 x 666 mm | E6010-0012 | 100.57 € |
| 13-High SS Rack, Locking Rod | 141 x 148 x 721 mm | E6010-0013 | 106.35 € |

***OFFER REF: Q26000966457**

must be stated on your PO (excl. eShop). Offer valid to 29.05.26 / while stock lasts.

TipOne® sterile refills.



Made from recycled material

Reuse your TipOne® racks with our refills!

PLASTIC SAVING REFILLS

- Fully enclosed sterile refill
- Saves up to 63% plastic compared to racked tips
- Space saving of up to 20% compared to racked tips
- No need to autoclave – I'm fully sterile
- Made in a factory run by sustainable energy
- Recyclable. Join our unique closed-loop recycling service!



| Description | Pack Size | Art. No. | OFFER* |
|---|-----------|------------|---------|
| 10/20µl XL Graduated TipOne® Pipette Tip, Sterile Refill | 10 x 96 | S1110-3710 | 21.75 € |
| 200µl Bevelled TipOne® Pipette Tip, Sterile Refill | 10 x 96 | S1111-1710 | 22.74 € |
| 200µl UltraPoint® Graduated TipOne® Pipette Tip, Sterile Refill | 10 x 96 | S1113-1710 | 25.05 € |
| 300µL, Graduated TipOne® Pipette Tip, Sterile Refill | 10 x 96 | S1110-9710 | 36.92 € |
| 1250µl XL Graduated TipOne® Pipette Tip, Sterile Refill | 10 x 96 | S1112-1730 | 43.90 € |

Win a TipOne® lego set!

Enter our online competition for the chance to win! Scan the QR code or visit the link

Competition closes:
29-05-2026



ENTER HERE

starlab.click/win-lego



Erlenmeyer flasks.

IDEAL WHEN SAFETY OR WEIGHT ARE PRIORITIES

These optically clear, polycarbonate flasks are safer, and visibly clearer than glass.

- Sterile (SAL 10⁻⁶), individually-wrapped
- Shatter resistant
- Dual-function cap featuring 0.22 µm PTFE membrane means you can confidently use the same flask for both aerobic cultures (vented) or anaerobic cultures (closed)
- No more foil or bungs that can breed contamination




REUSABLE! Withstands at least 10 autoclave cycles!

| | Description | Pack Size | Art. No. | OFFER* |
|---|--|-----------|------------|----------|
| FLAT BASE  | 125 ml Erlenmeyer Flask, Flat Base | 24 | E5000-1012 | 136.88 € |
| | 250 ml Erlenmeyer Flask, Flat Base | 12 | E5000-1025 | 90.32 € |
| | 500 ml Erlenmeyer Flask, Flat Base | 12 | E5000-1050 | 118.37 € |
| | 1000 ml Erlenmeyer Flask, Flat Base | 6 | E5000-1100 | 82.75 € |
| | 2000 ml Erlenmeyer Flask, Flat Base | 6 | E5000-1200 | 104.63 € |
| | 3000 ml Erlenmeyer Flask, Flat Base | 3 | E5000-1300 | 125.38 € |
| BAFFLED BASE  | 125 ml Erlenmeyer Flask, Baffled Base | 24 | E5001-1012 | 136.88 € |
| | 250 ml Erlenmeyer Flask, Baffled Base | 12 | E5001-1025 | 90.32 € |
| | 500 ml Erlenmeyer Flask, Baffled Base | 12 | E5001-1050 | 118.37 € |
| | 1000 ml Erlenmeyer Flask, Baffled Base | 6 | E5001-1100 | 82.75 € |
| | 2000 ml Erlenmeyer Flask, Baffled Base | 6 | E5001-1200 | 104.63 € |
| | 3000 ml Erlenmeyer Flask, Baffled Base | 3 | E5001-1300 | 125.38 € |

Reagent reservoirs.

Autoclavable and reusable, reducing plastic waste. Compatible with all multi-channel pipettes.

| | Description | Pack Size | Article No. | OFFER* |
|---|---|-----------|-------------|---------|
|  | 50/60 ml StarLab, PP Reservoir. Non-sterile. V-shaped bottom. | 10 | S4026-5806 | 12.96 € |

Made in Germany. "Local for local" production reduces carbon footprint.

***OFFER REF: Q26000966457**
must be stated on your PO (excl. eShop).
Offer valid to 29.05.26 / while stock lasts.