

**Trusted for consistency,
designed for growth.**
You are passionate
for science.

SPECIAL OFFER HIGHLIGHTS

Valid until 31.08.2026



SCAN TO
BUY ONLINE

**IN THIS EDITION:
FLASKS, DISHES & PLATES
BOTTLE TOP FILTRATION
ASPIRATOR
TUBES & VIALS FOR CELL CULTURE**



CytoOne® Dishes. FROM: 41,86 €*



CytoOne® TC-treated dishes feature stacking rings with vented design for thermal equalisation. Untreated lids reduce condensation: ideal for microscopy.

Size	Pack Size	Art. No.	OFFER*
TC-Treated Dishes			
35 x 10 mm	30 x 10	CC7682-3340	43,90 €
60 x 15 mm	30 x 10	CC7682-3359	41,86 €
100 x 20 mm	30 x 10	CC7682-3394	92,15 €
150 x 20 mm	12 x 5	CC7682-3614	54,36 €

Certified and sterile (SAL 10⁻⁹).

All CytoOne® products are certified RNase, DNase, DNA & Pyrogen free and Gamma sterilised (SAL 10⁻⁹). TC-Treated for optimum cell attachment and growth.

CytoOne® Plates. FROM: 40,99 €*



Complete 360° chimney well design prevents cross-contamination and ensures a uniform temperature transfer across all wells to prevent the edge effect. Integrated vents and stacking rims create space for continuous airflow between plates and reduces evaporation.

Wells	Pack Size	Art. No.	OFFER*
TC-Treated Plates			
6-Well	50	CC7682-7506	40,99 €
12-Well	50	CC7682-7512	54,65 €
24-Well	50	CC7682-7524	47,67 €
48-Well	50	CC7682-7548	59,01 €
96-Well	50	CC7682-7596	56,98 €

Certified and sterile (SAL 10⁻⁹).



CytoOne® Flasks. FROM: 64,26 €*

All flasks supplied in convenient 'mini-packs' in easy-open, resealable zip top bags that protect unused flasks.



Vented filter caps with a 0.2 µm hydrophobic membrane for gas exchange and to exclude bacteria and fungal contaminants.

Two-position plug seal caps allow you to choose secure manual venting or an airtight seal.

Size	Pack Size	Art. No.	OFFER*
TC Treated Flasks, Vented Filter Caps			
T-25	30 x 10	CC7682-4825	226,77 €
T-75	20 x 5	CC7682-4875	134,23 €
T-150	8 x 5	CC7682-4815	87,11 €
T-225	5 x 5	CC7682-4822	84,25 €

Certified and sterile (SAL 10⁻⁹).

Size	Pack Size	Art. No.	OFFER*
TC Treated Flasks, Plug Seal Caps			
T-25	30 x 10	CC7682-4325	132,23 €
T-75	20 x 5	CC7682-4175	90,54 €
T-150	8 x 5	CC7682-4415	71,11 €
T-225	5 x 5	CC7682-4225	64,26 €

Certified and sterile (SAL 10⁻⁹).

Centrifuge Tubes.

Safeguard your cell viability. Designed to support your most sensitive cell work. Graduated to full 15 ml and 50 ml volumes, with extra graduations on the conical base.

Extra moulded graduations!
15 ml Conical Tube: 0 – 1.0 ml
50 ml Conical Tube: 0 – 3.0 ml

Description	Pack Size	Art. No.	OFFER*
15 ml Conical, Loose	20 x 25	E1415-0100	55,81 €
15 ml Conical, Racked	20 x 25	E1415-0700	67,32 €
50 ml Conical, Loose	20 x 25	E1450-0100	72,68 €
50 ml Conical, Racked	20 x 25	E1450-0700	83,94 €
50 ml Skirted, Loose	20 x 25	E1450-0400	85,07 €

Certified and sterile (SAL 10⁻⁹).

MICROTUBE PESTLE for 1.5/2.0 ml tubes. Assorted colours. Certified. Pack 10. I1415-5390 7,04 €*



CULTURE TUBES
14 ml tubes with dual-position caps that can be left loose for aerobic work or sealed for anaerobic culture.

Polypropylene tube -190 °C / 3000 xg. 20 x 25-Tube bags I1485-0810 140,91 €*

Polystyrene tube -10 °C / 1400 xg 20 x 25-Tube bags I1485-2810 136,55 €*



Serological Pipettes.



Vol.	Qty	Art. No.	OFFER*
1.0 ml	500	E4860-0001	42,98 €
2.0 ml	500	E4860-0002	51,23 €
5. ml	200	E4860-0005	28,74 €
10 ml	200	E4860-0010	36,99 €
25 ml	150	E4860-0025	45,98 €
50 ml	100	E4860-0050	51,23 €
100 ml	50	E4860-1100	66,72 €

Certified RNase, DNase, DNA, Latex and Pyrogen free. Sterile (SAL 10⁻⁹). Individually wrapped.

Aspirator. Includes 4L PP Bottle & Hand Controller. Art. No. N2400-9000 OFFER: 1.466,62 €*

Aspiration Pipettes.

Volume	Qty	Art. No.	OFFER*
2.0 ml	500	E4861-0002	59,48 €
5.0 ml	200	E4861-0005	27,74 €
10 ml	200	E4861-0010	35,49 €

Certified RNase, DNase, DNA, Latex and Pyrogen free. Sterile (SAL 10⁻⁹). Individually wrapped.

Save time and effort during lab workflows with the Starlab Aspirator. The powerful pump can operate two hand controllers without loss of performance.



Replaceable rechargeable battery!



PIPETTE CONTROLLER. Pressure-responsive speed controls make the ErgoOne® FAST ideal for cell culture work such as resuspension of bacteria or cell pellets, nucleic acid precipitates and aspiration of a cell layer.

ErgoOne® FAST Art. No. S7166-0010 OFFER: 262,21 €

Cryogenic Vials.

Supplied in convenient, self-standing resealable pouches. Certified. Sterile (SAL 10⁻⁹).



	Vol:	1.0 ml	1.8 ml	1.8 ml	3.6 ml	3.6 ml	4.5 ml	4.5 ml
	Base:	Skirted	Conical	Skirted	Conical	Skirted	Conical	Skirted
EXTERNAL Threaded Tube Self-Sealing Cap, Sterile	E3090-6212	E3090-6321	E3090-6222	n/a	n/a	n/a	E3090-6242	
	10 x 50	10 x 50	10 x 50				6 x 50	
	108,00 €*	116,36 €*	100,84 €*				80,36 €*	
INTERNAL Threaded Tube, Silicone Seal in Cap, Sterile	E3110-6112	E3110-6121	E3110-6122	E3110-6131	E3110-6132	E3110-6141	E3110-6142	
	10 x 50	10 x 50	10 x 50	8 x 50	8 x 50	6 x 50	6 x 50	
	126,20 €*	141,44 €*	118,45 €*	103,71 €*	117,20 €*	108,21 €*	108,21 €*	

! *OFFER REF: Q26000992730 must be stated on your PO (excl. eShop). Offer valid to 31.08.26 / while stock lasts.

ErgoOne® Starter Packs.



ErgoOne® pipettes paired with Starlab's trusted TipOne® pipette tips create the perfect pipetting partnership!

3 x ErgoOne® Single-Channel Pipettes: 20 µl, 200 µl, 1000 µl
3 x racks of sterile TipOne® Tips
3 x Shelf-Adhesive Pipette Holders
1 x Laboratory Timer

STARTERPACK
382,70 €

4 x ErgoOne® Single-Channel Pipettes: 2.5 µl, 20 µl, 200 µl, 1000 µl,
3 x racks of sterile TipOne® Tips
4 x Shelf-Adhesive Pipette Holders
1 x Laboratory Timer

STARTERPACK42.5
511,70 €

4 x ErgoOne® Single-Channel Pipettes: 10 µl, 20 µl, 200 µl, 1000 µl,
3 x racks of sterile TipOne® Tips
4 x Shelf-Adhesive Pipette Holders
1 x Laboratory Timer

STARTERPACK410
511.70 €

Cell Scrapers.

Gently raise/lift the cells with Starlab cell scrapers. The silicone blades of 220 mm and 250 mm scrapers prevent cell lysing. The flexible ABS handles on all scrapers provides easy handling inside flasks. Individually wrapped. Certified. Sterile.



220 mm straight handle, 11 mm fixed blade
CC7600-0220
Pack Size: 100
SAL 10⁶ **141,93 €**



250 mm angled handle, 11 mm pivoting blade
CC7600-0250
Pack Size: 100
SAL 10⁶ **135,20 €**



300 mm straight handle, 20 mm pivoting blade
CC7600-0300
Pack Size: 150
SAL 10⁶ **114,72 €**

Win a TipOne® lego set!

Enter our online competition for the chance to win!

Scan the QR code or visit the link.

Competition closes: **31.08.2026**



starlab.click/win-lego




***OFFER REF: Q26000992730** must be stated on your PO (excl. eShop). Offer valid to 31.08.26 / while stock lasts.

Bottle Top Filtration.

Flexible system! Choose complete full assembly units for ready-to-use filtration, GL45 compatible bottle top funnels to pair with your own bottles or media bottles for versatile storage and preparation.

Superior filtration performance!



Bottle Top Filtration System					
	150 ml	250 ml	500 ml	1000 ml	
Full Assembly Unit					
	0,22 µm PES membrane Pack size: 12	CC8221-3226 60,64 €* 	CC8221-4226 67,74 €* 	CC8221-5226 85,49 €* 	CC8221-6226 136,66 €*
Funnel Tops only (GL45 compatible)					
	0,22 µm PES membrane Pack size: 24	CC8222-3226 69,32 €* 	CC8222-4226 77,07 €* 	CC8222-5226 97,48 €* 	CC8222-6226 155,54 €*
Media Bottles only					
	Bottle and Cap Pack size: 12	CC8221-3000 49,88 €* 	CC8221-4000 56,30 €* 	CC8221-5000 62,42 €* 	CC8221-6000 97,92 €*

All products are certified RNase, DNase and Pyrogen free, are sterile (SAL 10⁻⁶), and supplied individually wrapped.

THE PERFECT PARTNER

Starlab Aspirator. N2400-9000.
1.466,62 €*



Syringe Filters.

Designed for efficiency and safety. Low hold-up volume means you can recover more of your valuable solution. 33 mm diameter filters supplied in easy-open blister packs showing sterile validity date (SAL 10⁻⁶).

100 Filters: Double the usual pack size!



0.22 µm PES Sterile.
E4780-1226
Pack Size: 100
73,83 €*



0.22 µm PVDF Sterile.
E4780-1221
Pack Size: 100
73,83 €*



0.22 µm CA Sterile.
E4780-1223
Pack Size: 100
73,83 €*

Cell Strainers.

Cell Strainers provide a simple and reliable way to achieve uniform cell suspensions. Designed to fit perfectly into 50 ml centrifuge tubes. Gamma sterilised (SAL 10⁻⁶). Individually wrapped.



40 µm mesh. Sterile.
CC8111-0042
Pack Size: 50
67,73 €*



70 µm mesh. Sterile.
CC8111-0072
Pack Size: 50
67,73 €*



100 µm mesh. Sterile.
CC8111-0102
Pack Size: 50
67,73 €*