


# Rigid Frame PCR Plates.

- Polycarbonate frame for precise automated handling
- Ultra-thin PP wells for optimal heat transfer, high reaction efficiency
- Ridged wells prevent evaporation when sealed
- White printed matrix
- RNase, DNase, DNA & PCR inhibitor free



| Frame Colour                 | Well Colour | Cut Corner | Approx. Max. Well Capacity | Pack Size | Article No. | OFFER   |
|------------------------------|-------------|------------|----------------------------|-----------|-------------|---------|
| 96-WELL SEMI-SKIRTED         |             |            |                            |           |             |         |
| Natural                      | Natural     | A12        | 300 µl                     | 10        | E1404-0200  | 30,63 € |
| Blue                         | Natural     | A12        | 300 µl                     | 10        | E1404-0201  | 30,63 € |
| 96-WELL SKIRTED, LOW-PROFILE |             |            |                            |           |             |         |
| Natural                      | Natural     | H1         | 200 µl                     | 10        | E1404-0100  | 30,63 € |
| Blue                         | Natural     | H1         | 200 µl                     | 10        | E1404-0101  | 30,63 € |
| 384-WELL SKIRTED             |             |            |                            |           |             |         |
| Natural                      | Natural     | A24        | 40 µl                      | 10        | E1404-0300  | 37,33 € |
| Blue                         | Natural     | A24        | 40 µl                      | 10        | E1404-0301  | 37,33 € |

## Plate sealing.



| Description                                        | Art. No.   | OFFER*   |
|----------------------------------------------------|------------|----------|
| Self-Adhesive Plate Seal, Polypropylene, Clear     | E2796-0793 | 43,35 €  |
| Self-Adhesive Sealing Strips, Polypropylene, Clear | E2796-2850 | 38,82 €  |
| Aluminium sealing film                             | E2796-0792 | 72,82 €  |
| Clear Polypropylene Seal (PCR)                     | E2796-0793 | 43,35 €  |
| Peelable Heat Sealing Film                         | E2796-3020 | 244,33 € |
| Easi-Peel Heat Sealing Film PP                     | E2796-3060 | 164,29 € |
| Easy Pierce Film (PCR)                             | E2796-3090 | 164,29 € |
| Sealing Film for qPCR, self sealing                | E2796-9895 | 193,09 € |
| Sealing Foil, Aluminium                            | E2796-1100 | 89,07 €  |

## Black Friday week is coming. Be first to know!

A full week of savings across equipment and more is on the way.

Scan the code and we'll send you the discounts before they go live.



[starlab.click/blackfriday](https://starlab.click/blackfriday)

## Storage offers.



Cardboard storage boxes



Polpropylene storage boxes



Cryogenic storage boxes

Explore the full offer range online.

Scan for full details!




[starlab.click/offers](https://starlab.click/offers)

## Freezer racking

### SUSTAINABILITY JUST GOT COOLER

- Freezer racking helps optimise capacity by keeping freezers tidy and maximising storage space for efficient energy use
- High quality racking prevents the need for regular repair or exchange
- Manufactured in a factory powered by 100% sustainable energy
- Made from 85% recycled stainless steel
- Made in Germany. Short journey from factory to customer reduces carbon footprint



**\*OFFER REF: Q26001021191**

must be stated on your PO (excl. eShop).  
Offer valid to 16.12.26 / while stock lasts.

# Save on equipment, spend more on research. You are passionate for science.

## SPECIAL OFFER HIGHLIGHTS


Valid until 16.12.2026




IN THIS EDITION:  
LAB EQUIPMENT  
PIPETTES  
SAMPLE STORAGE  
PCR PLATES AND SEALS  
FREEZER RACKING



# Our lowest ever priced equipment for your lab.



| Description                                                     | Art. No.   | OFFER*    |
|-----------------------------------------------------------------|------------|-----------|
| MIXING, ROLLING, STIRRING AND SHAKING.                          |            |           |
| Magnetic Stirrer                                                | N2400-2010 | 134,26 €  |
| Magnetic Stirrer with Heating                                   | N2400-3010 | 260,02 €  |
| Rotator Mixer Multi-1, Programmable (incl. Platform N2400-5104) | N2400-5000 | 595,55 €  |
| Tube Roller                                                     | N2400-7010 | 433,00 €  |
| Rocking Shaker                                                  | N2400-8010 | 852,38 €  |
| 3D Rocking Shaker                                               | N2400-8020 | 836,59 €  |
| Orbital Shaker (a platform must be purchased separately)        | N2400-8030 | 749,48 €  |
| Microplate Mixer                                                | N2400-8040 | 302,23 €  |
| Plate Centrifuge                                                | N2631-0100 | 690,35 €  |
| Vortex Mixer, 100-240 V                                         | S8010-0000 | 307,53 €  |
| Vortex Mixer, Standard attachment                               | S8010-0010 | 6,16 €    |
| Vortex Mixer, Adapter for 1.5/2.0 ml tubes                      | S8010-0011 | 16,73 €   |
| Vortex Mixer, Adapter for microplates                           | S8010-0012 | 10,17 €   |
| Minicentrifuge                                                  | S8030-5000 | 209,40 €  |
| Minicentrifuge, tube rotor, 8 place                             | S8030-5010 | 19,46 €   |
| Minicentrifuge, rotor for PCR tubes                             | S8030-5011 | 19,46 €   |
| GUARDONE' WORKSTATION.                                          |            |           |
| GuardOne laminar flow 32in, Starlab                             | S8040-0000 | 3068,48 € |
| GuardOne laminar flow 48in, Starlab                             | S8040-0001 | 4089,74 € |
| GuardOne dead air 32in, Starlab                                 | S8040-1000 | 2582,36 € |
| GuardOne dead air 48in, Starlab                                 | S8040-1001 | 3275,72 € |
| LABEL PRINTER.                                                  |            |           |
| Brady® M210-LAB Printer                                         | N3000-2000 | 118,43 €  |



- Ergonomic and lightweight
- Minimal operating forces
- Volume lock
- Universal tip fit (excl. 5.0 ml)
- 3 year warranty



- Easy to use, intuitive menu
- Lightweight, ergonomic design
- Universal tip fit
- 12-month warranty

## Manual pipettes.

### Single-Channel Manual Pipettes


|                                        |                                       |                                           |                                              |
|----------------------------------------|---------------------------------------|-------------------------------------------|----------------------------------------------|
| S7100-0125<br>0.1 – 2.5 µl<br>142,49 € | S7100-0510<br>0.5 – 10 µl<br>142,49 € | S7100-0221<br>2 – 20 µl Micro<br>142,49 € | S7100-0220<br>2 – 20 µl Standard<br>142,49 € |
| S7100-1100<br>10 – 100 µl<br>142,49 €  | S7100-2200<br>20 – 200 µl<br>142,49 € | S7110-1000<br>100 – 1,000 µl<br>142,49 €  | S7150-5000<br>500 – 5,000 µl<br>142,49 €     |

### 8-Channel Pipettes

|                                       |                                       |                                       |
|---------------------------------------|---------------------------------------|---------------------------------------|
| S7108-0510<br>0.5 – 10 µl<br>388,52 € | S7108-1100<br>10 – 100 µl<br>388,52 € | S7108-3300<br>30 – 300 µl<br>388,52 € |
|---------------------------------------|---------------------------------------|---------------------------------------|

### 12-Channel Pipettes

|                                       |                                       |                                       |
|---------------------------------------|---------------------------------------|---------------------------------------|
| S7112-0510<br>0.5 – 10 µl<br>445,23 € | S7112-1100<br>10 – 100 µl<br>445,23 € | S7112-3300<br>30 – 300 µl<br>445,23 € |
|---------------------------------------|---------------------------------------|---------------------------------------|



**RepeatOne S®**  
Dispenses up to 48 precise volumes from a single sample.  
1-year warranty.

Art. No: G4026-1500  
255,48 €



**ErgoOne® carousel**  
Holds up to 8 pipettes.  
Art. No: S7200-0000  
88,93 € / each

**ErgoOne Wall Mount**  
Art. No: S7200-0001  
4,96 € / each

**Acrylic Stand**  
Holds up to 5 ErgoOne Pipettes  
Art. No: E3099-0000  
36,21 € / each



**ErgoOne® Fast pipette controller.**  
The ErgoOne® FAST is one of the lightest pipette controllers on the market.

Art. No: S7166-0010  
276,01 €

## Electronic pipettes.

### Single-Channel ErgoOne® E Pipettes

|                                       |                                       |                                       |                                         |                                          |
|---------------------------------------|---------------------------------------|---------------------------------------|-----------------------------------------|------------------------------------------|
| G9001-0010<br>0.5 – 10 µl<br>398,36 € | G9001-0020<br>2.0 – 20 µl<br>398,36 € | G9001-0200<br>10 – 200 µl<br>398,36 € | G9001-1000<br>50 – 1,000 µl<br>398,36 € | G9001-5000<br>250 – 5,000 µl<br>398,36 € |
|---------------------------------------|---------------------------------------|---------------------------------------|-----------------------------------------|------------------------------------------|

### 8-Channel ErgoOne® E Pipettes

|                                       |                                       |                                        |                                       |                                        |
|---------------------------------------|---------------------------------------|----------------------------------------|---------------------------------------|----------------------------------------|
| G9008-0010<br>0.5 – 10 µl<br>690,43 € | G9008-0020<br>1.0 – 20 µl<br>690,43 € | G9008-0100<br>5.0 – 100 µl<br>690,43 € | G9008-0300<br>15 – 300 µl<br>690,43 € | G9008-1250<br>50 – 1250 µl<br>690,43 € |
|---------------------------------------|---------------------------------------|----------------------------------------|---------------------------------------|----------------------------------------|

### 12-Channel ErgoOne® E Pipettes

|                                       |                                       |                                        |                                       |                                        |
|---------------------------------------|---------------------------------------|----------------------------------------|---------------------------------------|----------------------------------------|
| G9012-0010<br>0.5 – 10 µl<br>836,09 € | G9012-0020<br>1.0 – 20 µl<br>836,09 € | G9012-0100<br>5.0 – 100 µl<br>836,09 € | G9012-0300<br>15 – 300 µl<br>836,09 € | G9012-1250<br>50 – 1250 µl<br>836,09 € |
|---------------------------------------|---------------------------------------|----------------------------------------|---------------------------------------|----------------------------------------|

**\*OFFER REF: Q26001021191**  
must be stated on your PO (excl. eShop).  
Offer valid to 16.12.26 / while stock lasts.