

More than just tips. Starlab special offers.

Valid to 31 May 2022

IN THIS ISSUE

Equipment

Multifunctional devices for a variety of applications.

Consumables

Everything you need for your multi-functional lab devices.

Sample Storage

Storage solutions for a range of applications.

Pipettes

ErgoOne® pipettes for accuracy and precision all day.

Protection

StarGuard® nitrile gloves for daily protection.



Passionate for science



Contents of the Starlab Special Offers.

MORE THAN JUST TIPS...

Products to help you in a range of different applications.

EQUIPMENT




Mixer HC
page 4



Mini Fuge PLUS
page 8



Plate Centrifuge
page 10



Microplate Mixer
page 12



GuardOne® Workstation
page 6



Vortexers
page 7




Rotator Mixer
page 15



Magnetic Stirrers
page 14

GLOVES



Nitrile Gloves
pages 18 & 19

CONSUMABLES



Microcentrifuge Tubes
page 5



Preparation Tubes
page 5



Screw Cap Microtubes
page 9




PCR Plate Seals
page 11



PlateOne® Plates
page 13

PIPETTES



ErgoOne® Pipettes
page 16



ErgoOne® FAST Controller
page 17

SAMPLE STORAGE



Multi-Tube Racks
page 7



Freezer storage / cooling
page 15

The Mixer HC.

ONE SMALL DEVICE, SO MANY USES.

IT'S FOOTPRINT IS ROUGHLY THE SIZE OF AN iPad!

The Mixer HC is a smart bit of kit that does the work of several other pieces of equipment in the lab. Easy to use and to program, the Mixer HC will become an invaluable benchtop companion to help speed up your laboratory tasks. You will wonder how you ever managed without one!



LEARN MORE

starlab.click/mixer-hc

Mixer HC

Art. No. **S8012-0000** (Excludes Thermoblock)

OFFER: 1550 €



NINE THERMOBLOCKS TO COVER ALL APPLICATION NEEDS

Versatile! Use the Mixer HC with all popular tubes and plates.



24 x 0.5 ml
Microcentrifuge Tubes
(Includes Iso-Rack)
S8012-0010
OFFER: 315 €



24 x 1.5 ml
Microcentrifuge Tubes
(Includes Iso-Rack)
S8012-0011
OFFER: 315 €



24 x 2.0 ml Microcentrifuge
Tubes, Round Bottoms
(Includes Iso-Rack)
S8012-0012
OFFER: 315 €



24 x 2.0 ml
Microcentrifuge Tubes
Ø 12 mm, V-Bottoms
S8012-0013
OFFER: 385 €



8 x 5.0 ml
Preparation Tubes
S8012-0014
OFFER: 315 €



8 x 15 ml
Conical Centrifuge Tubes
S8012-0015
OFFER: 315 €




4 x 50 ml
Conical Centrifuge Tubes
S8012-0016
OFFER: 315 €



24 x 1.5/2.0 ml
Cryo Tubes
S8012-0017
OFFER: 315 €



Micro/Deepwell Plate,
with Lid
S8012-0018
OFFER: 385 €



Adapter Plate 96 x 0.2 ml
Tubes or a PCR plate
(Use with S8012-0018)
S8012-0019
OFFER: 135 €



ISO-Rack: 24 x 0.5 ml
Microcentrifuge Tubes
S8012-0020
OFFER: 45 €



ISO-Rack: 24 x 1.5/2.0 ml
Microcentrifuge Tubes
S8012-0021
OFFER: 45 €



Tubes for use with the Mixer HC.

WELL, JUST SOME OF THE VAST RANGE ...

USE THE MIXER HC FOR:

- Purification of plasmids, DNA, RNA
- cDNA synthesis
- Enzyme reactions (restriction digest, ligations, proteinase-k digestion)
- Denaturation of DNA, RNA, proteins
- Labelling of DNA, RNA, proteins
- Cultivation of bacteria in microcentrifuge or centrifuge tubes & deepwell plates
- Transformation of bacteria and plasmids
- Resuspension of pellets

THE MIXER HC: WHAT OUR CUSTOMERS ARE SAYING:

"Brilliant product for use during lots of different laboratory incubations and steps. Love the small footprint. The range of blocks means I can use 15 ml tubes, 1.0 ml tubes and plates on the one machine. Great price, and lots of uses in the laboratory."

SL, University of Cardiff



TOP TIP!
Convenient mini packs (5x100) decreases contamination risk

LEARN MORE

starlab.click/tubeone



LEARN MORE

starlab.click/preparationtubes

TubeOne® MICROCENTRIFUGE TUBES

TubeOne® stays sealed during boiling, freezing or centrifuging, and are still easy to open afterwards! Made from 100% virgin polypropylene. Certified.

Volume	Colour	Pack Size	Article No.	OFFER*
0.5 ml	Natural	5 x 200	S1605-0000	18 €
0.5 ml	Mixed ¹⁾	5 x 200	S1605-0099	23 €
1.5 ml	Natural	500	S1615-5500	8 €
1.5 ml	Natural	5 x 100	S1615-5550	9 €
1.5 ml	Mixed ¹⁾	5 x 100	S1615-5599	11 €
2.0 ml	Natural	500	S1620-2700	22 €
2.0 ml	Mixed ¹⁾	5 x 100	S1620-2799	22 €

¹⁾ Mixed packs contain the following colours: blue, green, red, violet and yellow.

5.0 ml PREPARATION TUBES

Need to prepare samples up to 5.0 ml? Don't waste excessive plastic and use 15 ml tubes, use our certified 5.0 ml preparation tubes!

Colour	Pack Size	Article No.	OFFER*
Natural	200	E1450-1100	21 €
Blue	200	E1450-1101	21 €
Green	200	E1450-1102	21 €
Red	200	E1450-1104	21 €
Yellow	200	E1450-1106	21 €
Black	200	E1450-1108	21 €

CERTIFIED



All-in for your applications. The GuardOne® Workstation.

Starlab's GuardOne® Workstation provides the ideal workspace for sensitive and critical applications requiring a clean environment or protection from potential contamination by bacteria, DNA or RNA. The GuardOne® Workstation minimizes cross-contamination risk to protect your valuable samples and reagents, thus helping to prevent both the time and cost involved by having to duplicate contaminated work.

- Safe: Efficient UV_{254 nm} decontamination with timer function
- Reduced waste: UV lamp lifespan of up to 9000 hrs
- Easy to use: Glove compatible touchscreen
- Smart: Consumables life-span monitoring
- 2-year warranty

GuardOne® Workstation Laminar Flow 32"

Art. No. S8040-0000

OFFER: 2850 €



BUY ONLINE

starlab.click/guardone

ALSO AVAILABLE



Laminar Flow Workstation 48"
Art. No. S8040-0001

OFFER: 3750 €



Dead Air Workstation 48"
Art. No. S8040-1001

OFFER: 3050 €



Dead Air Workstation 32"
Art. No. S8040-1000

OFFER: 2450 €



For an organised workstation.

A range of versatile multi-tube racks that are **ideal for use in a PCR workstation** or on the bench. Made from polypropylene which provides excellent chemical resistance, all racks are autoclavable, with a temperature use range of 90 °C to 121 °C.



TWISSTAR RACK

Four different coloured polypropylene sections that can be locked together to suit your needs. Holds an array of different tube sizes (0.2 ml to 50 ml), and a PCR plate. Autoclavable.



COMBI-MINI TUBE RACK

80-place polypropylene rack where each side holds different sized tubes: 4 x 50 ml, 12 x 15 ml, 32 x 1.5/2.0 ml and 32 x 0.5 ml. Autoclavable.

PCR LOCK-RACK

Polypropylene racks which lock together to make larger racks. Transparent lid fits over two sections locked together. Each section holds 24 x 1.5 / 2.0 ml tubes, 16 x 0.5 ml tubes and 60 x 0.2 ml tubes. Autoclavable.



LEARN MORE

starlab.click/samplestorage

Art. No. E2344-1099

OFFER: 29 €

Art. No. E2345-1099

Pack of 5. Mixed Neon

OFFER: 25 €

Art. No. E2344-2099

OFFER: 40 €

SIMPLE AND SMART!

Choose between the economical Vortex with speeds up to 2500 rpm, or the Vortex with optional IR motion detector and speeds up to 3000 rpm. Which ever you choose, you know either of these compact units will fit nicely in your GuardOne®.



Vortex

- Variable speed control up to 2500 rpm
- Continuous and touch modes
- Wide range of accessories
- 2-year warranty

N2400-6010

OFFER: 140 €

BUY ONLINE

starlab.click/vortex



Vortex IR

- Variable speed control up to 3000 rpm
- Continuous and sensor modes
- Sensor automatically starts vibration
- Wide range of accessories
- 2-year warranty

N2400-0001

OFFER: 190 €

BUY ONLINE

starlab.click/vortexir

Mini Fuge PLUS.

EVERY BENCH SHOULD HAVE ONE!

Starlab's Mini Fuge Plus is supplied with two rotors: an 8-place rotor for 0.2, 0.5, 1.5 and 2.0 ml microcentrifuge tubes and screw caps tubes, and a 32-place rotor for 0.2 ml PCR tubes (32 singles or 4 x 8-strips).

- 6000 rpm / 2000 x g
- Quick and easy rotor change
- Small footprint

Mini Fuge PLUS

Art. No. N2631-0017

OFFER: 145 €



USE ME WITH PCR TUBES & STRIPS

Go online to see our range of PCR tubes and strips.

PCR TUBES & STRIPS

starlab.click/pcr

USE ME FOR:

- The quick spin down of condensate, vortexed or thawed samples

USE ME WITH

- 8 x 0.5 / 1.5 / 2.0 ml Microcentrifuge Tubes (see page 5)
- 8 x Screw Cap Microtubes (see next page)
- 32 x 0.2 ml PCR Tubes
- 4 x 0.2 ml PCR Strips



BUY ONLINE

starlab.click/minifuge-plus

Screw Cap Microtubes

MIX AND MATCH, TUBES AND CAPS

TUBES

																			
0.5 ml Ribbed, Skirted Tube		0.5 ml Ribbed, Skirted Tube		1.5 ml Ribbed Skirted Tube		1.5 ml Ribbed Skirted Tube		1.5 ml Graduated Conical Tube		1.5 ml Graduated Conical Tube		2.0 ml Graduated Conical Tube		2.0 ml Graduated Conical Tube		2.0 ml Graduated Skirted Tube		2.0 ml Graduated Skirted Tube	
Natural		Amber		Natural		Amber		Natural		Amber		Natural		Amber		Natural		Amber	
E1405-2120		E1405-2127		E1415-2220		E1415-2227		E1415-2230		E1415-2237		E1420-2330		E1420-2337		E1420-2340		E1420-2347	
OFFER: 14 €		OFFER: 18 €		OFFER: 14 €		OFFER: 14 €		OFFER: 14 €		OFFER: 15 €		OFFER: 13 €		OFFER: 15 €		OFFER: 14 €		OFFER: 13 €	

- Routinely tested to 20.000 x g
- All tubes and caps autoclavable
- Made from 100 % medical grade polypropylene
- Certified RNase/DNase and DNA free.
- Tested Endotoxin free

STANDARD CAPS

								
Blue	Green	Orange	Red	Violet	Yellow	Amber	Black	White
E1480-0101	E1480-0102	E1480-0103	E1480-0104	E1480-0105	E1480-0106	E1480-0107	E1480-0108	E1480-0109
OFFER: 27 €	OFFER: 25 €	OFFER: 25 €	OFFER: 27 €	OFFER: 25 €	OFFER: 25 €	OFFER: 25 €	OFFER: 27 €	OFFER: 25 €

TETHERED CAPS

					
Natural	Blue	Green	Violet	Amber	Black
E1480-0200	E1480-0201	E1480-0202	E1480-0205	E1480-0207	E1480-0208
OFFER: 25 €	OFFER: 27 €	OFFER: 27 €	OFFER: 27 €	OFFER: 25 €	OFFER: 25 €

EASY GRIP CAPS

							
Natural/White	White/Blue	White/Green	White/Orange	White/Red	White/Violet	White/White	Mixed
E1480-0300	E1480-0301	E1480-0302	E1480-0303	E1480-0304	E1480-0305	E1480-0309	E1480-0399
OFFER: 25 €	OFFER: 27 €	OFFER: 25 €	OFFER: 25 €	OFFER: 25 €	OFFER: 27 €	OFFER: 27 €	OFFER: 25 €

Plate Centrifuge.

YOU SPIN ME ROUND!

Sealed plates are loaded vertically and surface tension keeps samples in place so there is no danger of spilling samples. The Plate Centrifuge accepts all plates up to 25 mm high and is less than quarter the size of most plate centrifuges.

FEATURES

- 2500 rpm / 500 x g
- Takes any plate up to 25 mm high
- Rapid start to stop

USE ME FOR:

- Safe and quick spin down of droplets and condensation
- Spinning before and after thermal cycling to increase PCR yield.

Plate Centrifuge

Art. No. N2631-0008

OFFER: 515 €



BUY ONLINE

starlab.click/platecentrifuge



STARLAB PCR PLATES

Discover our range of high quality PCR plates: skirted, semi-skirted & non-skirted, including white plates for qPCR. Our range includes both 96- and 384-well plates to fit a wide range of thermal cyclers, including FAST.

LEARN MORE

starlab.click/pcrplates

USE ME WITH 1/2 PLATES, STRIPS & TUBES TOO

An adapter plate (N2631-0018) is available separately which allows you to use the plate centrifuge with 48-well plates or PCR tubes or strips.

PCR Plate Seals.




Aluminium StarSeal (PCR)
Certified RNase, DNase, DNA
& Pyrogen free

Art. No. E2796-9792
100 self-adhesive seals
OFFER: 35 €




**Aluminium Sealing Film
60 µm (PCR)**

Art. No. E2796-0792
100 self-adhesive seals
OFFER: 55 €



Clear Polyolefin StarSeal (PCR)
Certified RNase, DNase, DNA
& Pyrogen free

Art. No. E2796-9793
100 self-adhesive seals
OFFER: 35 €



**X-Clear Advanced Polyolefin
StarSeal (qPCR)**
Certified RNase, DNase, DNA
& Pyrogen free

Art. No. E2796-9795
100 self-adhesive seals
OFFER: 93 €



**Polypropylene Sealing
Film Strips (PCR)**
Certified RNase, DNase &
Nucleid Acid free

Art. No. E2796-2850
200 Strips
OFFER: 27 €




Polyester Sealing Film (PCR)
Certified RNase, DNase & Nucleid
Acid free

Art. No. E2796-0100
100 self-adhesive seals
OFFER: 38 €



Aluminium Sealing Film (PCR)
Certified RNase & DNase free

Art. No. E2796-1100
100 self-adhesive seals
OFFER: 59 €



Microplate Mixer.

MIX IT. SHAKE IT.

Starlab's Microplate Mixer has a stepless, variable vibration speed and is suitable for use with all skirted plates with a SBS footprint. The small footprint saves space on your bench too.

FEATURES

- 0 - 1500 rpm
- Takes any skirted plate with SBS footprint
- Stepless speed control

USE ME FOR:

- Mixing and shaking any skirted plate
- Resuspending pellets

Microplate Mixer



Art. No. N2400-8040

OFFER: 265 €



BUY ONLINE

starlab.click/microplatemixer



96 Square Wells, Conical Bottoms

Polypropylene deepwell plates with moulded matrix. Certified RNase, DNase & Pyrogen free. Tested endotoxin free.

1.2 ml Deepwell

10 plates

Art. No. E2896-1820

OFFER: 34 €

2.2 ml Deepwell

10 plates

Art. No. E2896-1810

OFFER: 54 €

PlateOne®

Starlab's PlateOne® Microplates and Deepwell plates with SBS footprint are all suitable for use with the Microplate Mixer.

PlateOne® plates are made from high-quality virgin polypropylene, the plates feature a highly-visible printed matrix and are certified RNase, DNase, DNA and Pyrogen free. Sterile options also available.

PlateOne® Plates.

PlateOne® MICROPLATES

96 Round Wells, Round Bottoms (320 µl)



Non-sterile	4 x 20 plates	S1830-9600	111 €
Sterile	4 x 20 plates	S1830-9610	129 €

96 Round Wells, Conical Bottoms (320 µl)



Non-sterile	4 x 20 plates	S1833-9600	111 €
Sterile	4 x 20 plates	S1833-9610	129 €

96 Round Wells, Flat Bottoms (370 µl)



Non-sterile	4 x 20 plates	S1837-9600	111 €
Sterile	4 x 20 plates	S1837-9610	129 €

PlateOne® DEEPWELL PLATES

96 Round Wells, Round Bottoms (0.5 ml)



Non-sterile (Bulk)	4 x 8 plates	S1896-5800	105 €
Non-sterile	32 plates	S1896-5000	94 €
Sterile (Ind. Wrap)	32 plates	S1896-5110	121 €

96 Round Wells, Round Bottoms (1.0 ml)



Non-sterile (Bulk)	4 x 5 plates	S1896-1800	66 €
Non-sterile	20 plates	S1896-1000	63 €
Sterile (Ind. Wrap)	20 plates	S1896-1110	76 €

96 Square Wells, V-Shaped Bottoms (2.0 ml)



Non-sterile (Bulk)	4 x 5 plates	S1896-2800	75 €
Non-sterile	20 plates	S1896-2000	68 €
Sterile (Ind. Wrap)	20 plates	S1896-2110	87 €

384 Square Wells, V-Shaped Bottoms (200 µl)



Non-sterile (Bulk)	4 x 5 plates	S1884-2400	157 €
Non-sterile	20 plates	S1884-2000	141 €
Sterile (Ind. Wrap)	20 plates	S1884-2410	181 €



BUY ONLINE

starlab.click/plateone

Last but not least ...

MORE EQUIPMENT FOR YOUR DAILY TASKS...



Mini Dry Bath Incubator

Art. No. N2400-4020

A very small yet powerful unit suitable for a range of applications. The digital microprocessor controls the preset temperature of the metal block. **Metal blocks sold separately**, see website for ordering info.

OFFER: 155 €



BUY ONLINE

starlab.click/mini-dry-bath



Magnetic Stirrer

Art. No. N2400-2010

An economical magnetic stirrer with a variable speed range of 100 – 1500 rpm and a stirring volume of up to 3 litres. Supplied with a 40 mm stir bar. Easy to clean and small footprint.

OFFER: 120 €



BUY ONLINE

starlab.click/stirrers



Magnetic Stirrer

WITH HEATING Art. No. N2400-3010



Fast temperature ramping saves you time! This stirrer features a speed range of 200 – 1500 rpm and a temperature range from ambient to 280 °C (in 1 °C increments). Supplied with a 40 mm stir bar.

OFFER: 235 €



BUY ONLINE

starlab.click/stirrers



Rotator Mixer Multi 1

Art. No. N2400-5000

A compact, multifunctional unit that provides thorough mixing of microtubes with consistent results. Features three different motions (vertical, reciprocal, vibration) that can be used separately or combined.

OFFER: 495 €



BUY ONLINE

starlab.click/rotatormixer

Sample storage/ cooling solutions.




StarChill PCR, Purple/Pink

Provides you with an instant, visual indication of sample temperature. Temperature colour change at 7° C. 96-well rack with lid for use with 0.2 ml tubes or strips, or a 0.2 ml 96-well PCR plate.

Art. No. E2396-2004

Pack of 2

OFFER: 57 €



StarStore 100

100-place storage box with fixed, 10x10 dividers for tubes up to 13 mm, eg. 1.5/2.0 ml microcentrifuge or screw cap tubes. Autoclavable. Temperature use: -90°C to 121°C. WxDxH: 140 x 155 x 55 mm.

Art. No. I2310-5848

Pack of 5, Mixed Neon

OFFER: 25 €



Polycarbonate Cryogenic Storage Boxes

White base with transparent, lift-off lid with printed reference grid. Temperature use from -196°C to 121°C. Standard 130x130 footprint.

Art. No. I2381-2305

81-place

OFFER: 36 € each

Art. No. I2510-5005

100-place

OFFER: 35 € each

Starlab ErgoOne® manual pipettes.

ErgoOne®

Single and multi-channel
manual pipettes

from 120 €




ORDER ONLINE

starlab.click/manualpipettes

**COMFORTABLE,
RELIABLE PIPETTING:
ALL DAY, EVERY DAY**

- Ergonomic and lightweight
- Volume lock
- Universal tip fit



Single-Channel Pipettes

S7100-0125 0,1 - 2,5 µl OFFER: 120 €	S7100-0510 0,5 - 10 µl OFFER: 120 €	S7100-0221 2 - 20 µl Micro OFFER: 120 €	S7100-0220 2 - 20 µl Standard OFFER: 120 €
S7100-1100 10 - 100 µl OFFER: 120 €	S7100-2200 20 - 200 µl OFFER: 120 €	S7110-1000 100 - 1000 µl OFFER: 120 €	S7150-5000 500 - 5000 µl OFFER: 120 €

8-Channel Pipettes

S7108-0510 0,5 - 10 µl OFFER: 415 €	S7108-1100 10 - 100 µl OFFER: 415 €	S7108-3300 30 - 300 µl OFFER: 415 €
--	--	--

12-Channel Pipettes

S7112-0510 0,5 - 10 µl OFFER: 510 €	S7112-1100 10 - 100 µl OFFER: 510 €	S7112-3300 30 - 300 µl OFFER: 510 €
--	--	--

Starlab ErgoOne® FAST pipette controller.

ERGONOMY MEETS MODERNITY

ErgoOne® FAST is one of the lightest pipette controllers on the market. It is designed for comfort and balances perfectly in the hand for long, fatigue-free pipetting.

Use with all plastic or glass serological pipettes (1.0 to 100 ml). The ErgoOne® FAST can fill a 25 ml pipette in less than four seconds!



ErgoOne® Fast Pipette Controller

Art. No. S7166-0010

OFFER: 245 €

ORDER ONLINE

starlab.click/pipettecontroller



Choose your personal protection. StarGuard® Gloves.


COMFORT

for enhanced nitrile
comfort and dexterity



RECOMMENDED FOR USE WHEN:

- Work demands dexterity and tactile sensitivity
- You require a lightweight glove with exceptional fit
- You want to reduce hand fatigue caused by gloves
- You need to combine comfort and chemical protection
- Testing or handling food stuffs



Description	Pack Size	Art. No.	OFFER
StarGuard® Comfort, Nitrile Gloves, XS	10 × 250	SG-C-XS	245 €
StarGuard® Comfort, Nitrile Gloves, S	10 × 250	SG-C-S	245 €
StarGuard® Comfort, Nitrile Gloves, M	10 × 250	SG-C-M	245 €
StarGuard® Comfort, Nitrile Gloves, L	10 × 250	SG-C-L	245 €
StarGuard® Comfort, Nitrile Gloves, XL ¹⁾	10 × 230	SG-C-XL	245 €

¹⁾ 230 gloves per box XL



AQL
0.65



PROTECT

for nitrile strength and durability



RECOMMENDED FOR USE WHEN:

- Combines strength and durability for everyday protection
- Comfort, tactile sensitivity and consistent grip
- Excellent chemical splash protection
- Longer glove length (250 mm) covers the wrist

Description	Pack Size	Art. No.	OFFER
StarGuard® Protect, Nitrile Gloves, XS	10 × 100	SG-P-XS	155 €
StarGuard® Protect, Nitrile Gloves, S	10 × 100	SG-P-S	155 €
StarGuard® Protect, Nitrile Gloves, M	10 × 100	SG-P-M	155 €
StarGuard® Protect, Nitrile Gloves, L	10 × 100	SG-P-L	155 €
StarGuard® Protect, Nitrile Gloves, XL	10 × 100	SG-P-XL	155 €



BUY ONLINE

starlab.click/starguard

SENSITIVE

comfortable nitrile protection that is clinically proven to reduce the potential of contact dermatitis



RECOMMENDED FOR USE WHEN:

- Low Dermatitis Potential - ideal for users with sensitive skin
- Contains no contaminants to interfere with experiments
- Improved tactile sensitivity
- Excellent grip in wet conditions

Description	Pack Size	Art. No.	OFFER
StarGuard® Sensitive, Nitrile Gloves, XS	10 × 200	SG-N-XS	195 €
StarGuard® Sensitive, Nitrile Gloves, S	10 × 200	SG-N-S	195 €
StarGuard® Sensitive, Nitrile Gloves, M	10 × 200	SG-N-M	195 €
StarGuard® Sensitive, Nitrile Gloves, L	10 × 200	SG-N-L	195 €
StarGuard® Sensitive, Nitrile Gloves, XL	10 × 200	SG-N-XL	195 €



Starlab International GmbH

Neuer Höltigbaum 38
22143 Hamburg
Germany

T: +49 (0)40 675 99 39 0
F: +49 (0)40 675 99 39 20
info@starlab.de
www.starlabgroup.com

Starlab GmbH

Neuer Höltigbaum 38
22143 Hamburg
Germany

T: +49 (0)40 675 99 39 0
F: +49 (0)40 675 99 39 20
info@starlab.de
www.starlab.de

Starlab (UK), Ltd

5 Tanners Drive
Milton Keynes MK14 5BU
United Kingdom

T: +44 (0)1908 283800
F: +44 (0)1908 283802
infoline@starlab.co.uk
www.starlab.co.uk

Starlab FRANCE SARL

30 Rue Jean Rostand
91400 Orsay
France

T: +33 (0)1 60 13 71 70
F: +33 (0)1 69 41 48 65
info@starlab.fr
www.starlab.fr

Starlab S.r.l.

Via Pinturicchio 1
20133 Milano
Italy

T: +39 (0)2 7020 1040
F: +39 (0)2 7020 1033
info@starlab.it
www.starlab.it



STARGUARD® GLOVES

Back on stock: our range of StarGuard® nitrile gloves offer daily protection in the laboratory without sacrificing comfort.

For more information on our StarGuard® gloves and details on Chemical Breakthrough Times, please see our website.



VISIT ONLINE SHOP

starlab.click/starguard

