

This is the sound of the Starlab special offers.

valid to 15 December 2022

IN THIS ISSUE

Equipment & Pipettes

Starlab is best known for high-quality TipOne® pipette tips, but we have a whole orchestra of equipment to cover your daily needs in the lab. Mixers, shakers, vortexers: you name it, we've got you covered.

Explore the complete symphony of lab equipment we have on offer, and you'll be lab-dancing along in no time.

www.starlab.click/soundofstarlab

Passionate for science.



Listen to the sound of Starlab:
We have created a cool bit of music -
generated solely with the sounds from
Starlab equipment.

 starlab


Contents of the Starlab Special Offers.



Mixer HC
page 4




GuardOne® Workstation
page 5



Tube Roller
page 7



Plate Centrifuge
page 7



Microplate Mixer
page 6



Dry Baths
page 8



Vortexers
page 10




Rotator Mixer
page 6



Magnetic Stirrers
page 9

PIPETTES



ErgoOne® Pipettes
page 13



ErgoOne® FAST Controller
page 12




ErgoOne®E Pipettes
page 14

GLOVES



Nitrile and Latex Gloves
page 15



NEW: THE STARLAB MINICENTRIFUGE



SMOOTH OPERATOR

With its distinctive curved design that will look fantastic on your bench, the Starlab Minicentrifuge is not all about style, it also has substance! This is a compact and solid unit that delivers a quiet, smooth operation and unique features.

page 11

Shake it. Heat it.

The Mixer HC.

MIXING, HEATING, COOLING & TIMING
Combined in one unit to help speed up your daily tasks

EASY TO USE
Digital display for precise adjustment. Clear menu navigation and simple parameter setting

PRECISE TEMPERATURE CONTROL
Up to 99°C or down to 13°C below room temperature

MIXING
Up to 1500 rpm

TIMING
1 min to 99:59 hours

SMALL, COMPACT UNIT
You will not believe how small it is!

FUNCTIONAL
INTERVAL MIX: alternate between mixing & resting
PROGRAM MODE: create your own protocols
SHORT MIX: instant mixing at the push of a button

AUTOMATIC THERMOBLOCK RECOGNITION
Provides quick and easy changing of thermoblocks between protocols. No programming required!



Mixer HC

The Mixer HC will fit perfectly into any application where you need to heat, cool or mix, or perform a combination of the three. The tiny footprint fits on any bench to provide a convenient, all-in-one solution that will simplify your workflows.

Art. No. S8012-0000

OFFER 1.550 €*

*Offer price, additional discounts not applicable. Valid only until 15.12.2022. Thermoblock not included.



VISIT ONLINE SHOP

starlab.click/mixerhc



NINE THERMOBLOCKS AVAILABLE

Versatile! Use the Mixer HC with all popular tubes and plates.

All-in for your applications. The GuardOne® Workstation.

GuardOne® Workstation

Starlab's GuardOne® Workstation provides the ideal workspace for sensitive and critical applications requiring a clean environment or protection from potential contamination by bacteria, DNA or RNA.


From 2.450 €*

*Offer price, additional discounts not applicable. Valid only until 15.12.2022.



ORDER ONLINE

starlab.click/guardone



AVAILABLE MODELS



Laminar Flow Workstation 32"
Art. No. S8040-0000
OFFER: 2.850 €



Laminar Flow Workstation 48"
Art. No. S8040-0001
OFFER: 3.750 €



Dead Air Workstation 32"
Art. No. S8040-1000
OFFER: 2.450 €



Dead Air Workstation 48"
Art. No. S8040-1001
OFFER: 3.050 €

Rotator Mixer 1.

MIX IT. SHAKE IT.

A compact, multifunctional unit that provides thorough mixing of microtubes with consistent results. Features three different motions (vertical, reciprocal, vibration) that can be used separately or combined.

FEATURES

- Controlled sample inverting
- Programmable
- Three different motions

Art. No. N2400-5000

OFFER: 525 €*




BUY ONLINE

starlab.click/rotatormixer

Microplate Mixer.

MIX IT. SHAKE IT.

Starlab's Microplate Mixer has a stepless, variable vibration speed and is suitable for use with all skirted plates with a SBS footprint. The small footprint saves space on your bench too.



FEATURES

- 0 - 1500 rpm
- Takes any skirted plate with SBS footprint
- Stepless speed control

USE ME FOR:

- Mixing and shaking any skirted plate
- Resuspending pellets

Microplate Mixer

Art. No. N2400-8040

OFFER: 287 €*

*Offer price, additional discounts not applicable.
Valid only until 15.12.2022.



BUY ONLINE

starlab.click/microplatemixer

Plate Centrifuge.

YOU SPIN ME ROUND!

Sealed plates are loaded vertically and surface tension keeps samples in place so there is no danger of spilling samples. The Plate Centrifuge accepts all plates up to 25 mm high and is less than quarter the size of most plate centrifuges.



FEATURES

- 2500 rpm / 500 x g
- Takes any plate up to 25 mm high
- Rapid start to stop

USE ME FOR:

- Safe and quick spin down of droplets and condensation
- Spinning before and after thermal cycling to increase PCR yield.

Art. No. N2631-0008

OFFER: 575 €*



BUY ONLINE

starlab.click/platecentrifuge

USE ME WITH 1/2 PLATES, STRIPS & TUBES TOO

An adapter plate (N2631-0018) is available separately which allows you to use the plate centrifuge with 48-well plates or PCR tubes or strips.

Tube Roller.

IT'S A ROCK AND ROLL STAR

Simultaneous rolling and tilting movements ensures thorough mixing of your samples. The unit features variable speed control and a sturdy and robust construction. The tube roller is built for continuous operation and is suitable for applications such as hybridisation.

- Simultaneous rolling and tilting motion
- Continuous operation
- Speed control up to 70 rpm
- Six rollers holds load up to 4 kg
- Easy to use

Art. No. N2400-7010


OFFER: 355 €*

*Offer price, additional discounts not applicable.
Valid only until 15.12.2022.



BUY ONLINE

starlab.click/tuberoller



Mini Dry Bath & Dry Bath System.

TINY, BUT POWERFUL

A very small yet powerful unit suitable for a range of applications. The digital microprocessor controls the preset temperature of the metal block. Metal blocks sold separately, see website for ordering info.



- Microprocessor controlled temperature and timer
- Digital display
- Very small footprint

MINI DRY BATH INCUBATOR

ART. NO. N2400-4020

OFFER: 185 €*



BUY ONLINE

starlab.click/mini-dry-bath

DRY BATH. WET BATH.

Available as single or dual block units, these systems offer great value and convenience for the user. The range of blocks make these units suitable for a wide variety of applications. Both systems can also function as a water bath.



- Microprocessor controlled temperature and timer
- Digital display
- Wide range of metal blocks

DRY BATH INCUBATOR

ART. NO. N2400-4001
(SINGLE BLOCK)

OFFER: 335 €*

ART. NO. N2400-4002
(DUAL BLOCK)

OFFER: 425 €*

*Offer price, additional discounts not applicable. Valid only until 15.12.2022.



BUY ONLINE

starlab.click/drybaths



Stir it up!

MAGNETIC STIRRERS FOR YOUR DAILY TASKS...



MAGNETIC STIRRER
ART. NO. N2400-2010

An economical magnetic stirrer with a variable speed range of 100 – 1500 rpm and a stirring volume of up to 3 litres. Supplied with a 40 mm stir bar. Easy to clean and small footprint.

OFFER: 130 €*



BUY ONLINE

starlab.click/stirrers



MAGNETIC STIRRER
WITH HEATING
ART. NO. N2400-3010

Fast temperature ramping saves you time! This stirrer features a speed range of 200 – 1500 rpm and a temperature range from ambient to 280°C (in 1°C increments). Supplied with a 40 mm stir bar.

OFFER: 250 €*

*Offer price, additional discounts not applicable.
Valid only until 15.12.2022.



BUY ONLINE

starlab.click/stirrers

Vortexers

SIMPLE AND SMART!

Choose between the economical Vortex with speeds up to 2500 rpm, or the Vortex with optional IR motion detector and speeds up to 3000 rpm. Which ever you choose, you know either of these compact units will fit nicely in your GuardOne®.



- Variable speed control up to 2500 rpm
- Continuous and touch modes
- Wide range of accessories
- 2-year warranty

ART. NO. N2400-6010

OFFER: 145 €*



BUY ONLINE

starlab.click/vortex



- Variable speed control up to 3000 rpm
- Continuous and sensor modes
- Sensor automatically starts vibration
- Wide range of accessories
- 2-year warranty


ART. NO. N2400-0001

OFFER: 205 €*



BUY ONLINE

starlab.click/vortexir



New: The Minicentrifuge.

AHEAD OF THE CURVE!

With its distinctive curved design that will look fantastic on your bench, the Starlab Minicentrifuge is not all about style, it also has substance! This is a compact and solid unit that delivers a quiet, smooth operation and unique features.

- Quick start on lid closure
- Magnetic lid lock
- Smooth, quiet operation
- One-touch STOP for fast braking and automatic lid release

ART. NO. S8030-5000

PRICE: 245 €

MAGNETIC LID LOCK
aids easy,
one-handed
operation

QUICK START
just close the lid
to start spinning!

TOOL FREE
quick and easy
rotor exchange

ONE-TOUCH STOP
for fast
braking and
lid release

OH SO QUIET
noise level
55 dBA



Starlab ErgoOne® FAST pipette controller.

ERGONOMY MEETS MODERNITY

ErgoOne® FAST is one of the lightest pipette controllers on the market. It is designed for comfort and balances perfectly in the hand for long, fatigue-free pipetting.

Use with all plastic or glass serological pipettes (1.0 to 100 ml). The ErgoOne® FAST can fill a 25 ml pipette in less than four seconds!



ErgoOne® Fast Pipette Controller

Art. No. S7166-0010

OFFER: 245 €*

*Offer price, additional discounts not applicable.
Valid only until 15.12.2022.



ORDER ONLINE

starlab.click/pipettecontroller

Starlab ErgoOne® manual pipettes.

STYLE AND COMFORT: YOU CAN HAVE BOTH!

Reliability is key in the laboratory and your pipette must be trusted to perform well, day in, day out. ErgoOne® manual pipettes give you the accuracy and precision you can trust, with a stylish new look you will love.

ErgoOne®

Single and multi-channel manual pipettes

from 125 €*

*Offer price, additional discounts not applicable.
Valid only until 15.12.2022.



ORDER ONLINE

starlab.click/manualpipettes



Single-Channel Pipettes

S7100-0125 0,1 - 2,5 µl OFFER: 125 €	S7100-0510 0,5 - 10 µl OFFER: 125 €	S7100-0221 2 - 20 µl Micro OFFER: 125 €	S7100-0220 2 - 20 µl Standard OFFER: 125 €
S7100-1100 10 - 100 µl OFFER: 125 €	S7100-2200 20 - 200 µl OFFER: 125 €	S7110-1000 100 - 1.000 µl OFFER: 125 €	S7150-5000 500 - 5.000 µl OFFER: 125 €

8-Channel Pipettes

S7108-0510 0,5 - 10 µl OFFER: 425 €	S7108-1100 10 - 100 µl OFFER: 425 €	S7108-3300 30 - 300 µl OFFER: 425 €
--	--	--

12-Channel Pipettes

S7112-0510 0,5 - 10 µl OFFER: 515 €	S7112-1100 10 - 100 µl OFFER: 515 €	S7112-3300 30 - 300 µl OFFER: 515 €
--	--	--

- Ergonomic and lightweight
- Volume lock
- Universal tip fit

ErgoOne® E electronic pipettes.

RELIABLE & ERGONOMIC

ErgoOne® E single and multi-channel pipettes provide reliable pipetting every day. These micro-processor controlled, battery-operated, piston-stroke members of the ErgoOne® family have been optimised for ergonomic operation and ease of use.

Five operating modes: standard pipetting, reverse pipetting, mixing and dispensing, plus a unique Electrophoresis mode to help you load gels easily.

Single-Channel Pipettes

G9001-0010 0,5 - 10 µl OFFER: 295 €	G9001-0020 2,0 - 20 µl OFFER: 295 €	G9001-0200 10 - 200 µl OFFER: 295 €	G9001-1000 50 - 1.000 µl OFFER: 295 €
G9001-5000 250 - 5.000 µl OFFER: 295 €			

8-Channel Pipettes

G9008-0010 0,5 - 10 µl OFFER: 545 €	G9008-0020 1,0 - 20 µl OFFER: 545 €	G9008-0100 5,0 - 100 µl OFFER: 545 €	G9008-0300 15 - 300 µl OFFER: 545 €
---	---	--	---

12-Channel Pipettes

G9012-0010 0,5 - 10 µl OFFER: 655 €	G9012-0020 1,0 - 20 µl OFFER: 655 €	G9012-0100 5,0 - 100 µl OFFER: 655 €	G9012-0300 15 - 300 µl OFFER: 655 €
---	---	--	---

ErgoOne® E Single and multi-channel electronic pipettes

from **295 €***

*Offer price, additional discounts not applicable.
Valid only until 15.12.2022.



ORDER ONLINE

starlab.click/electronicpipettes



Choose your personal protection.

COMFORT

for enhanced nitrile comfort and dexterity



Description	Pack Size	Art. No.	OFFER
StarGuard® Comfort, Nitrile Gloves, XS	10 × 250	SG-C-XS	199 €*
StarGuard® Comfort, Nitrile Gloves, S	10 × 250	SG-C-S	199 €*
StarGuard® Comfort, Nitrile Gloves, M	10 × 250	SG-C-M	199 €*
StarGuard® Comfort, Nitrile Gloves, L	10 × 250	SG-C-L	199 €*
StarGuard® Comfort, Nitrile Gloves, XL ¹⁾	10 × 230	SG-C-XL	199 €*

¹⁾ 230 gloves per box XL

SENSITIVE

comfortable nitrile protection that is clinically proven to reduce the potential of contact dermatitis



Description	Pack Size	Art. No.	OFFER
StarGuard® Sensitive, Nitrile Gloves, XS	10 × 200	SG-N-XS	175 €*
StarGuard® Sensitive, Nitrile Gloves, S	10 × 200	SG-N-S	175 €*
StarGuard® Sensitive, Nitrile Gloves, M	10 × 200	SG-N-M	175 €*
StarGuard® Sensitive, Nitrile Gloves, L	10 × 200	SG-N-L	175 €*
StarGuard® Sensitive, Nitrile Gloves, XL	10 × 200	SG-N-XL	175 €*

PROTECT

for nitrile strength and durability



Description	Pack Size	Art. No.	OFFER
StarGuard® Protect, Nitrile Gloves, XS	10 × 100	SG-P-XS	125 €*
StarGuard® Protect, Nitrile Gloves, S	10 × 100	SG-P-S	125 €*
StarGuard® Protect, Nitrile Gloves, M	10 × 100	SG-P-M	125 €*
StarGuard® Protect, Nitrile Gloves, L	10 × 100	SG-P-L	125 €*
StarGuard® Protect, Nitrile Gloves, XL	10 × 100	SG-P-XL	125 €*

TOUCH

if you prefer the secure grip and tactile sensitivity of latex



Description	Pack Size	Art. No.	OFFER
StarGuard® Touch, Latex Gloves, XS	10 × 100	SG-T-XS	75 €*
StarGuard® Touch, Latex Gloves, S	10 × 100	SG-T-S	75 €*
StarGuard® Touch, Latex Gloves, M	10 × 100	SG-T-M	75 €*
StarGuard® Touch, Latex Gloves, L	10 × 100	SG-T-L	75 €*
StarGuard® Touch, Latex Gloves, XL	10 × 90	SG-T-XL	75 €*



BUY ONLINE

starlab.click/starguard



Starlab International GmbH

Neuer Höltigbaum 38
22143 Hamburg
Germany

T: +49 (0)40 675 99 39 0
F: +49 (0)40 675 99 39 20
info@starlab.de
www.starlabgroup.com

Starlab GmbH

Neuer Höltigbaum 38
22143 Hamburg
Germany

T: +49 (0)40 675 99 39 0
F: +49 (0)40 675 99 39 20
info@starlab.de
www.starlab.de

Starlab (UK), Ltd

5 Tanners Drive
Milton Keynes MK14 5BU
United Kingdom

T: +44 (0)1908 283800
F: +44 (0)1908 283802
infoline@starlab.co.uk
www.starlab.co.uk

Starlab FRANCE SARL

30 Rue Jean Rostand
91400 Orsay
France

T: +33 (0)1 60 13 71 70
F: +33 (0)1 69 41 48 65
info@starlab.fr
www.starlab.fr

Starlab S.r.l.

Via Bracco 6
20159 Milano
Italy

T: +39 (0)2 7020 1040
F: +39 (0)2 7020 1033
info@starlab.it
www.starlab.it

STARLAB SPECIAL OFFERS

Find the full range of our products on offer on our website. If you haven't registered already, set up an account and order online where the discount is applied automatically.



VISIT ONLINE SHOP


starlab.click/specialoffers



Listen to the Sound of Starlab!

WATCH VIDEO

starlab.click/soundofstarlab



Scan QR code for more information & ordering