


**Stock up now because
your research can't wait.
You are passionate
for science.**

SPECIAL OFFERS

Pipettes & Lab Equipment

Valid until 13.12.2024.



Passionate for science.



GuardOne® Workstation.



from: £2,985*

DESIGNED WITH THE ENVIRONMENT IN MIND

Long-life consumables reduces waste. Minimal packaging volume reduces transportation needs. Minimal power consumption.

AVOID CROSS-CONTAMINATION & THE RISK OF DUPLICATING WORK

Protect your valuable samples and reagents¹⁾ to help prevent both the extra time, cost and waste involved when having to duplicate contaminated work. Starlab's GuardOne® Workstation provides the ideal clean and controlled environment. GuardOne® is perfect for your PCR / qPCR applications or any other application that requires a controlled, clean environment.

GuardOne® has been designed with convenience in mind. All functions are controlled via the glove compatible, easy to use touchscreen. Tool-free HEPA filter exchange and UV lamp exchange in just 2 minutes!

- One-touch QuickStart / QuickStop settings
- Efficient UV_{245 nm} decontamination with timer function
- Efficient, class H14 HEPA filter with pre-filter provides ISO class 5 working environment²⁾
- Real-time airflow measurement²⁾
- Life-span monitoring of consumables
- Virtually maintenance-free
- Easy to clean
- Easy, tool-free 15-minute assembly



Dead Air Workstation 32"	Dead Air Workstation 48"	Laminar Flow Workstation 32"	Laminar Flow Workstation 48"
S8040-1000	S8040-1001	S8040-0000	S8040-0001
£3,553.45	£4,560.39	£4,539.14	£5,741.06
£2,665*	£3,420*	£3,488*	£4,480*

starlab.click/guardone

¹⁾GuardOne® is not designed to protect the user.
²⁾Laminar flow models only.



***OFFER REF: Q24000671068**
must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

Starlab Vortex Mixer.



WATCH VIDEO

GREAT FEATURES. GOOD VIBRATIONS!

The new Starlab Vortex Mixer is ideal for the gentle mixing of samples right up to vigorous vortex mixing and the resuspending of cells or chemical pellets.

- Unique, glove compatible touch screen
- Precise vortexing speeds (300 to 3,500 rpm)
- One-touch pre-set options for speed / time
- Touch or continuous operation
- End of run audible alarm
- Power-saving sleep mode
- Optional adapters available separately to mix up to 15 microcentrifuge tubes or one microplate (up to 1,500 rpm)

TOUCH MODE
Minimal pressure required to maintain vibration. Mix single or multiple tubes up to 3,500 rpm.

GLOVE COMPATIBLE TOUCHSCREEN
for precise speed and time settings

SMOOTH, QUIET OPERATION



POWER-SAVING SLEEP MODE

SPEED RANGE
300 – 3,500 RPM

ANTI-SLIP FEET
provides stability on the bench during mixing

Vortex Mixer.

Art. No. S8010-0000

£350

OFFER: £306*

 starlab.click/shop-vortex



***OFFER REF: Q24000671068**
must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

The Starlab Aspirator.

USE TWO HAND CONTROLLERS WITHOUT ANY LOSS OF PERFORMANCE

- Powerful pump that can operate two hand controllers without any loss of performance
- Auto ON-OFF pump switch
- Low noise emission for a quieter working environment
- Supplied with everything you need to get started
- Versatile! A whole range of accessories also available
- The perfect partner for use with Starlab's bottle top filtration system

N2400-9000 includes everything you need to get started: Aspirator, 4 Litre Polypropylene Bottle with Protection Filter and Inlet Pipe, Hand Controller, including 2.5 m Tubing, a Wall Mount, an Adapter for Pipettes, and an Adapter for Pipette Tips.



DESIGNED FOR THE SAFE AND EFFICIENT ASPIRATION OF FLUIDS

- Vacuum level and suction power is adjustable in 8 steps (150 to 850 mbar) for safe and reproducible working
- The chemically-resistant pump automatically switches on/off to reduce noise level in the workplace
- You can also use the Aspirator as a mobile vacuum source for sterile filtration and solid phase extraction
- Additional hand controllers and a range of other accessories are available. See website for details

Aspirator

Art. No. N2400-9000

£2154.98

OFFER: £1,341*

 starlab.click/aspirator

Discover our range of Aspiration Pipettes



starlab.click/aspirationpipettes

FULLY AUTOCLAVABLE
ergonomically designed hand controller



***OFFER REF: Q24000671068**
must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.


The Minicentrifuge.

AHEAD OF THE CURVE!

Ideal for all your quick spin-downs, the Starlab Minicentrifuge is not all about style, it also has substance! This is a compact and solid unit that delivers a quiet, smooth operation and a maximum speed of 6,500 rpm.

- Quick start on lid closure
- Magnetic lid lock
- One-touch STOP for fast braking and automatic lid release

SUPPLIED WITH TWO ROTORS



Tube Rotor holds up to eight microcentrifuge tubes or eight conical or skirted screw cap microtubes



PCR Strip Rotor holds two, 0.2 ml 8-strips or 32 x 0.2 ml single tubes.



starlab.click/minicentrifuge



QUICK START
just close the lid
to start spinning!

**QUICK & EASY
TOOL FREE
ROTOR EXCHANGE**

**ONE-TOUCH
STOP** for fast
braking and
lid release

**MAGNETIC
LID LOCK**
aids easy,
one-handed
operation



WATCH VIDEO

Minicentrifuge.

Art. No. S8030-5000

~~£233.97~~

OFFER: £190*



***OFFER REF: Q24000671068**
must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

The versatile Mixer HC.

 starlab.click/mixer-hc

ONE SMALL DEVICE, SO MANY USES.

The Mixer HC is a smart bit of kit that does the work of several other pieces of equipment in the lab. Easy to use and to program, the Mixer HC will become an invaluable benchtop companion to help speed up your laboratory tasks. You will wonder how you ever managed without one!

MIXING
up 1,500 rpm

HEATING
up 99 °C

COOLING
to 13 °C
below room
temperature



TIMING
1 min to
99:59 hours

FUNCTIONAL
Alternate between
mixing and resting.
Create your own
protocols.
Or instant mixing
at the push of a
button



WATCH VIDEO

Mixer HC

Art. No. **S8012-0000** (Excludes Thermoblock)

£2,335.38

OFFER: £1,635*

NINE EXCHANGEABLE THERMOBLOCKS

Use the Mixer HC with all popular tubes and plates.

			
24 x 0.5 ml Microcentrifuge Tubes S8012-0010 ¹⁾ £449.93 £376*	24 x 1.5 ml Microcentrifuge Tubes S8012-0011 ¹⁾ £449.93 £376*	24 x 2.0 ml Microcentrifuge Tubes, Round Bottoms S8012-0012 ¹⁾ £449.93 £376*	24 x 2.0 ml Microcentrifuge Tubes Ø 12 mm, V-Bottoms S8012-0013 £620.88 £466*
			
8 x 5.0 ml Preparation Tubes S8012-0014 £449.93 £376*	8 x 15 ml Conical Centrifuge Tubes S8012-0015 £449.93 £376*	4 x 50 ml Conical Centrifuge Tubes S8012-0016 £449.93 £376*	24 x 1.5/2.0 ml Cryo Tubes S8012-0017 £449.93 £376*
			
Micro/Deepwell Plate, with Lid S8012-0018 £620.88 £466*	Adapter Plate 96 x 0.2 ml Tubes or one PCR plate (Use with S8012-0018) S8012-0019 £169.35 £128*	ISO-Rack: 24 x 0.5 ml Microcentrifuge Tubes S8012-0020 £69.87 £54*	ISO-Rack: 24 x 1.5/2.0 ml Microcentrifuge Tubes S8012-0021 £69.87 £54*


¹⁾ Includes ISO Rack. Tubes not included.

***OFFER REF: Q24000671068**
must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

Shaking & Mixing.

ESSENTIAL EQUIPMENT FOR DAILY TASKS

Orbital Shaker.



- Speed range: 100–500 rpm
- LCD displays for speed & time
- Continuous or timed mixing
- Supplied without platform

Orbital Shaker
Art. No: N2400-8030
 (excludes platform)
~~£1,383.82~~
OFFER: £830*



Platform with clamping rolls
 (for use with Orbital, 3D and Rocking Shakers)


Art. No: N2400-8011
~~£382.00~~
OFFER: £306*



Platform with rubber mat and bungee cords
 (for use with Orbital, 3D and Rocking Shakers)

Art. No: N2400-8021
~~£382.00~~
OFFER: £306*

3D / Rocking Shakers.



- See saw rocking or 3D mixing
- Speed range: 10 –70 rpm
- LCD displays for speed & time
- Continuous or timed mixing.
- Supplied with platform with bungee cords

Rocking Shaker
Art. No: N2400-8010
~~£1,577.12~~
OFFER: £945*

3D Rocking Shaker
Art. No: N2400-8020
~~£1,577.12~~
OFFER: £945*



Rotator Mixer.

Art. No. N2400-5000
~~£886.75~~
OFFER: £621*

- Controlled sample inverting
- Speed range: 1 – 100 rpm
- Programmable
- Three different motions
- Variable speed and timer
- Low profile and small footprint



Vertical Rotation
 360° vertical motion for time period of 0–250 seconds, or non-stop.



Reciprocal Motion
 Platform direction changes in turns (turning angle 1–90°) for 0–250 seconds, or non-stop.



Vibration
 Same as reciprocal motion, plus additional vibration (turning angle 1–5° for 1–5 seconds.



Pause
 Platform direction changes in turns, incl. a temporary motion stop of 1–5 seconds (in conjunction with reciprocal motion only).

starlab.click/shakers

starlab.click/rotatormixer



***OFFER REF: Q24000671068**
 must be stated on your PO (excl. eShop).
 Offer valid to 13.12.2024.

Mixing & Stirring.

ESSENTIAL EQUIPMENT FOR DAILY TASKS



Use for mixing and shaking any skirted plate with a SBS footprint, and for resuspending pellets.

- Speed range: 0 – 1500 rpm
- Stepless speed control
- Range of accessories available

Microplate Mixer.

Art. No. N2400-8040

£571.16

OFFER: £343*



starlab.click/microplate-mixer



An economical magnetic stirrer.

- Speed range: 100 – 1500 rpm
- Stirring volume up to 3 litres
- Supplied with 1 x 40 mm stir bar

Magnetic Stirrer.

Art. No. N2400-2010

£244.11

OFFER: £146*



starlab.click/stirrers



Fast temperature ramping saves you time!

- Speed range: 200 – 1500 rpm.
- Temperature range:
ambient to 280 °C (1°C increments)
- Stirring volume up to 3 litres
- Supplied with 1 x 40 mm stir bar

Magnetic Stirrer with Heating.

Art. No. N2400-3010

£430.32

OFFER: £280*



starlab.click/stirrers



***OFFER REF: Q24000671068**

must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

Rolling.



Built for continuous operation. Suitable for applications such as hybridisation.

- Simultaneous rolling and tilting motion
- Speed control up to 70 rpm
- Six rollers holds load up to 4 kg

Tube Roller.

Art. No. N2400-7010

£696.33

OFFER: £487*

starlab.click/tuberoller

Pipette Stands.



ErgoOne® Universal Carousel Stand for up to 8 pipettes (UV Resistant)

For use with: ErgoOne®, Biohit®, Finn®, Gilson®, Rainin® single or multi-channel pipettes

S7200-0000
~~£147.22~~
£90*



ErgoOne® Self-Adhesive Shelf Holder for 1 pipette (UV Resistant)

For ErgoOne® single or multi-channel pipettes

E3099-0001
~~£7.51~~
£4*



ErgoOne® Wall Mount for 1 pipette (UV Resistant)


For use with: ErgoOne®, Eppendorf®, Finn®, Gilson®, Rainin®, Sartorius® single or multi-channel pipettes with finger hook

S7200-0001
~~£8.50~~
£6*



ErgoOne® E Individual Stand for 1 single-channel pipette

10 µl to 1000 µl G9999-1002 £68.58 £62*	5000 µl G9999-1003 £62.08 £56*
---	--



NEW! ErgoOne® E Shelf-Holder
for 1 pipette

For ErgoOne® E single or multi-channel pipettes, all volumes

G9999-1009
~~£19.77~~
£18*



NEW! ErgoOne® E Individual Stand
for 1 pipette

For ErgoOne® E single or multi-channel pipettes, all volumes

G9999-1011
~~£98.57~~
£88*



Universal Acrylic Pipette Stand for up to 5 pipettes (UV Resistant)

For use with: ErgoOne®, Biohit®, Finn®, Eppendorf®, Gilson®, Rainin® single or multi-channel pipettes

E3099-0000
~~£65.31~~
£32*



ErgoOne® Bench Support for 1 pipette

For ErgoOne® single-channel pipettes

E3099-0002
~~£16.86~~
£8*
Pack of 5.



ErgoOne® E Charging Stand for up to 3 pipettes

For ErgoOne® E single-channel pipettes (excluding 5 ml)


G9999-1001
~~£273.21~~
£186*



ErgoOne® E Individual Stand for 1 multi-channel pipette

For ErgoOne® E multi-channel pipettes (excluding 1250 µl)

G9999-1004
~~£36.80~~
£34*



NEW! ErgoOne® E Wall Mount
for 1 pipette

For ErgoOne® E single or multi-channel pipettes, all volumes

G9999-1010
~~£43.29~~
£32*

starlab.click/pipetestands

Pipettes not included.



***OFFER REF: Q24000671068**
must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

ErgoOne® FAST Pipette Controller.

FLEXIBLE, ACCURATE, SAFE TRANSFER



starlab.click/ergoonefast

ERGONOMY MEETS MODERNITY

Every aspect of the ErgoOne® FAST has been developed with ergonomics in mind.

ErgoOne® FAST

Art. No. S7166-0010

~~£382.64~~

OFFER: £240*



WATCH VIDEO

- Lightweight and ergonomic
- Designed to balance perfectly in the hand for maximum comfort
- Smooth operation via pressure sensitive control buttons
- Fast! Fills a 25 ml pipette in less than 4 seconds
- Long-life battery provides 25 x ~2000 aspiration cycles on a full charge
- Battery life indicator
- Battery maintains constant pipetting speed that does not diminish as battery charge declines
- Can be recharged during use
- Use with all glass or plastic pipettes, 1.0 ml to 100 ml

Discover our range of Serological Pipettes



starlab.click/serologicals



Pressure responsive speed controls



Autoclavable aspiration cone and silicone pipette holder ensures a secure seal for your liquids while the filter protects your pipette



Traffic light battery life indicator



Easy to replace battery



Benchtop support allows you to place the controller on the bench or hood floor without risk of contamination

ErgoOne® FAST Benchtop Support

Art. No: P7166-6306

~~£7.81~~

OFFER: £6*



***OFFER REF: Q24000671068**

must be stated on your PO (excl. eShop). Offer valid to 13.12.2024.

ErgoOne® E Electronic Pipettes.

INCREASE EFFICIENCY IN YOUR LAB

- Easy to use intuitive menu structure
- Lightweight, ergonomic design
- Adjustable finger rest for individual comfort
- Significantly reduced tip loading and ejection forces
- Universal tip fit

ErgoOne® E
from: **£384***

ErgoOne® E OPERATING MODES

Pipetting | Reverse Pipetting | Mixing
Dispense | Electrophoresis¹⁾

Single-Channel ErgoOne® E Pipettes

G9001-0010	G9001-0020	G9001-0200	G9001-1000	G9001-5000
0.5 – 10 µl	2.0 – 20 µl	10 – 200 µl	50 – 1,000 µl	250 – 5,000 µl
£498.75	£498.75	£498.75	£498.75	£498.75
£384*	£384*	£384*	£384*	£384*

8-Channel ErgoOne® E Pipettes

G9008-0010	G9008-0020	G9008-0100	G9008-0300	G9008-1250
0.5 – 10 µl	1.0 – 20 µl	5.0 – 100 µl	15 – 300 µl	50 – 1250 µl
£902.16	£902.16	£902.16	£902.42	£897.74
£705	£705	£705	£705	£705

12-Channel ErgoOne® E Pipettes

G9012-0010	G9012-0020	G9012-0100	G9012-0300	G9012-1250
0.5 – 10 µl	1.0 – 20 µl	5.0 – 100 µl	15 – 300 µl	50 – 1250 µl
£1076.97	£1076.97	£1076.97	£1076.97	£1047.91
£841	£841	£841	£841	£841

¹⁾ GEL mode is not available on ErgoOne® E models: 1000 µl, 1250 µl, 5000 µl.

MULTI-CHANNEL COMBO PACKS! 10 µl to 300 µl multi-channel pipettes supplied with an individual stand (G9999-1004) and one rack of TipOne® pipette tips. 1250 µl models supplied with a pipette shelf holder (G9999-1009) and two racks of TipOne® pipette tips.



NEW!
50 – 1250 µl
multi-channel
pipettes



See page 9
for stands

starlab.click/electronicpipettes



***OFFER REF: Q24000671068**
must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

ErgoOne® Manual Pipettes.



See page 9
for stands

ACCURACY AND PRECISION YOU CAN TRUST!

- Ergonomic and lightweight
- Minimal operating forces
- Volume lock
- Universal tip fit (excl. 5.0 ml)
- Colour coded

ErgoOne®

from: **£135***

 starlab.click/ergoone



WATCH VIDEO

Single-Channel Manual Pipettes

0.1 – 2.5 µl (red) S7100-0125 £229.00 £135*	0.5 – 10 µl (red) S7100-0510 £229.00 £135*	2 – 20 µl Micro (red) S7100-0221 £229.00 £135*	2 – 20 µl (yellow) S7100-0220 £229.00 £135*
10 – 100 µl (yellow) S7100-1100 £229.00 £135*	20 – 200 µl (yellow) S7100-2200 £229.00 £135*	100 – 1,000 µl (blue) S7110-1000 £229.00 £135*	500 – 5,000 µl (violet) S7150-5000 £229.00 £135*

8-Channel Manual Pipettes

0.5 – 10 µl (red) S7108-0510 £654.89 £430*	10 – 100 µl (yellow) S7108-1100 £654.89 £430*	30 – 300 µl (green) S7108-3300 £654.89 £430*
--	---	--

12-Channel Manual Pipettes

0.5 – 10 µl (red) S7112-0510 £780.61 £510*	10 – 100 µl (yellow) S7112-1100 £780.61 £510*	30 – 300 µl (green) S7112-3300 £780.61 £510*
--	---	--



***OFFER REF: Q24000671068**
must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

ErgoOne® Starter Packs.



Reliability is key so your pipette must be trusted to provide accurate and precise pipetting, day in, day out. ErgoOne® manual pipettes provide you with both the reliability you need and an ergonomic design that prevents fatigue over long periods to help keep your pipetting consistent. ErgoOne® pipettes have tighter specifications than ISO8655, and feature a volume lock to ensure the correct volumes in your reaction.

CONTENTS	STARTERPACK	STARTERPACK42.5	STARTERPACK410
0.1 – 2.5 µl ErgoOne® Single-Ch Pipette S7100-0125		●	
0.5 – 10 µl ErgoOne® Single-Ch Pipette S7100-0510			●
2.0 – 20 µl ErgoOne® Single-Ch Pipette, Standard Cone S7100-0220	●	●	●
20 – 200 µl ErgoOne® Single-Ch Pipette S7100-2200	●	●	●
100 – 1000 µl ErgoOne® Single-Ch Pipette S7110-1000	●	●	●
TipOne® Pipette Tips Three racks of sterile TipOne® tips per pack	●	●	●
Self-Adhesive Pipette Holder E3099-0001	●	●	●
Laboratory Timer One timer per pack	●	●	●
	OFFER* £675 £400*	£875 £530*	£875 £530*

ErgoOne® pipettes and Starlab's trusted TipOne® pipette tips create the perfect pipetting partnership.

 starlab.click/starterpacks

***OFFER REF: Q24000671068**
must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

Brady® Printer & Labels.



PRINTER FEATURES

- 133 symbols
- Barcoding options
- Auto label function cuts according to the size of the vial
- Automated serialisation
- Save up to 12 label styles
- Time and date stamp
- Three-line LCD display for clear visibility
- Auto shut-off saves battery life

Brady M210-LAB PRINTER

Blister pack containing Printer, Quick Start Guide, N5127-7425 label cartridge ready installed. Batteries not included. 2-year warranty.

Art. No: N3000-2000

List: ~~£171.43~~

OFFER: £104*

starlab.click/brady

HIGH QUALITY LABEL MATERIALS DESIGNED FOR THE LAB ENVIRONMENT

Label Cartridge Information & Ordering						Label Applications			
Label material	Colour	Width	Article No.	List	OFFER	Material specifications / applications	Vials	Slides	Plates
B-595 Vinyl	Various ¹⁾	9.5 mm	N3795-439#¹⁾	£32.39	£24*	Freezer: Min. -40 °C.	0.5 ml		✓
	White	9.5 mm	N3795-4399	£32.39	£24*		0.5 ml		✓
B-461 Polyester	White	9.5 mm	N3795-4619	£41.30	£29*	Freezer: Min. -80 °C. Liquid Nitrogen ²⁾ : -196 °C. Autoclave: 121 °C. Hot Water Baths: 100 °C. Barcoding (19.1 mm labels only) Low profile, extremely thin labels.	0.5 ml		✓
	White	12.7 mm	N5127-4619	£47.08	£34*		1.5 ml		✓
	White	19.1 mm	N7519-4619	£51.29	£39*		2.0 – 50 ml		✓
B-488 Polyester	White	9.5 mm	N3795-4889	£35.59	£27*	Freezer: Min. -80 °C. Autoclave: 121 °C. Barcoding (19.1 mm labels only)	0.5 ml	✓	✓
	White	12.7 mm	N5127-4889	£41.92	£32*		1.5 ml	✓	✓
	White	19.1 mm	N7519-4889	£50.62	£36*		2.0 – 50 ml	✓	✓
B-499 Nylon	White	9.5 mm	N3795-4999	£38.22	£26*	Freezer: Min. -80 °C. Liquid Nitrogen ²⁾ : -196 °C. Autoclave 121 °C.	0.5 ml	✓	✓
	White	12.7 mm	N5127-4999	£41.31	£31*		1.5 ml	✓	✓
	White	19.1 mm	N7519-4999	£46.00	£35*		2.0 – 50 ml	✓	✓
B-7425 Polypropylene	White	9.5 mm	N3795-7425	£45.24	£26*	Freezer: Min. -80 °C. Liquid Nitrogen ²⁾ : -196 °C. Autoclave: 121 °C. Hot Water Baths: 100 °C. Barcoding (19.1 mm labels only)	0.5 ml	✓	✓
	White	12.7 mm	N5127-7425	£50.76	£30*		1.5 ml	✓	✓
	White	19.1 mm	N7519-7425	£59.10	£35*		2.0 – 50 ml	✓	✓
B-423 Polyester	White	19.1 mm	N7519-4239	£41.58	£33*	Freezer: Min. -80 °C. Autoclave: 121 °C. Hot Water Bath: 100 °C Barcoding	2.0 – 50 ml	✓	✓
B-430 Polyester	Clear	19.1 mm	N7519-4300	£44.01	£33*	Freezer: Min. -80 °C. Autoclave: 121 °C. Hot Water Bath: 100 °C Barcoding	2.0 – 50 ml	✓	✓
B-403 Paper Water Dissolvable	White	19.1 mm	N7519-4039	£39.19	£30*	For temporary labelling/barcoding. Dissolves within 30 seconds of contact with warm water.	Erlenmeyers or any reusable container that requires frequent re-labelling		

¹⁾ To order colour required, replace # as last digit of Article No. with: 1 for blue, 2 for green, 4 for red, 6 for yellow.

²⁾ Acceptable for use in liquid nitrogen only if the ends of the label overlap.

All labels 6.4 m long, except B-499 Nylon which are 4.8 m long. All suitable for use with Brady printers: LabPal™, BMP™21, BMP™21-PLUS and M210-LAB.



***OFFER REF: Q24000671068**

must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

Erlenmeyer Flasks.

DESIGNED TO BE REUSED

Ideal for the cultivation of bacteria, yeast and fungi, as well as plant and animal cells in suspension.

PRODUCT FEATURES:

- Optically clear polycarbonate flasks
- Moulded graduations for easy volume measurement
- Leakproof and shatterproof
- Dual-function, vented polypropylene cap provides sterile air exchange or a tight seal
- Use for both aerobic and anaerobic cultures
- 0.22 µm PTFE cap membrane
- Individually-wrapped
- Sterility level (SAL 10⁻⁶) fulfils the highest demands
- Certified RNase, DNase, DNA and Pyrogen free
- Autoclavable and reusable up to 10 times



FLAT BASE



Description	Pack Size	Art. No.	List	OFFER
125 ml Erlenmeyer Flask	24	E5000-1012	£200.19	£156*
250 ml Erlenmeyer Flask	12	E5000-1025	£132.07	£103*
500 ml Erlenmeyer Flask	12	E5000-1050	£173.33	£135*
1000 ml Erlenmeyer Flask	6	E5000-1100	£121.13	£95*
2000 ml Erlenmeyer Flask	6	E5000-1200	£152.75	£119*
3000 ml Erlenmeyer Flask	3	E5000-1300	£182.87	£143*

BAFFLED BASE



125 ml Erlenmeyer Flask	24	E5001-1012	£200.19	£156*
250 ml Erlenmeyer Flask	12	E5001-1025	£132.07	£103*
500 ml Erlenmeyer Flask	12	E5001-1050	£173.33	£135*
1000 ml Erlenmeyer Flask	6	E5001-1100	£121.13	£95*
2000 ml Erlenmeyer Flask	6	E5001-1200	£152.75	£120*
3000 ml Erlenmeyer Flask	3	E5001-1300	£182.87	£143*

 starlab.click/erlenmeyers



DUAL-FUNCTION CAPS

Polypropylene caps with a 0.22 µm PTFE membrane, ideal for both aerobic and anaerobic cultures. Remove the top section to allow sterile exchange, and then replace when ready. Avoid contamination risk and the hassle of switching caps and stock only one flask type for all cultures.

Lab Timer.



4-CHANNEL TIMER WITH CLOCK

- Set up to 4 different timings
- Alarm auto switch-off saves batteries
- Bench stand and magnetic strip
- Large, easy to read display

Timer/Clock

Art. No. **N9900-0001**

List: **£21.93**

OFFER: £16*

 starlab.click/timer



***OFFER REF: Q24000671068**

must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

Vortex.



CLEARANCE OFFER
Available only while
stock lasts!

SIMPLE, YET SMART

- Variable speed control up to 2500 rpm
- Operating modes: continuous and touch
- Anti-slip feet
- Small footprint

Vortex


Art. No. N2400-6110

£258.20

OFFER: £137*

 starlab.click/vortex

-86 °C Freezer.



91 litre capacity
Ergonomic handle
USB port for
temperature logger

UNDERBENCH FREEZER

- Use as your personal -80 °C storage space
- No filter, maintenance free!
- 48-hour battery back-up
- Supplied with two adjustable shelves

Freezer

Art. No. N2700-1000

£6,636.17

OFFER: £5,442*

 starlab.click/underbench-freezer

5.0 ml Tips.



**The perfect fit
for your 5.0 ml
ErgoOne® E**

Also compatible with
Brand and Finnpiquette

See page 11
for ErgoOne® E
offer prices



LARGE VOLUME TIPS

- Slim design, perfect for pipetting from tall and narrow vessels
- Certified RNase, DNA and ATP free
- Supplied in an easy to handle round rack

5.0 ml Tips for ErgoOne® E

Art. No. E1070-2607 (1 x 28-tip rack)

£43.38

OFFER: £16.50*

 starlab.click/5ml-tip



T: +44 (0)1908 283800
infoline@starlab.co.uk
www.starlab.co.uk

Passionate for science.



***OFFER REF: Q24000671068**
must be stated on your PO (excl. eShop).
Offer valid to 13.12.2024.

