

SPECIAL *Offers*

January – March 2020

ErgoOne®



www.starlab.click/manualpipettes

ErgoOne® Pipettes

Universal tip fit. Fully Autoclavable. Volume Lock.



Single-Channels

S7100-0125 0.1–2.5 µl List: £242.00	S7100-0510 0.5–10 µl List: £242.00	S7100-0221 2–20 µl Micro List: £242.00	S7100-0220 2–20 µl Standard List: £242.00
S7100-1100 10–100 µl List: £242.00	S7100-2200 20–200 µl List: £242.00	S7110-1000 100–1000 µl List: £242.00	S7150-5000 500–5000 µl List: £242.00

SINGLE
CHANNELS
£108.00



8-CHANNELS

£360.00

12-CHANNELS

£430.00

8-Channels

S7108-0510 0.5–10 µl List: £593.85	S7108-1100 10–100 µl List: £593.85	S7108-3300 30–300 µl List: £593.85
---	---	---

12-Channels

S7112-0510 0.5–10 µl List: £675.00	S7112-1100 10–100 µl List: £675.00	S7112-3300 30–300 µl List: £675.00
---	---	---

IMPORTANT: OFFER REF: 18329ALL must be stated on your Purchase Order. If the Ref is not quoted against the relevant product line, discount will NOT be applied and CREDIT WILL NOT BE GIVEN RETROSPECTIVELY. A promotion code is not required for Starlab eShop orders. All prices exclude VAT. Offer prices valid until 31 March 2020.

Passionate for Science.



www.starlab.co.uk

ErgoOne® Starter Packs

ErgoOne® and TipOne® – the perfect pipetting partnership

PRICES FROM
£215.00



www.starlab.click/manualpipettes

ErgoOne® Starter Packs		
Description	Cat. No.	OFFER £
MICRO STARTER PACK 2 x ErgoOne® single-channel pipettes (2.5 µl, 10 µl); 2 x Self-Adhesive Pipette Shelf Holders; 2 x Racks Sterile TipOne® Tips; 1 x Lab Timer; 1 x Pen / Screen Stylus.	STARTMICRO	£215.00
STARTER PACK 3 x ErgoOne® single-channel pipettes (20 µl*, 200 µl, 1000 µl); 3 x Self-Adhesive Pipette Shelf Holders; 3 x Racks Sterile TipOne® Tips; 1 x Lab Timer; 1 x Pen / Screen Stylus.	STARTERPACK	£320.00
*standard cone S7100-0220		

ErgoOne® Starter Packs		
Description	Cat. No.	OFFER £
STARTER PACK 4 (2.5) 4 x ErgoOne® single-channel pipettes (2.5 µl, 20 µl*, 200 µl, & 1000 µl); 4 x Self-Adhesive Pipette Shelf Holders; 4 x Racks Sterile TipOne® Tips; 1 x Lab Timer; 1 x Pen / Screen Stylus.	STARTERPACK42.5	£425.00
STARTER PACK 4 (10) 4 x ErgoOne® single-channel pipettes (10 µl, 20 µl*, 200 µl, & 1000 µl); 4 x Self-Adhesive Pipette Shelf Holders; 4 x Racks Sterile TipOne® Tips; 1 x Lab Timer; 1 x Pen / Screen Stylus.	STARTERPACK410	£425.00
*standard cone S7100-0220		

TipOne® Sterile Refills

The economical and eco-friendly alternative to racked sterile tips

- ▶ Certified and sterile, individually-wrapped refills
- ▶ Completely enclosed design ensures contamination prevention
- ▶ Up to 65% less waste compared to racked sterile tips
- ▶ 100% recyclable - all components made from recyclable polypropylene
- ▶ Sterile refills available for both filter and standard tips
- ▶ Quick and easy to use



www.starlab.click/refills



Easy-open outer wrapper



Place complete refill into empty rack



Remove the lid and start pipetting

IMPORTANT: OFFER REF: 18329ALL must be stated on your Purchase Order. If the Ref is not quoted against the relevant product line, discount will NOT be applied and CREDIT WILL NOT BE GIVEN RETROSPECTIVELY. A promotion code is not required for Starlab eShop orders. All prices exclude VAT. Offer prices valid until 31 March 2020.

Thermomixer Mixer HC

LIST
£1753.64
£1425*

Mixing, Heating, Cooling, Timing: in one convenient unit

The Starlab Thermomixer, Mixer HC, is a multifunctional unit small enough to sit on your bench. Suitable for a wide range of applications, the Mixer HC is the all-round lab assistant that provides mixing, heating, cooling and timing in one convenient unit. Every bench should have one!

Mixer features

- ▶ Very small footprint
- ▶ Automatic thermoblock recognition
- ▶ 2-year warranty



www.starlab.click/mixerhc



Cat. No: S8012-0000

*Excludes Thermoblock

Mixer HC Thermoblocks		
Description	Cat. No.	OFFER £
For 24 x 0.5 ml Microcentrifuge Tubes ¹⁾	S8012-0010	£299.00
For 24 x 1.5 ml Microcentrifuge Tubes ¹⁾	S8012-0011	£299.00
For 24 x 2.0 ml Microcentrifuge Tubes ¹⁾	S8012-0012	£299.00
For 8 x 15 ml Centrifuge Tubes	S8012-0015	£299.00
For 4 x 50 ml Centrifuge Tubes	S8012-0016	£299.00

¹⁾ Price includes Iso Rack.



The Iso-Rack provides the quick transfer of multiple tubes in one go



Nine Thermoblocks are available to cover all popular tubes and plates

TubeOne® Tubes

PRICES FROM
£6.45



- ▶ Conservatively rated to 20,000 x g
- ▶ Temperature range: -80°C to 121°C (Autoclavable)
- ▶ Certified RNase, DNase, DNA and Pyrogen free

TubeOne® Microcentrifuge Tubes			
Description	Pack Size	Cat. No.	OFFER £
1.5 ml TubeOne® Tubes	500	S1615-5500	£6.45
2.0 ml TubeOne® Tubes	200	S1620-2700	£20.00



www.starlab.click/tubeone

Centrifuge Tubes

500 TUBES FROM
£37.00



- ▶ No polystyrene! Racked tubes supplied in cardboard trays
- ▶ USP Class VI, non-cytotoxic and sterile (SAL 10⁻⁶)
- ▶ Certified RNase, DNase, DNA and Pyrogen free

Centrifuge Tubes			
Description	Pack Size	Cat. No.	OFFER £
15 ml Conical, Loose (Sterile)	10 x 50	E1415-0100	£37.00
15 ml Conical, Racked (Sterile)	20 x 25	E1415-0700	£42.50
50 ml Conical, Loose (Sterile)	10 x 50	E1450-0100	£47.50
50 ml Conical, Racked (Sterile)	20 x 25	E1450-0700	£54.50
50 ml Skirted, Loose (Sterile)	20 x 25	E1450-0400	£54.50

www.starlab.click/centrifuge tubes

IMPORTANT: OFFER REF: 18329ALL must be stated on your Purchase Order. If the Ref is not quoted against the relevant product line, discount will NOT be applied and CREDIT WILL NOT BE GIVEN RETROSPECTIVELY. A promotion code is not required for Starlab eShop orders. All prices exclude VAT. Offer prices valid until 31 March 2020.

PCR Plates for Popular Cyclers

100 PLATES FROM
£130.00

Universal PCR Plates for use with a variety of machines




0.1 ml Non-Skirted Universal
compatible with most 0.1 ml cyclers (excl. FAST)

E1403-0200-C (100 Plates)
OFFER: £130.00




0.2 ml Non-Skirted Universal
compatible with most 0.2 ml cyclers (excl. FAST)

E1403-1200-C (100 Plates)
OFFER: £130.00



0.2 ml Semi-Skirted, Straight Edges, Universal
compatible with most 0.2 ml cyclers (excl. FAST)


I1402-9800-C (100 Plates)
OFFER: £130.00



0.2 ml Skirted Plate
also compatible with Eppendorf Gradient, Nexus family, Mastercycler Pro & Realplex


E1403-5200-C (100 Plates)
OFFER: £130.00

PCR Plates for ABI Systems



0.1 ml FAST Plate
StepOne Plus, Veriti 0, 1ml, 7500 FAST, 9800 FAST, Quant Studio 12K Flex, ViiA7

E1403-7700-C (100 Plates)
OFFER: £130.00



0.2 ml Semi-Skirted TAQ-MAN® style: use with most 0.2ml ABI, Techne & Thermo Hybaid cyclers (excl. FAST)


E1403-8200-C (100 Plates)
OFFER £130.00

PCR Plates for Roche® Lightcycler



96-Well Plate
for LightCycler 96 or LightCycler 480

I1402-9909-C (100 Plates)
OFFER: £140.00



384-Well Plate
for LightCycler 480

E1042-9909-C (100 Plates)
OFFER: £215.00


DOES OUR PLATE FIT YOUR CYCLER? See our compatibility chart ►

www.starlab.click/pcr-compatibility



PCR Plate Sealing

Silicone Sealing Mat (Autoclavable)



Reusable. Ultra Flexible. Piercable. Can be cut to size.

E1403-0000 (5 Mats)
(NOTE: not suitable for qPCR)

LIST
£34.50
£21.00

www.starlab.click/platesealing

X-Clear, Advanced Polyolefin StarSeal (qPCR)
E2796-9795 (100 Seals)

LIST
£137.65
£95.00

Clear, Polyolefin StarSeal (PCR)
E2796-9793 (100 Seals)

LIST
£49.11
£31.00

Aluminium StarSeal (PCR)
E2796-9792 (100 Seals)

LIST
£49.11
£31.00

IMPORTANT: OFFER REF: 18329ALL must be stated on your Purchase Order. If the Ref is not quoted against the relevant product line, discount will NOT be applied and CREDIT WILL NOT BE GIVEN RETROSPECTIVELY. A promotion code is not required for Starlab eShop orders. All prices exclude VAT. Offer prices valid until 31 March 2020.

CytoOne® Cell Culture Ware



For optimum cell growth and handling

Stringent testing! Culture tests are performed to assure proper cell attachment and growth. Every flask is pressure tested for leak-free assurance. Flasks, plates and dishes are tested for wettability and flatness, and visually inspected to eliminate defects and scratches. CytoOne® plates, flasks and dishes are certified and gamma sterilised (SAL 10⁻⁵).



CytoOne® TC-Treated Plates

- ▶ 360° chimney wells prevent cross-contamination
- ▶ Uniform temperature across all wells
- ▶ Vented skirt reduces condensation
- ▶ Optical clarity for distortion free microscopy

CytoOne® Plates, TC-Treated			
Description	Pack Size	Cat. No.	OFFER £
6-Well Plate, TC-Treated	50	CC7682-7506	£29.50
12-Well Plate, TC-Treated	50	CC7682-7512	£39.00
24-Well Plate, TC-Treated	50	CC7682-7524	£34.50
48-Well Plate, TC-Treated	50	CC7682-7548	£42.50
96-Well Plate, TC-Treated	50	CC7682-7596	£41.00



PRICES FROM
£29.50

CytoOne® Flasks with Vented Filter Caps

- ▶ 100% pressure tested
- ▶ Stacking rims with vented design for thermal equalisation
- ▶ Anti-tilt skirts and base bars improve stability
- ▶ Wide necks provide easy access

CytoOne® Flasks, Vented Filter Caps			
Description	Pack Size	Cat. No.	OFFER £
T-25 Flask, TC-Treated	30 × 10	CC7682-4825	£122.00
T-75 Flask, TC-Treated	20 × 5	CC7682-4875	£71.00
T-150 Flask, TC-Treated	8 × 5	CC7682-4815	£57.50
T-225 Flask, TC-Treated	5 × 5	CC7682-4822	£56.00
T-25 Flask, Non-Treated	30 × 10	CC7672-4825	£163.00
T-75 Flask, Non-Treated	20 × 5	CC7672-4875	£96.50
T-150 Flask, Non-Treated	8 × 5	CC7672-4815	£62.50
T-225 Flask, Non-Treated	5 × 5	CC7672-4822	£60.00



PRICES FROM
£56.00

TC-Treated flasks with vented filter caps

CytoOne® TC-Treated Dishes

- ▶ Optimal adherence and growth of adherent cells
- ▶ Vented design for thermal equalisation
- ▶ Easy-open sleeves in convenient small pack sizes
- ▶ Untreated lids to minimise condensation

CytoOne® Dishes, TC-Treated			
Description	Pack Size	Cat. No.	OFFER £
35 × 10 mm Dish, TC-Treated	10 × 30	CC7682-3340	£31.50
60 × 15 mm Dish, TC-Treated	10 × 30	CC7682-3359	£30.00
100 × 20 mm Dish, TC-Treated	10 × 30	CC7682-3394	£66.00
150 × 20 mm Dish, TC-Treated	5 × 12	CC7682-3614	£39.00



PRICES FROM
£30.00

www.starlab.click/cytoone

IMPORTANT: OFFER REF: 18329ALL must be stated on your Purchase Order. If the Ref is not quoted against the relevant product line, discount will NOT be applied and CREDIT WILL NOT BE GIVEN RETROSPECTIVELY. A promotion code is not required for Starlab eShop orders. All prices exclude VAT. Offer prices valid until 31 March 2020.

What's NEW? Microcentrifuge 24

LIST PRICE
£1516.50



Cat. No:
S8030-0000



www.starlab.click/microcentrifuge24

Small space, Big results!

Starlab's new Microcentrifuge 24 is the ideal centrifuge if you are looking for a versatile, personal microcentrifuge for all your routine molecular biology applications.


Product features

- ▶ Convenient and easy to use
- ▶ Small footprint
- ▶ Automatic lid opening at the end of the run
- ▶ Motorised lid closure
- ▶ "One-Touch" short-spin function
- ▶ Easy to clean keypad

Capacity

- ▶ Up to 21,300 g (15,060 rpm)
- ▶ Up to 24 x 2.0ml tubes
- ▶ Up to 32 x 0.2 ml PCR tubes
- ▶ Various adapters available


DNA/RNA EXTRACTION



DNA GEL EXTRACTION



DNA, PLASMID, PCR PURIFICATION




WIN A PIZZA LUNCH FOR YOUR LAB!



Complete the **STARLAB WORDSEARCH** for a chance to win a Pizza Lunch for your lab.

Scan the QR code now, or see the Special Offers section of our website to download the puzzle and enter the prize draw.

No purchase necessary. T&Cs apply.



STARLAB (UK) Ltd 5 Tanners Drive, Milton Keynes, MK14 5BU
Tel: 01908 283800 Fax: 01908 283802 infoline@starlab.co.uk

Passionate for Science.



www.starlab.co.uk