



Extra moulded graduations!
 15 ml Tube: 0 – 1.0 ml
 50 ml Tube: 0 – 3.0 ml



Centrifuge Tubes

STARLAB's new range of centrifuge tubes are designed for the centrifugation of cells, but can be used for standard centrifugation, preparation and storage applications. USP Class VI polymers and certified non-toxicity ensure the viability and integrity of cells. Maximum centrifugation rate: 12,000 x g.

Product features

- ▶ Clear, sterile polypropylene tubes with green HDPE screw caps
- ▶ White writing surface and clearly printed graduations
- ▶ Longer length screw caps with sealing ring prevents any leakage
- ▶ Perfect for centrifugation of cells: USP Class VI, non-cytotoxic
- ▶ Sterility level (SAL 10⁻⁶) fulfils the highest demands, meeting EN 556-1
- ▶ Certified RNase, DNase, DNA and Pyrogen free
- ▶ Racked tubes supplied in eco-friendly recyclable cardboard racks



ORDERING INFORMATION				
Description	Pack Size	Cat. No.	List £	OFFER £
15 ml Tubes				
15 ml Centrifuge Tubes, Conical, Loose (Sterile)	10 x 50	E1415-0100	£105.00	£37.50
15 ml Centrifuge Tubes, Conical, Racked (Sterile)	20 x 25	E1415-0700	£110.00	£42.95
50 ml Tubes				
50 ml Centrifuge Tubes, Conical, Loose (Sterile)	10 x 50	E1450-0100	£115.00	£47.50
50 ml Centrifuge Tubes, Conical, Racked (Sterile)	20 x 25	E1450-0700	£120.00	£53.50

IMPORTANT ORDERING INFORMATION

OFFER REF: 17934ALL must be stated on any Purchase Order that is emailed, faxed or posted. If the Ref is not quoted against the relevant product line, discount will not be applied and credit cannot be given retrospectively. A promotion code is not required for orders via **STARLAB eShop** or **eProcurement sites** that show the discounted pricing. Prices exclude VAT and valid until 31 December 2019.

