

END OF YEAR BUDGET MONEY TO SPEND...?

TREAT YOUR LAB TO SOME SHINY NEW EQUIPMENT & PIPETTES...

THERMOMIXER HC




S8102-0000
(with S8012-0011.
Thermoblocks sold separately.)

Mixer HC			
Cat. No.	Description	List Price	OFFER Price
S8012-0000	Mixer HC (excludes exchangeable Thermoblock)	£1,670.13	£1,335.00
10% OFF ACCESSORIES			
S8012-0010	Thermoblock for 24 x 0.5 ml microcentrifuge tubes*	£341.97	£307.77
S8012-0011	Thermoblock for 24 x 1.5 ml microcentrifuge tubes**	£341.97	£307.77
S8012-0012	Thermoblock for 24 x 2.0 ml microcentrifuge tubes**	£341.97	£307.77
S8012-0013	Thermoblock for 24 x micro tubes, diameter 12 mm	£429.38	£386.44
S8012-0014	Thermoblock for 8 x 5.0 ml preparation tubes	£341.97	£307.77
S8012-0015	Thermoblock for 8 x 15 ml conical centrifuge tubes	£341.97	£307.77
S8012-0016	Thermoblock for 4 x 50 ml conical centrifuge tubes	£341.97	£307.77
S8012-0017	Thermoblock for 24 x 1.5 ml – 2.0 ml cryogenic tubes	£341.97	£307.77
S8012-0018	Thermoblock for microplates / deepwell plates, with lid	£429.38	£386.44
S8012-0019	0.2ml Adaptor for use with S8012-0018	£117.35	£105.61

Thermoblocks sold separately. *Incl. ISO-Rack for 0.5ml tubes. ** Incl, ISO-Rack for 1.5/2.0ml tubes.

¹⁾ Excludes S8102-0013 and S8012-0018

MICROPLATE MIXER




Microplate Mixer			
Cat. No.	Description	List Price	OFFER Price
N2400-8040	Microplate Mixer	£457.23	£254.90
10% OFF ACCESSORIES			
N2400-8041	Attachment for one microtiter plate	£13.15	£11.84
N2400-8042	Attachment for two microtiter plates	£28.56	£25.70
N2400-8043	Microplate Mixer Adapter Plate (for use with N2400-6012 & N2400-6013)	£13.15	£11.84
N2400-6012	Tube Adapter for 18 x 1.5/2.0ml Microcentrifuge Tubes up to 10mm Ø (Can only be used in Combination with N2400-8043)	£45.72	£41.15
N2400-6013	Tube Adapter for 8 x 15 ml Centrifuge Tubes up to 16mm Ø (Can only be used in Combination with N2400-8043)	£45.72	£41.15

ROTATOR MIXER



Rotator Mixer			
Cat. No.	Description	List Price	OFFER Price
N2400-5000	Rotator Mixer	£643.75	£475.00

TUBE ROLLER



Tube Roller			
Cat. No.	Description	List Price	OFFER Price
N2400-7010	Tube Roller	£557.45	£335.00

IMPORTANT: OFFER REF: 18272ALL must be quoted on your Purchase Order against each item line if your PO is being emailed, faxed or posted. eShop orders: Offer Ref must be not in the "Special Requests" box at Checkout. No Offer Ref, no discount! Discount cannot be applied retrospectively. All prices excl. VAT. List prices valid to 31-Dec-19; **OFFER PRICES valid until 17-December-2019.**

Passionate for Science.



www.starlab.co.uk

VORTEX



VORTEX IR



ASPIRATOR



MINI FUGE



PLATE CENTRIFUGE



PCR CABINETS



FREEZER



Vortex			
Cat. No.	Description	List Price	OFFER Price
N2400-6110	Vortex with UK Plug	£206.69	£128.25
10% OFF ACCESSORIES			
N2400-6011	Replacement Universal Cup for Tubes up to 30mm Ø	£13.15	£11.84
N2400-6012	Tube Adapter for 18 x 1.5/2.0ml Microcentrifuge Tubes up to 10mm Ø (Can only be used in Combination with N2400-6015)	£45.72	£41.15
N2400-6013	Tube Adapter for 8 x 15 ml Centrifuge Tubes up to 16mm Ø) (Can only be used in Combination with N2400-6015)	£45.72	£41.15
N2400-6014	Adapter for 96-Well Microtitre Plates	£90.20	£81.18
N2400-6015	Vortex Adapter Plate (for use with N2400-6012 & N2400-6013)	£32.57	£29.31

Vortex IR			
Cat. No.	Description	List Price	OFFER Price
N2400-0001	Vortex IR	£280.21	£181.65
10% OFF ACCESSORIES			
N2400-1001	Attachment for 19 x 1.5ml Microcentrifuge Tubes	£58.80	£52.92
N2400-1002	Attachment for 5 x 15ml Centrifuge Tubes	£58.80	£52.92
N2400-1003	Attachment for 4 x 50ml Centrifuge Tubes	£58.80	£52.92
N2400-1004	Attachment for one Microtitre Plate	£58.80	£52.92
N2400-1005	Universal Rubber Attachment	£32.75	£29.48
N2400-1006	Customisable Soft Foam Attachment	£47.84	£43.06

Aspirator			
Cat. No.	Description	List Price	OFFER Price
N2400-9000	Aspirator (incl. 4L PP Bottle, Hand Controller and Wall Mount)	£1,560.45	£1,190.00
N2400-9001	Hand Controller for Aspirator (incl. adapter for seros, adapter for tips, wall mount, 2.5m tubing)	£134.64	£121.18
N2400-9006	Connection Kit with Inlet Pipe for Additional Hand Controller	£31.21	£28.09
N2400-9009	Quick Coupling (Bottle to Pump)	£87.26	£78.53


Microcentrifuge			
Cat. No.	Description	List Price	OFFER Price
N2631-0017	Mini Fuge PLUS	£253.42	£145.32

Plate Centrifuge			
Cat. No.	Description	List Price	OFFER Price
N2631-0008	Plate Centrifuge	£626.34	£482.20
N2631-0018	Adapter for 48-Well Plates and 0.2 ml Strips and Tubes. Pack of 2.	£44.21	£39.79

AirClean PCR Cabinets			
Cat. No.	Description	List Price	OFFER Price
N2540-8400	32" Dead Air PCR Cabinet with UV Light	£2,786.15	£2,300.00
N2540-8500	48" Dead Air PCR Cabinet with UV Light	£3,517.07	£2,700.00
N2530-8200	32" Laminar Flow PCR Cabinet with UV Light	£3,415.84	£2,550.00
N2530-8300	48" Laminar Flow PCR Cabinet with UV Light	£4,621.55	£3,450.00


-86°C Underbench Freezer			
Cat. No.	Description	List Price	OFFER Price
N2700-1000	Underbench Freezer	£5,005.00	£4,322.00
N2700-1001	U100 Rack for 50mm High Boxes, 74 places	£1,200.00	£1,080.00
N2700-1002	U100 Rack for 75mm High Boxes, 45 places	£1,200.00	£1,080.00
N2700-1003	U100 Rack for 100mm High Boxes, 33 places	£1,200.00	£1,080.00

ROCKING/3D ROCKING SHAKERS




ORBITAL SHAKER

N2400-8030
(with N2400-8011)



DRY BATH SYSTEMS



N2400-0001

N2400-0002
(with N2400-4108)

MINI DRY BATH



N2400-4020
(with N2400-4023)

BATH BEADS



MAGNETIC STIRRERS



Rocking & 3D Rocking Shakers

Cat. No.	Description	List Price	OFFER Price
N2400-8010	Rocking Shaker (incl. platform N2400-8021)	£1,262.51	£741.50
N2400-8020	3D Rocking Shaker (incl. platform N2400-8021)	£1,262.51	£741.50
N2400-8011	Optional Universal Platform with Clamping Rolls	£305.81	£177.75

Orbital Shaker

Cat. No.	Description	List Price	OFFER Price
N2400-8030	Orbital Shaker*	£1,107.77	£648.85
N2400-8011	Universal Platform with Clamping Rolls	£305.81	£177.75
N2400-8021	Universal Platform with Rubber Mat and Bungee Cords	£305.81	£177.75

*IMPORTANT: N2400-8030 is NOT supplied with a platform. A platform must be purchased separately.

Dry Bath Systems

Cat. No.	Description	List Price	OFFER Price
N2400-4001	Dry Bath System (Single Block)*	£343.92	£258.20
N2400-4002	Dry Bath System (Dual Block)*	£497.76	£324.50

10% OFF ALL METAL BLOCKS

N2400-4100	Metal Block for 64 x 0.2ml Tubes	£93.83	£84.45
N2400-4101	Metal Block for 20 x 1.5/2.0ml Tubes	£93.83	£84.45
N2400-4102	Metal Block for 12 x 1.5ml Tubes	£93.83	£84.45
N2400-4103	Metal Block for 4 x 50ml Tubes	£93.83	£84.45
N2400-4104	Metal Block for 20 x 1.5ml Tubes (V-Bottom)	£93.83	£84.45
N2400-4105	Metal Block for 20 x 0.5ml Tubes	£93.83	£84.45
N2400-4106	Metal Block for 20 Tubes up to 13mm diameter	£93.83	£84.45
N2400-4107	Metal Block for 96-Well Plate (for use only with N2400-4001)	£156.38	£140.74
N2400-4108	Metal Block for 96-Well Plate (for use only with N2400-4002)	£212.78	£191.50

*Metal blocks not included and must be purchased separately.

Mini Dry Bath

Cat. No.	Description	List Price	OFFER Price
N2400-4020	Mini Dry Bath Incubator	£236.04	£140.00

10% OFF ALL METAL BLOCKS

N2400-4021	Mini Dry Bath Metal Block for 32 x 0.2ml PCR Tubes	£48.17	£43.35
N2400-4022	Mini Dry Bath Metal Block for 12 x 0.5ml Microcentrifuge Tubes	£48.17	£43.35
N2400-4023	Mini Dry Bath Metal Block for 12 x 1.5ml Microcentrifuge Tubes	£48.17	£43.35
N2400-4024	Mini Dry Bath Metal Block for 12 x 1.5/2.0ml Microcentrifuge Tubes	£48.17	£43.35
N2400-4025	Mini Dry Bath Metal Block for 6 x 15ml Centrifuge Tubes	£48.17	£43.35
N2400-4026	Mini Dry Bath Metal Block for 2 x 50ml Centrifuge Tubes	£48.17	£43.35

*Metal blocks not included and must be purchased separately.

Bath Beads

Cat. No.	Description	List Price	OFFER Price
N2400-4900	Bath Beads 0.5L	£138.54	£96.00
N2400-4901	Bath Beads 4.0L	£540.75	£420.22
N2400-4902	Bath Beads 12L	£1,699.50	£1,134.80

Magnetic Stirrers

Cat. No.	Description	List Price	OFFER Price
N2400-2010	Magnetic Stirrer	£195.42	£112.75
N2400-3010	Magnetic Stirrer with Heating	£344.49	£220.95

PIPETTE Offers FOR YOUR END OF YEAR BUDGET SPEND



- ▶ Fully autoclavable
- ▶ Universal tip fit
- ▶ Colour coded
- ▶ Volume lock
- ▶ Ergonomic comfort
- ▶ Great accuracy
- ▶ Great precision
- ▶ Great price



ErgoOne® Starter Packs			OFFER £
MICRO STARTER PACK: TWO ErgoOne® single-channel pipettes (2.5µl, 10µl). TWO Self-Adhesive Pipette Shelf Holders. TWO Racks Sterile TipOne® Tips. ONE Lab Timer. ONE Pen / Screen Stylus.	STARTMICRO	£446.34	£210.00
STARTER PACK: THREE ErgoOne® single-channel pipettes (20µl Standard, 200µl, 1000µl). THREE Self-Adhesive Pipette Shelf Holders. THREE Racks Sterile TipOne® Tips. ONE Lab Timer. ONE Pen / Screen Stylus.	STARTERPACK	£639.21	£315.00
STARTER PACK 4 (2.5): FOUR ErgoOne® single-channel pipettes (2.5µl, 20µl Standard, 200µl, & 1000µl). FOUR Self-Adhesive Pipette Shelf Holders. FOUR Racks Sterile TipOne® Tips. ONE Lab Timer. ONE Pen / Screen Stylus.	STARTERPACK42.5	£861.07	£420.00
STARTER PACK 4 (10): FOUR ErgoOne® single-channel pipettes (10µl, 20µl Standard, 200µl, & 1000µl). FOUR Self-Adhesive Pipette Shelf Holders. FOUR Racks Sterile TipOne® Tips. ONE Lab Timer. ONE Pen / Screen Stylus.	STARTERPACK410	£861.07	£420.00

8-CHANNEL

S7108-0510 0.5–10µl List: £555.00	S7108-1100 10–100µl List: £555.00	S7108-3300 30–300µl List: £555.00
--	--	--

8-CH
£355
EACH

12-CHANNEL

S7112-0510 0.5–10µl List: £656.00	S7112-1100 10–100µl List: £656.00	S7112-3300 30–300µl List: £656.00
--	--	--

12-CH
£420
EACH



SINGLE-CHANNEL

- S7100-0125**
0.1–2.5µl
List: £212.00 each
▶ **OFFER: £105**
- S7100-0510**
0.5–10µl
List: £212.00 each
▶ **OFFER: £105**
- S7100-0221**
2–20µl Micro
List: £212.00 each
▶ **OFFER: £105**
- S7100-0220**
2–20µl Standard
List: £212.00 each
▶ **OFFER: £105**
- S7100-1100**
10–100µl
List: £197.00 each
▶ **OFFER: £105**
- S7100-2200**
20–200µl
List: £197.00 each
▶ **OFFER: £105**
- S7110-1000**
100–1000µl
List: £197.00 each
▶ **OFFER: £105**
- S7150-5000**
500–5000µl
List: £216.00 each
▶ **OFFER: £105**

ONLY
£105
EACH

ELECTRONIC PIPETTES

ErgoOne® E Electronic Pipettes			OFFER £
Single-Channel			
0.5–10µl	G9001-0010	£421.79	£295.00
2.0–20µl	G9001-0020	£421.79	£295.00
10–200µl	G9001-0200	£421.79	£295.00
50–1000µl	G9001-1000	£421.79	£295.00
500–5000µl	G9001-5000	£448.82	£300.00



ErgoOne® E Electronic Pipettes			OFFER £
8-Channel			
0.5–10µl	G9008-0010	£702.98	£510.00
1.0–20µl	G9008-0020	£702.98	£510.00
5.0–100µl	G9008-0100	£702.98	£510.00
15–300µl	G9008-0300	£702.98	£510.00
12-Channel			
0.5–10µl	G9012-0010	£811.13	£610.00
1.0–20µl	G9012-0020	£811.13	£610.00
5.0–100µl	G9012-0100	£811.13	£610.00
15–300µl	G9012-0300	£811.13	£610.00



S7166-0010
ErgoOne® FAST
Pipette Controller
List Price: £285.00 each
▶ **OFFER: £197.50**

IMPORTANT: OFFER REF: 18272ALL must be quoted on your Purchase Order against each item line if your PO is being emailed, faxed or posted. eShop orders: Offer Ref must be not in the "Special Requests" box at Checkout. No Offer Ref, no discount! Discount cannot be applied retrospectively. All prices excl. VAT. List prices valid to 31-Dec-19; OFFER PRICES valid until 17-December-2019.

Passionate for Science.



www.starlab.co.uk