

LAB START UP Offer

2019

SETTING UP A NEW LAB?
JUST RECEIVED YOUR GRANT?



LET US
HELP YOU
EQUIP
YOUR LAB



Discounted prices on **EQUIPMENT & PIPETTES**

- ▶ Spend **£2000** and get **£200** consumables free
- ▶ Spend **£3000** and get **£300** consumables free
- ▶ Spend **£5000** and get **£500** consumables free

Terms apply.

Passionate for Science.



www.starlab.co.uk

THERMOMIXER

www.starlab.click/mixerhc

Mixer HC			
Cat. No.	Description	List £	OFFER £
S8012-0000	Mixer HC ¹⁾	£1,670.13	£1380.13
10% OFF THERMOBLOCKS			
S8012-0010	For 24 x 0.5 ml microcentrifuge tubes ²⁾	£341.97	£307.77
S8012-0011	For 24 x 1.5 ml microcentrifuge tubes ³⁾	£341.97	£307.77
S8012-0012	For 24 x 2.0 ml microcentrifuge tubes ³⁾	£341.97	£307.77
S8012-0013	For 24 x micro tubes, diameter 12 mm	£429.38	£386.44
S8012-0014	For 8 x 5.0 ml preparation tubes	£341.97	£307.77
S8012-0015	For 8 x 15 ml conical centrifuge tubes	£341.97	£307.77
S8012-0016	For 4 x 50 ml conical centrifuge tubes	£341.97	£307.77
S8012-0017	For 24 x 1.5 ml – 2.0 ml cryogenic tubes	£341.97	£307.77
S8012-0018	For microplates / deepwell plates, includes lid	£429.38	£386.44
S8012-0019	Adapter plate 96 x 0.2ml (for use with S8012-0018)	£117.35	£105.62
S8012-0020	ISO-Rack for 24 x 0.5 ml Microcentrifuge Tubes	£49.51	£44.56
S8012-0021	ISO-Rack for 24 x 1.5/2.0 ml Microcentrifuge Tubes	£48.75	£43.88

¹⁾ Thermoblock is not included ²⁾ Includes 1 x S8012-0020 ISO-Rack for 0.5 ml tubes.
³⁾ Includes 1 x S8012-0020 ISO-Rack for 1.5/2.0 ml tubes.

ASPIRATOR

www.starlab.click/aspirator

Aspirator			
Cat. No.	Description	List £	OFFER £
N2400-9000	Aspirator (incl. 4L PP Bottle, Hand Controller and Wall Mount)	£1,560.45	£1272.57
10% OFF SELECTED ACCESSORIES			
N2400-9001	Hand Controller (incl. adapter for seros, adapter for tips, wall mount, 2.5m tubing)	£134.64	£121.18
N2400-9006	Connection Kit with Inlet Pipe for Additional Hand Controller	£31.21	£28.09
N2400-9009	Quick Coupling (Bottle to Pump)	£87.26	£78.54

MIXERS

www.starlab.click/rotatormixer

Rotator Mixer			
Cat. No.	Description	List £	OFFER £
N2400-5000	Rotator Mixer	£643.75	£508.69

www.starlab.click/microplate-mixer

Microplate Mixer			
Cat. No.	Description	List £	OFFER £
N2400-8040	Microplate Mixer	£457.23	£274.51
10% OFF ACCESSORIES			
N2400-8041	Attachment for one microtiter plate	£13.15	£11.84
N2400-8042	Attachment for two microtiter plates	£28.56	£25.70
N2400-8043	Microplate Mixer Adapter Plate (for use with N2400-6012 & N2400-6013)	£13.15	£11.84
N2400-6012	Tube Adapter for 18 x 1.5/2.0ml Microcentrifuge Tubes up to 10mm ¹⁾	£45.72	£41.15
N2400-6013	Tube Adapter for 8 x 15 ml Centrifuge Tubes up to 16mm Ø ¹⁾	£45.72	£41.15

¹⁾ Can only be used in Combination with N2400-8043)

www.starlab.click/tuberoller

Tube Roller			
Cat. No.	Description	List £	OFFER £
N2400-7010	Tube Roller	£557.45	£344.38

VORTEXERS

www.starlab.click/vortex

Vortex			
Cat. No.	Description	List £	OFFER £
N2400-6110	Vortex with UK Plug	£206.69	£138.09
10% OFF ACCESSORIES			
N2400-6011	Replacement Universal Cup for Tubes up to 30mm Ø	£13.15	£11.84
N2400-6012	Tube Adapter for 18 x 1.5/2.0ml Microcentrifuge Tubes up to 10mm ¹⁾	£45.72	£41.15
N2400-6013	Tube Adapter for 8 x 15 ml Centrifuge Tubes up to 16mm Ø ¹⁾	£45.72	£41.15
N2400-6014	Adapter for 96-Well Microtitre Plates	£90.20	£81.18
N2400-6015	Vortex Adapter Plate (for use with N2400-6012 & N2400-6013)	£32.57	£29.31

¹⁾ Can only be used in Combination with N2400-6015

www.starlab.click/vortexir

Vortex IR			
Cat. No.	Description	List £	OFFER £
N2400-0001	Vortex IR	£280.21	£195.65
10% OFF ACCESSORIES			
N2400-1001	Attachment for 19 x 1.5ml Microcentrifuge Tubes	£58.80	£52.92
N2400-1002	Attachment for 5 x 15ml Centrifuge Tubes	£58.80	£52.92
N2400-1003	Attachment for 4 x 50ml Centrifuge Tubes	£58.80	£52.92
N2400-1004	Attachment for one Microtitre Plate	£58.80	£52.92
N2400-1005	Universal Rubber Attachment	£32.75	£29.48
N2400-1006	Customisable Soft Foam Attachment	£47.84	£43.06

SHAKERS

www.starlab.click/orbitalshaker

Orbital Shaker			
Cat. No.	Description	List £	OFFER £
N2400-8030	Orbital Shaker ¹⁾	£1,107.77	£698.75
N2400-8011	Universal Platform with Clamping Rolls	£305.81	£191.32
N2400-8021	Universal Platform with Rubber Mat and Bungee Cords	£305.81	£191.32

¹⁾ IMPORTANT NOT supplied with a platform. A platform must be purchased separately.

www.starlab.click/shakers

Rocking & 3D Rocking Shakers			
Cat. No.	Description	List £	OFFER £
N2400-8010	Rocking Shaker (includes platform N2400-8021)	£1,262.51	£798.59
N2400-8020	3D Rocking Shaker (includes platform N2400-8021)	£1,262.51	£798.59
N2400-8011	Optional Universal Platform with Clamping Rolls	£305.81	£191.32

STIRRERS

www.starlab.click/magnetic-stirrer

Magnetic Stirrers			
Cat. No.	Description	List £	OFFER £
N2400-2010	Magnetic Stirrer	£195.42	£121.46
N2400-3010	Magnetic Stirrer with Heating	£344.49	£237.91

DRY BATH INCUBATORS

www.starlab.click/dry-bath-systems

Dry Bath Systems			
Cat. No.	Description	List £	OFFER £
N2400-4001	Dry Bath System (Single Block) ¹⁾	£343.92	£278.08
N2400-4002	Dry Bath System (Dual Block) ¹⁾	£497.76	£349.39
10% OFF METAL BLOCKS			
N2400-4100	Metal Block for 64 x 0.2ml Tubes	£93.83	£84.45
N2400-4101	Metal Block for 20 x 1.5/2.0ml Tubes	£93.83	£84.45
N2400-4102	Metal Block for 12 x 1.5ml Tubes	£93.83	£84.45
N2400-4103	Metal Block for 4 x 50ml Tubes	£93.83	£84.45
N2400-4104	Metal Block for 20 x 1.5ml Tubes (V-Bottom)	£93.83	£84.45
N2400-4105	Metal Block for 20 x 0.5ml Tubes	£93.83	£84.45
N2400-4106	Metal Block for 20 Tubes up to 13mm diameter	£93.83	£84.45
N2400-4107	Metal Block for 96-Well Plate (for use only with N2400-4001)	£156.38	£140.74
N2400-4108	Metal Block for 96-Well Plate (for use only with N2400-4002)	£112.78	£191.50

¹⁾ Metal blocks not included and must be purchased separately.

www.starlab.click/mini-dry-bath

Mini Dry Bath			
Cat. No.	Description	List £	OFFER £
N2400-4020	Mini Dry Bath Incubator	£236.04	£149.74
10% OFF METAL BLOCKS			
N2400-4021	Mini Dry Bath Metal Block for 32 x 0.2 ml PCR Tubes	£48.17	£43.35
N2400-4022	Mini Dry Bath Metal Block for 12 x 0.5 ml Microcentrifuge Tubes	£48.17	£43.35
N2400-4023	Mini Dry Bath Metal Block for 12 x 1.5 ml Microcentrifuge Tubes	£48.17	£43.35
N2400-4024	Mini Dry Bath Metal Block for 12 x 1.5/2.0 ml Microcentrifuge Tubes	£48.17	£43.35
N2400-4025	Mini Dry Bath Metal Block for 6 x 15 ml Centrifuge Tubes	£48.17	£43.35
N2400-4026	Mini Dry Bath Metal Block for 2 x 50 ml Centrifuge Tubes	£48.17	£43.35

*Metal blocks not included and must be purchased separately.

BATH BEADS

www.starlab.click/bathbeads

Bath Beads			
Cat. No.	Description	List £	OFFER £
N2400-4900	Bath Beads 0.5L	£138.54	£103.38
N2400-4901	Bath Beads 4.0L	£540.75	£452.57
N2400-4902	Bath Beads 12L	£1,699.50	£1,222.15

UNDERBENCH FREEZER

www.starlab.click/underbench-freezer

-86°C Underbench Freezer			
Cat. No.	Description	List Price	OFFER Price
N2700-1000	Underbench Freezer	£5,005.00	£4,654.74
10% OFF ACCESSORIES			
N2700-1001	U100 Rack for 50mm High Boxes, 74 places	£1,200.00	£1,108.92
N2700-1002	U100 Rack for 75mm High Boxes, 45 places	£1,200.00	£1,108.92
N2700-1003	U100 Rack for 100mm High Boxes, 33 places	£1,200.00	£1,107.52

CENTRIFUGES

www.starlab.click/minifuge-plus

Microcentrifuge			
Cat. No.	Description	List Price	OFFER Price
N2631-0017	Mini Fuge PLUS	£253.42	£156.51

www.starlab.click/platecentrifuge

Plate Centrifuge			
Cat. No.	Description	List Price	OFFER Price
N2631-0008	Plate Centrifuge	£626.34	£519.28
10% OFF ACCESSORIES			
N2631-0018	Adapter for 48-Well Plates and 0.2 ml Strips and Tubes (Pack of 2)	£44.21	£39.79

PCR CABINETS

www.starlab.click/airclean

AirClean PCR Cabinets			
Cat. No.	Description	List Price	OFFER Price
N2540-8400	32" Dead Air PCR Cabinet with UV Light	£2,786.15	£2,405.31
N2540-8500	48" Dead Air PCR Cabinet with UV Light	£3,517.07	£2,813.80
N2530-8200	32" Laminar Flow PCR Cabinet with UV Light	£3,415.84	£2,618.91
N2530-8300	48" Laminar Flow PCR Cabinet with UV Light	£4,621.55	£3,646.85

TIMER

www.starlab.click/timer

4-Channel Timer			
Cat. No.	Description	List £	OFFER £
N9900-0001	Four Channel Timer / Clock	£17.27	£12.00

MANUAL PIPETTES

www.starlab.click/manualpipettes

ErgoOne® Single-Channel Pipettes			
Cat. No.	Description	List £	OFFER £
S7100-0125	0.1-2.5µl ErgoOne Single-Channel	£212.00	£120.00
S7100-0220	2-20µl ErgoOne Single-Channel ¹⁾	£212.00	£120.00
S7100-0221	2-20µl ErgoOne Single Channel ²⁾	£212.00	£120.00
S7100-0510	0.5-10µl ErgoOne Single Channel	£212.00	£120.00
S7100-1100	10-100µl ErgoOne Single Channel	£197.00	£120.00
S7100-2200	20-200µl ErgoOne Single Channel	£197.00	£120.00
S7110-1000	100-1000µl ErgoOne Single Channel	£197.00	£120.00
S7150-5000	0.5-5ml ErgoOne Single Channel	£216.00	£120.00

¹⁾ with micro cone. ²⁾ with standard cone.

ErgoOne® Multi-Channel Pipettes			
Cat. No.	Description	List £	OFFER £
8-CHANNEL			
S7108-0510	0.5-10µl ErgoOne 8-Channel	£555.00	£375.00
S7108-1100	10-100µl ErgoOne 8-Channel	£555.00	£375.00
S7108-3300	30-300µl ErgoOne 8-Channel	£555.00	£375.00
12-CHANNEL			
S7112-0510	0.5-10µl ErgoOne 12-Channel	£656.00	£445.00
S7112-1100	10-100µl ErgoOne 12-Channel	£656.00	£445.00
S7112-3300	30-300µl ErgoOne 12-Channel	£656.00	£445.00

PIPETTE CONTROLLER

www.starlab.click/pipettecontroller

ErgoOne® Single-Channel Pipettes			
Cat. No.	Description	List £	OFFER £
S7166-0010	ErgoOne FAST Pipette Controller	£285.00	£227.28

ELECTRONIC PIPETTES

www.starlab.click/electronicpipettes

ErgoOne® E Single-Channel Pipettes			
Cat. No.	Description	List £	OFFER £
G9001-0010	ErgoOne E Single Channel 0.5-10µl	£421.79	£300.00
G9001-0020	ErgoOne E Single Channel, 2-20µl	£421.79	£300.00
G9001-0200	ErgoOne E Single Channel, 20-200µl	£421.79	£300.00
G9001-1000	ErgoOne E Single Channel, 100-1000µl	£421.79	£300.00
G9001-5000	ErgoOne E Single Channel, 0.5-5ml	£448.82	£300.00

ErgoOne® E Multi-Channel Pipettes			
Cat. No.	Description	List £	OFFER £
8-CHANNEL			
G9008-0010	ErgoOne E 8-Channel, 0.5-10µl	£702.98	£550.00
G9008-0020	ErgoOne E 8-Channel, 1-20µl	£702.98	£550.00
G9008-0100	ErgoOne E 8-Channel, 5-100µl	£702.98	£550.00
G9008-0300	ErgoOne E 8-Channel, 15-300µl	£702.98	£550.00
12-CHANNEL			
G9012-0010	ErgoOne E 12-Channel, 0.5-10µl	£811.13	£650.00
G9012-0020	ErgoOne E 12-Channel, 1-20µl	£811.13	£650.00
G9012-0100	ErgoOne E 12-Channel, 5-100µl	£811.13	£650.00
G9012-0300	ErgoOne E 12-Channel, 15-300µl	£811.13	£650.00

HOW IT WORKS

- Select the equipment/pipettes you need at the discounted pricing outlined in this flyer.
- Decide which free consumables you would like to add to your package.
- Contact your Account Manager who will arrange the quote for you, or if needed, they can arrange a meeting with you to discuss your requirements
- Your Lab Start Up Quote will be emailed to you.
- You place your Order (quoting the Quote Reference No. you have been given on the PO)
- Your shiny new equipment/pipettes and free consumables are delivered to your lab

TERMS OF OFFER

Discounted pricing shown on this order can only be used for Lab Start Up Quotes with a minimum order value of £2000; Pricing shown on this offer is not valid for ordering individual items.

A quote will be raised for all required items; all items on the quote must be ordered at the same time on one Purchase Order; The Quote Reference No. must be quoted on PO at time of ordering.

Free consumables are at customer's choice; value of the free consumables is excluding VAT and based on list price, not on any discounted price that may already be agreed with the customer. Consumables is defined as anything in our product portfolio that is not part of the Pipettes, Pipette Servicing or Smart Instruments categories of product.

Pricing shown excludes VAT. Offer valid until 31 December 2019.