

Starlab Special Offers



Screw Cap Tubes with caps

Art. No.	Description	Pack Size	List Price	Your Price
E1405-0301	0.5 ml Plain Skirted Tube with EasyGrip Screw Cap White/Blue (Sterile)	5 x 50	£59.38	£39.19
E1405-0302	0.5 ml Plain Skirted Tube with EasyGrip Screw Cap White/Green (Sterile)	5 x 50	£59.38	£38.00
E1405-0304	0.5 ml Plain Skirted Tube with EasyGrip Screw Cap White/Red (Sterile)	5 x 50	£59.38	£38.00
E1405-0306	0.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Yellow (Sterile)	5 x 50	£59.38	£38.00
E1405-2141	0.5 ml Plain Skirted Tube, Natural Standard Screw Cap (Sterile)	5 x 50	£59.32	£39.15
E1415-0301	1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Blue (Sterile)	5 x 50	£59.38	£38.00
E1415-0302	1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Green (Sterile)	5 x 50	£59.38	£38.00
E1415-0304	1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Red (Sterile)	5 x 50	£59.38	£38.00
E1415-0306	1.5ml Graduated Conical Tube with EasyGrip Screw Cap White/Yellow (Sterile)	5 x 50	£59.38	£38.00
E1415-0311	1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Blue (Sterile)	5 x 50	£59.38	£38.00
E1415-0312	1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Green (Sterile)	5 x 50	£59.38	£38.00
E1415-0314	1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Red (Sterile)	5 x 50	£59.38	£38.00
E1415-0316	1.5ml Plain Skirted Tube with EasyGrip Screw Cap White/Yellow (Sterile)	5 x 50	£63.49	£38.09
E1415-2231	1.5 ml Graduated Conical Tube, Natural Standard Screw Cap (Sterile)	5 x 50	£54.60	£36.04
E1415-2241	1.5 ml Plain Skirted Tube, Natural Standard Screw Cap (Sterile)	5 x 50	£54.60	£36.04
E1415-2243	1.5 ml Plain, Skirted Tube, Natural EasyGrip Cap (Sterile)	5 x 50	£59.38	£38.00
E1420-0301	2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Blue (Sterile)	5 x 50	£80.98	£53.45
E1420-0302	2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Green (Sterile)	5 x 50	£80.98	£51.83
E1420-0304	2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Red (Sterile)	5 x 50	£80.98	£51.83
E1420-0306	2.0ml Graduated Conical Tube with EasyGrip Screw Cap White/Yellow (Sterile)	5 x 50	£80.98	£40.49
E1420-0311	2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Blue (Sterile)	5 x 50	£59.38	£38.00
E1420-0312	2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Green (Sterile)	5 x 50	£59.38	£38.00
E1420-0314	2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Red (Sterile)	5 x 50	£59.38	£38.00
E1420-0316	2.0ml Graduated Skirted Tube with EasyGrip Screw Cap White/Yellow (Sterile)	5 x 50	£59.38	£38.00
E1420-2331	2.0 ml Graduated Conical Tube, Natural Standard Screw Cap (Sterile)	5 x 50	£55.55	£36.66
E1420-2341	2.0 ml Graduated Skirted Tube, Natural Standard Screw Cap (Sterile)	5 x 50	£55.42	£36.02
E1420-2343	2.0 ml Graduated Skirted Tube, Natural EasyGrip Cap (Sterile)	5 x 50	£59.38	£38.00

Please use Q26000927108 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Starlab Special Offers



Screw Cap Tubes only

Art. No.	Description	Pack Size	List Price	Your Price
E1405-2130	0.5 ml Screw Cap Tube, Plain, Conical, Natural	500	£30.10	£19.87
E1405-2140	0.5 ml Screw Cap Tube, Plain, Skirted, Natural	500	£29.40	£18.23
E1405-2147	0.5 ml Screw Cap Tube, Plain, Skirted, Amber	500	£35.11	£22.47
E1415-2230	1.5 ml Screw Cap Tube, Graduated, Conical, Natural	500	£26.64	£15.72
E1415-2240	1.5 ml Screw Cap Tube, Plain, Skirted, Natural	500	£26.46	£15.35
E1420-2330	2.0 ml Screw Cap Tube, Graduated, Conical, Natural	500	£26.92	£15.62
E1420-2340	2.0 ml Screw Cap Tube, Graduated, Skirted, Natural	500	£26.60	£14.90



Ultra High Recovery Microcentrifuge Tubes

Art. No.	Description	Pack Size	List Price	Your Price
E1405-2600	0.5 ml Ultra High Recovery Microcentrifuge Tube	500	£27.64	£18.24
E1415-2600	1.5 ml Ultra High Recovery Microcentrifuge Tube	250	£17.36	£11.46
I1420-2600	2.0 ml Ultra High Recovery Microcentrifuge Tube	250	£20.45	£13.50



Centrifuge Tubes

Art. No.	Description	Pack Size	List Price	Your Price
E1415-0100	15 ml Centrifuge Tube, Conical, Loose (Sterile)	500	£146.05	£52.58
E1415-0700	15 ml Centrifuge Tube, Conical, Racked (Sterile)	500	£151.87	£60.75
E1450-0100	50 ml Centrifuge Tube, Conical, Loose (Sterile)	500	£155.42	£68.38
E1450-0400	50 ml Centrifuge Tube, Skirted, Loose (Sterile)	500	£146.06	£75.95
E1450-0700	50 ml Centrifuge Tube, Conical, Racked (Sterile)	500	£162.58	£76.88

Please use Q26000927108 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Starlab Special Offers



Screw Caps only


Art. No.	Description	Pack Size	List Price	Your Price
E1480-0100	Standard Caps for Screw Cap Microtubes, natural	500	£52.91	£29.63
E1480-0101	Standard Cap for Screw Cap Microtubes, Blue	500	£56.10	£31.42
E1480-0102	Standard Cap for Screw Cap Microtubes, Green	500	£56.10	£32.54
E1480-0103	Standard Cap for Screw Cap Microtubes, Orange	500	£56.10	£32.54
E1480-0104	Standard Caps for Screw Cap Microtubes, red	500	£56.10	£31.42
E1480-0105	Standard Cap for Screw Cap Microtubes, Violet	500	£56.10	£31.42
E1480-0106	Standard Cap for Screw Cap Microtubes, Yellow	500	£56.10	£31.42
E1480-0107	Standard Cap for Screw Cap Microtubes, Amber	500	£56.10	£31.42
E1480-0108	Standard Cap for Screw Cap Microtubes, Opaque Black	500	£56.93	£38.71
E1480-0109	Standard Cap for Screw Cap Microtubes, Opaque White	500	£53.46	£30.47
E1480-0300	EasyGrip Cap, Natural / White	500	£53.15	£30.83
E1480-0301	Easy Grip Cap, White / Blue	500	£55.54	£32.21
E1480-0302	Easy Grip Cap, White / Green	500	£55.54	£32.21
E1480-0303	Easy Grip Cap, White / Orange	500	£55.54	£32.21
E1480-0304	Easy Grip Cap, White / Red	500	£55.54	£32.21
E1480-0305	Easy Grip Cap, White / Violet	500	£55.54	£32.21
E1480-0306	Easy Grip Cap, White / Yellow	500	£55.54	£32.21
E1480-0309	Easy Grip Cap, White / White	500	£55.54	£32.21
E1480-0399	Easy Grip Caps, Mixed	500	£55.54	£32.21

Cryogenic Vials - Cap Inserts

Art. No.	Description	Pack Size	List Price	Your Price
E1494-6001	Coloured Cap Inserts for Cryogenic Vials, Blue	500	£30.63	£19.60
E1494-6002	Coloured Cap Inserts for Cryogenic Vials, Green	500	£30.63	£19.60
E1494-6003	Coloured Cap Inserts for Cryogenic Vials, Orange	500	£30.63	£19.60
E1494-6004	Coloured Cap Inserts for Cryogenic Vials, Red	500	£30.63	£19.60
E1494-6005	Coloured Cap Inserts for Cryogenic Vials, Violet	500	£30.63	£19.60
E1494-6006	Coloured Cap Inserts for Cryogenic Vials, Yellow	500	£30.63	£19.60
E1494-6007	Coloured Cap Inserts for Cryogenic Vials, Brown	500	£30.63	£19.60
E1494-6009	Coloured Cap Inserts for Cryogenic Vials, White	500	£30.63	£19.60
E1494-6099	Coloured Cap Inserts for Cryogenic Vials, Mixed	500	£30.63	£19.60
E1494-6104	Coloured Cap Inserts for Cryogenic Vials, Pink	500	£30.63	£19.60
E1494-6108	Coloured Cap Inserts for Cryogenic Vials, Grey	500	£30.63	£19.60


Please use Q26000927108 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Starlab Special Offers



Cardboard Storage Boxes

Art. No.	Description	Pack Size	List Price	Your Price
E2023-1615	Cardboard Storage Box for 15 ml Centrifuge Tubes, 36 places	1	£10.77	£4.96
E2023-1650	Cardboard Storage Box for 50 ml Centrifuge Tubes, 16 places	1	£10.77	£4.96
E2023-3100	Cardboard Storage Box, 133x133 Footprint, White, 100-Place Grid	1	£4.53	£3.26
E2023-8181	Cardboard Storage Box, 133x133 Footprint, White, 81-Place Grid	1	£4.53	£3.08



Cryogenic Storage Boxes

Art. No.	Description	Pack Size	List Price	Your Price
E2350-4001	50 Place Lock-Rack, Blue	4	£107.68	£48.46
E2350-4004	50 Place Lock-Rack, Red	4	£107.68	£48.46
E2350-4005	50 Place Lock-Rack, Violet	4	£107.68	£48.46
I2510-5005	100-place Polycarbonate Storage Box	5	£54.33	£27.17



Tube Racks

Art. No.	Description	Pack Size	List Price	Your Price
E2380-5000	StarRack 80 polypropylene fraction collection rack, Natural	5	£22.21	£14.44
E2380-5048	StarRack 80 polypropylene fraction collection rack, Neon Mixed	5	£23.00	£15.00
E2396-0010	StarRack 96 double-sided polypropylene rack, Natural	Each	£13.57	£8.82
E2396-1515	Racks for 15ml centrifuge tubes (15-place)	1	£29.59	£25.73
E2396-1550	Racks for 50 ml centrifuge tubes (8-place)	1	£38.38	£24.95
I2396-5048	StarRack 96 double-sided polypropylene rack, Neon Mixed	5	£69.99	£45.49



Cryogenic Vials

Art. No.	Description	Pack Size	List Price	Your Price
E3090-6212	1.0 ml Cryovial Vials with External Thread, Self-Sealing Cap, Skirted (Sterile)	10 x 50	£161.86	£113.30
E3090-6222	1.8 ml Cryovial Vials with External Thread, Self Sealing Seal Cap, Skirted (Sterile)	10 x 50	£161.39	£98.45
E3090-6242	4.5 ml Cryovial Vials with External Thread, Self-Sealing Cap, Skirted (Sterile)	6 x 50	£123.00	£79.23
E3090-6321	1.8 ml Cryovial Vials with External Thread, Self-Sealing Cap Conical (Sterile)	10 x 50	£168.91	£113.17
E3110-6112	1.0 ml Cryovial with Internal Thread, Silicone Seal Cap, Skirted (Sterile)	10 x 50	£189.74	£132.82
E3110-6121	1.8 ml Cryovial with Internal Thread, Silicone Seal Cap Conical (Sterile)	10 x 50	£208.84	£146.18
E3110-6122	1.8 ml Cryovial with Internal Thread, Silicone Seal Cap, Skirted (Sterile)	10 x 50	£192.06	£119.07
E3110-6131	3.6 ml Cryovial with Internal Thread, Silicone Seal Cap, Conical (Sterile)	8 x 50	£163.14	£104.41
E3110-6132	3.6 ml Cryovial with Internal Thread, Silicone Seal Cap, Skirted (Sterile)	8 x 50	£170.61	£119.43
E3110-6141	4.5 ml Cryovial with Internal Thread, Silicone Seal Cap, Conical (Sterile)	6 x 50	£164.17	£108.35
E3110-6142	4.5 ml Cryovial with Internal Thread, Silicone Seal Cap, Skirted (Sterile)	6 x 50	£164.17	£114.92



Please use Q26000927108 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Starlab Special Offers




Polypropylene Storage Boxes

Art. No.	Description	Pack Size	List Price	Your Price
E9025-1650	Centrifuge Tube Storage / Transport Box for 50ml Tubes, Blue	set of 2	£36.80	£23.92
E9025-3615	Centrifuge Tube Storage / Transport Box for 15ml Tubes, Green	set of 2	£36.81	£23.93
I2310-5800	StarStore 100-Place Storage Box, Natural	5	£36.26	£24.30
I2310-5848	StarStore 100-Place Storage Box, Neon Mixed*	5	£35.25	£22.91
I2381-5040	StarStore 81, Natural	5	£36.23	£22.46
I2381-5041	StarStore 81, Neon Blue	5	£37.97	£23.54
I2381-5042	StarStore 81, Neon Green	5	£37.97	£23.54
I2381-5043	StarStore 81, Neon Orange	5	£37.97	£23.54
I2381-5044	StarStore 81, Neon Pink	5	£37.97	£23.54
I2381-5046	StarStore 81, Neon Yellow	5	£37.97	£23.54
I2381-5048	StarStore 81, Neon Mixed	5	£37.97	£24.68



Startags

Art. No.	Description	Pack Size	List Price	Your Price
E9185-0001	StarTags Laboratory Labels (9.5 mm Diameter spot), blue	2000	£28.36	£15.60
E9187-1009	StarTags Laboratory Labels, Large (32.5 mm x 12.5 mm)	1000	£28.36	£15.60
E9188-1009	StarTags Laboratory Labels, Medium (24 mm x 12.5 mm)	1000	£28.36	£15.60
E9188-8699	StarTag Laboratory Labels, Small (20 mm x 6 mm)	1000	£28.36	£15.60



Cryogenic Storage

Art. No.	Description	Pack Size	List Price	Your Price
I2381-2305	81-place Polycarbonate Storage Box	5	£47.98	£23.99



TubeOne®

Art. No.	Description	Pack Size	List Price	Your Price
S1605-0000	0.5 ml TubeOne® Microcentrifuge Tubes, Natural (non-sterile)	5 x 200	£43.56	£26.14
S1615-5500	1.5 ml TubeOne® Microcentrifuge Tubes, Natural (non-sterile)	500	£17.67	£8.30
S1615-5510	1.5 ml TubeOne® Microcentrifuge Tubes, Natural (Sterile)	5 x 100	£36.82	£22.09
S1615-5550	1.5 ml TubeOne® Microcentrifuge Tubes, Natural (non-sterile)	5 x 100	£19.34	£9.85
S1620-2700	2.0ml TubeOne® Microcentrifuge Tubes, Natural	500	£62.01	£31.21

Please use Q26000927108 when ordering. Offer valid to 27.03.2026 / while stock lasts.

Starlab Special Offers



Manual Pipettes

Art. No.	Description	Pack Size	List Price	Your Price
S7100-0125	ErgoOne® Single Channel Pipette, 0.1 - 2.5µl	1	£246.50	£137.50
S7100-0220	ErgoOne® Single Channel Pipette (Standard Cone), 2 - 20µl	1	£246.50	£137.50
S7100-0221	ErgoOne® Single Channel Pipette (Micro Cone), 2 - 20µl	1	£246.50	£137.50
S7100-0510	ErgoOne® Single Channel Pipette, 0.5 - 10µl	1	£246.50	£137.50
S7100-1100	ErgoOne® Single Channel Pipette, 10 - 100µl	1	£246.50	£137.50
S7100-2200	ErgoOne® Single Channel Pipette, 20 - 200µl	1	£246.50	£137.50
S7108-0510	ErgoOne® 8-Channel Pipette, 0.5 - 10µl	1	£704.92	£385.00
S7108-1100	ErgoOne® 8-Channel Pipette, 10 - 100µl	1	£704.92	£385.00
S7108-3300	ErgoOne® 8-Channel Pipette, 30 - 300µl	1	£704.92	£385.00
S7110-1000	ErgoOne® Single Channel Pipette, 100 - 1000µl	1	£246.50	£137.50
S7112-0510	ErgoOne® 12-Channel Pipette, 0.5 - 10µl	1	£840.25	£470.00
S7112-1100	ErgoOne® 12-Channel Pipette, 10 - 100µl	1	£840.25	£470.00
S7112-3300	ErgoOne® 12-Channel Pipette, 30 - 300µl	1	£840.25	£470.00
S7150-5000	ErgoOne® Single Channel Pipette, 0.5 - 5ml	1	£246.50	£137.50

Equipment

Art. No.	Description	Pack Size	List Price	Your Price
S8010-0000	Starlab Vortex Mixer	1	£373.84	£280.38
S8030-5000	Starlab Minicentrifuge	1	£249.91	£187.43

Please use Q26000927108 when ordering. Offer valid to 27.03.2026 / while stock lasts.