

# TipOne® sterile refills.



Made from recycled material

Reuse your TipOne® racks with our refills!



## PLASTIC SAVING REFILLS

- Fully enclosed sterile refill
- Saves up to 63% plastic compared to racked tips
- Space saving of up to 20% compared to racked tips
- No need to autoclave – I'm fully sterile
- Made in a factory run by sustainable energy
- Recyclable. Join our unique closed-loop recycling service!

Description	Pack Size	Art. No.	OFFER*
10/20µl XL Graduated TipOne® Pipette Tip, Sterile Refill	10 x 96	S1110-3710	£23.91
200µl Bevelled TipOne® Pipette Tip, Sterile Refill	10 x 96	S1111-1710	£23.45
200µl UltraPoint® Graduated TipOne® Pipette Tip, Sterile Refill	10 x 96	S1113-1710	£27.56
300µL, Graduated TipOne® Pipette Tip, Sterile Refill	10 x 96	S1110-9710	£34.50
1250µl XL Graduated TipOne® Pipette Tip, Sterile Refill	10 x 96	S1112-1730	£42.88

## Win a TipOne® lego set!

Enter our online competition for the chance to win! Scan the QR code or visit the link

Competition closes:  
29-05-2026



ENTER HERE

starlab.click/win-lego



## Erlenmeyer flasks.

### IDEAL WHEN SAFETY OR WEIGHT ARE PRIORITIES

These optically clear, polycarbonate flasks are safer, and visibly clearer than glass.

- Sterile (SAL 10<sup>-6</sup>), individually-wrapped
- Shatter resistant
- Dual-function cap featuring 0.22 µm PTFE membrane means you can confidently use the same flask for both aerobic cultures (vented) or anaerobic cultures (closed)
- No more foil or bungs that can breed contamination



REUSABLE! Withstands at least 10 autoclave cycles!

	Description	Pack Size	Art. No.	OFFER*
FLAT BASE	125 ml Erlenmeyer Flask, Flat Base	24	E5000-1012	£134.71
	250 ml Erlenmeyer Flask, Flat Base	12	E5000-1025	£88.87
	500 ml Erlenmeyer Flask, Flat Base	12	E5000-1050	£116.64
	1000 ml Erlenmeyer Flask, Flat Base	6	E5000-1100	£81.51
	2000 ml Erlenmeyer Flask, Flat Base	6	E5000-1200	£102.79
BAFFLED BASE	3000 ml Erlenmeyer Flask, Flat Base	3	E5000-1300	£123.06
	125 ml Erlenmeyer Flask, Baffled Base	24	E5001-1012	£134.71
	250 ml Erlenmeyer Flask, Baffled Base	12	E5001-1025	£88.87
	500 ml Erlenmeyer Flask, Baffled Base	12	E5001-1050	£116.64
	1000 ml Erlenmeyer Flask, Baffled Base	6	E5001-1100	£81.51
	2000 ml Erlenmeyer Flask, Baffled Base	6	E5001-1200	£102.79
	3000 ml Erlenmeyer Flask, Baffled Base	3	E5001-1300	£123.06

## Reagent reservoirs.

Autoclavable and reusable, reducing plastic waste. Compatible with all multi-channel pipettes.

Description	Pack Size	Article No.	OFFER*
50/60 ml StarTub, PP Reservoir. Non-sterile. V-shaped bottom.	10	S4026-5806	£14.44

Made in Germany. "Local for local" production reduces carbon footprint.

**\*OFFER REF: Q26000965642**  
must be stated on your PO (excl. eShop).  
Offer valid to 29.05.26 / while stock lasts.

Built for your lab,  
made to reduce impact.  
You are passionate  
for science.

## SPECIAL OFFER HIGHLIGHTS

Valid until 29.05.2026



SCAN TO  
BUY ONLINE

IN THIS EDITION:  
PIPETTE TIPS  
MANUAL PIPETTES  
ERLENMEYER FLASKS  
FREEZER RACKING



# Repeat dispensing.



Made in Germany. "Local for local" production reduces carbon footprint.

Reusable saving on plastic waste.

## RepeatOne S®

Art. No. S4026-1500

**OFFER: £228.28**

Dispenser tips not included.

### YOUR FIRST CHOICE FOR CHALLENGING LIQUIDS

With up to 48 dispenses from the same sample, you can use the RepeatOne® S to dispense precise volumes of liquid in a series, without the need to aspirate in between each step. 1-year warranty.

## TipOne® Repeat dispenser tips.

Tip volume	Pack Size	Article No.	OFFER*	Tip Volume	Pack Size	Article No.	OFFER*
<b>Non-Sterile</b>				<b>Sterile, Individually Wrapped</b>			
0.1 ml	100	S4751-0010	£47.03	0.1 ml	100	S4761-0010	£80.31
0.5 ml	100	S4751-0050	£36.73	0.5 ml	100	S4761-0050	£71.45
1.0 ml	100	S4751-0100	£35.02	1.0 ml	100	S4761-0100	£71.45
2.5 ml	100	S4751-0250	£35.02	2.5 ml	100	S4761-0250	£71.45
5.0 ml	100	S4751-0500	£35.87	5.0 ml	100	S4761-0500	£71.45
10 ml	100	S4751-1000	£37.58	10 ml	100	S4761-1000	£71.45
25 ml <sup>1)</sup>	100	S4751-2500	£62.63	25 ml <sup>1)</sup>	100	S4761-2500	£102.47
50 ml <sup>1)</sup>	100	S4751-5000	£75.24	50 ml <sup>1)</sup>	100	S4761-5000	£108.95

<sup>1)</sup> 25 ml and 50 ml dispenser tips are supplied with four adapters.  
Note: TipOne® Repeat Dispenser Tips are not autoclavable.

# ErgoOne® manual pipettes.

### MADE FOR DURABILITY!

ErgoOne® manual pipettes provide accuracy and precision you can trust, alongside ergonomic comfort and a stylish look you will love. ErgoOne pipettes have a greased system, which means they rarely need their sealing system changed. Less parts = less waste! 3-year warranty.



Made in Germany. "Local for local" production reduces carbon footprint.

### Single-Channel Pipettes

<b>S7100-0125</b> 0.1–2.5µl <b>£135</b>	<b>S7100-0510</b> 0.5–10µl <b>£135</b>	<b>S7100-0221</b> 2–20µl Micro <b>£135</b>	<b>S7100-0220</b> 2–20µl Standard <b>£135</b>
<b>S7100-1100</b> 10–100µl <b>£135</b>	<b>S7100-2200</b> 20–200µl <b>£135</b>	<b>S7110-1000</b> 100–1,000µl <b>£135</b>	<b>S7150-5000</b> 500–5,000µl <b>£135</b>

### 8-Channel Pipettes

<b>S7108-0510</b> 0.5–10µl <b>£402</b>	<b>S7108-1100</b> 10–100µl <b>£402</b>	<b>S7108-3300</b> 30–300µl <b>£402</b>
--	--	--



**ErgoOne® carousel**  
Holds up to 8 pipettes. For use with ErgoOne®, Biohit®, Finn®, Gilson® and Rainin®, single and multi-channel pipettes. UV resistant.

Art. No: S7200-0000  
**£94.35**



**ErgoOne® Fast pipette controller.**  
The ErgoOne® FAST is one of the lightest pipette controllers on the market. Designed for your comfort, it balances perfectly in the hand.

Art. No: S7166-0010  
**£253.39**

# Freezer racking

### SUSTAINABILITY JUST GOT COOLER

- Optimise space in your freezer, reduce energy use and cut your running costs
- Manufactured using 100% sustainable energy
- Manufactured in Germany — shorter supply chain links reduce carbon footprint
- Racks specifically designed for Eppendorf, Thermo, PHCbi, Haier and Nordic Labs freezers

**MADE FROM 85% RECYCLED MATERIAL**

### RACKS WITH SIDE ACCESS



Boxes not included.

For use with **standard storage boxes up to 136 x 136 x 53 mm.**

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No	OFFER*
9-Place	3 deep x 3 high	139 x 426 x 169 mm	E6000-0033	£104.64
12-Place	3 deep x 4 high	139 x 426 x 224 mm	E6000-0034	£119.65
15-Place	3 deep x 5 high	139 x 426 x 297 mm	E6000-0035	£132.08
12-Place	4 deep x 3 high	139 x 565 x 169 mm	E6000-0043	£118.38
16-Place	4 deep x 4 high	139 x 565 x 224 mm	E6000-0044	£118.38
20-Place	4 deep x 5 high	139 x 565 x 279 mm	E6000-0045	£135.48
24-Place	4 deep x 6 high	139 x 565 x 334 mm	E6000-0046	£155.12



Boxes not included.



Boxes not included.

### RACKS WITH PULL OUT SHELVES

For use with all **standard storage boxes up to 136 x 136 x 53 mm.**

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No	OFFER*
9-Place	3 deep x 3 high	140 x 450 x 173 mm	E6001-0033	£143.25
12-Place	3 deep x 4 high	140 x 450 x 231 mm	E6001-0034	£148.43
15-Place	3 deep x 5 high	140 x 450 x 288 mm	E6001-0035	£153.63
12-Place	4 deep x 3 high	140 x 569 x 174 mm	E6001-0043	£153.63
16-Place	4 deep x 4 high	140 x 569 x 231 mm	E6001-0044	£158.79
20-Place	4 deep x 5 high	140 x 569 x 281 mm	E6001-0045	£163.96
24-Place	4 deep x 6 high	140 x 569 x 337 mm	E6001-0046	£176.41
28-Place	4 deep x 7 high	140 x 569 x 403 mm	E6001-0047	£187.81

For use with Starlab's **StarStore 100 storage boxes only.**

Capacity	Format	Ext. Dimensions (WxDxH)	Art. No	OFFER*
9-Place	3 deep x 3 high	147 x 500 x 179 mm	E6001-0133	£161.64
12-Place	3 deep x 4 high	147 x 500 x 239 mm	E6001-0134	£159.16
16-Place	4 deep x 4 high	147 x 655 x 239 mm	E6001-0144	£172.56

Racking suitable for use with StarStore 100 storage boxes (WxDxH) 140 x 155 x 55 mm.

### RACKING FOR CHEST FREEZERS

For use with all **standard storage boxes up to 136 x 136 x 53 mm.**

Capacity	Ext. Dimensions (WxDxH)	Art. No	OFFER*
9-High SS Rack, Locking Rod	141 x 148 x 501 mm	E6010-0009	£89.36
10-High SS Rack, Locking Rod	141 x 148 x 556 mm	E6010-0010	£89.36
11-High SS Rack, Locking Rod	141 x 148 x 611 mm	E6010-0011	£94.49
12-High SS Rack, Locking Rod	141 x 148 x 666 mm	E6010-0012	£99.89
13-High SS Rack, Locking Rod	141 x 148 x 721 mm	E6010-0013	£105.61

**\*OFFER REF: Q26000965642**

must be stated on your PO (excl. eShop). Offer valid to 29.05.26 / while stock lasts.