

**Trusted for consistency,  
designed for growth.**  
**You are passionate  
for science.**

**SPECIAL OFFER HIGHLIGHTS**

Valid until 31.08.2026



SCAN TO  
BUY ONLINE

**IN THIS EDITION:  
FLASKS, DISHES & PLATES  
BOTTLE TOP FILTRATION  
ASPIRATOR  
TUBES & VIALS FOR CELL CULTURE**



# CytoOne® Dishes.

FROM: **£39.00\***



CytoOne® TC-treated dishes feature stacking rings with vented design for thermal equalisation. Untreated lids reduce condensation: ideal for microscopy.

Size	Pack Size	Art. No.	OFFER*
TC-Treated Dishes			
35 x 10 mm	30 x 10	<b>CC7682-3340</b>	<b>£41.00</b>
60 x 15 mm	30 x 10	<b>CC7682-3359</b>	<b>£39.00</b>
100 x 20 mm	30 x 10	<b>CC7682-3394</b>	<b>£87.00</b>
150 x 20 mm	12 x 5	<b>CC7682-3614</b>	<b>£51.00</b>

Certified and sterile (SAL 10<sup>-9</sup>).

All CytoOne® products are certified RNase, DNase, DNA & Pyrogen free and Gamma sterilised (SAL 10<sup>-5</sup>). TC-Treated for optimum cell attachment and growth.

# CytoOne® Plates.

FROM: **£37.50\***



Complete 360° chimney well design prevents cross-contamination and ensures a uniform temperature transfer across all wells to prevent the edge effect. Integrated vents and stacking rims create space for continuous airflow between plates and reduces evaporation.

Wells	Pack Size	Art. No.	OFFER*
TC-Treated Plates			
6-Well	50	<b>CC7682-7506</b>	<b>£37.50</b>
12-Well	50	<b>CC7682-7512</b>	<b>£50.50</b>
24-Well	50	<b>CC7682-7524</b>	<b>£45.00</b>
48-Well	50	<b>CC7682-7548</b>	<b>£55.00</b>
96-Well	50	<b>CC7682-7596</b>	<b>£52.00</b>

Certified and sterile (SAL 10<sup>-9</sup>).

# CytoOne® Flasks.

FROM: **£59.50\***

All flasks supplied in convenient 'mini-packs' in easy-open, resealable zip top bags that protect unused flasks.



Vented filter caps with a 0.2 µm hydrophobic membrane for gas exchange and to exclude bacteria and fungal contaminants.

Size	Pack Size	Art. No.	OFFER*
TC Treated Flasks, Vented Filter Caps			
T-25	30 x 10	<b>CC7682-4825</b>	<b>£210.00</b>
T-75	20 x 5	<b>CC7682-4875</b>	<b>£125.00</b>
T-150	8 x 5	<b>CC7682-4815</b>	<b>£81.00</b>
T-225	5 x 5	<b>CC7682-4822</b>	<b>£78.00</b>

Certified and sterile (SAL 10<sup>-9</sup>).



Two-position plug seal caps allow you to choose secure manual venting or an airtight seal.

Size	Pack Size	Art. No.	OFFER*
TC Treated Flasks, Plug Seal Caps			
T-25	30 x 10	<b>CC7682-4325</b>	<b>£123.00</b>
T-75	20 x 5	<b>CC7682-4175</b>	<b>£83.50</b>
T-150	8 x 5	<b>CC7682-4415</b>	<b>£66.00</b>
T-225	5 x 5	<b>CC7682-4225</b>	<b>£59.50</b>

Certified and sterile (SAL 10<sup>-9</sup>).

# Centrifuge Tubes.

Safeguard your cell viability. Designed to support your most sensitive cell work. Graduated to full 15 ml and 50 ml volumes, with extra graduations on the conical base.

**Extra moulded graduations!**  
15 ml Conical Tube: 0 – 1.0 ml  
50 ml Conical Tube: 0 – 3.0 ml



Description	Pack Size	Art. No.	OFFER*
15 ml Conical, Loose	20 x 25	<b>E1415-0100</b>	<b>£52.00</b>
15 ml Conical, Racked	20 x 25	<b>E1415-0700</b>	<b>£65.00</b>
50 ml Conical, Loose	20 x 25	<b>E1450-0100</b>	<b>£70.00</b>
50 ml Conical, Racked	20 x 25	<b>E1450-0700</b>	<b>£83.00</b>
50 ml Skirted, Loose	20 x 25	<b>E1450-0400</b>	<b>£77.00</b>

Certified and sterile (SAL 10<sup>-9</sup>).

**MICROTUBE PESTLE**  
for 1.5/2.0ml tubes.  
Assorted colours.  
Certified. Pack 10.  
**I1415-5390**  
**£7.00\***



**CULTURE TUBES**  
14 ml tubes with dual-position caps that can be left loose for aerobic work or sealed for anaerobic culture.

**Polypropylene tube**  
-190 °C / 3000 xg.  
20 x 25-Tube bags  
**I1485-0810**  
**£130.00\***

**Polystyrene tube**  
-10 °C / 1400 xg  
20 x 25-Tube bags  
**I1485-2810**  
**£127.00\***



**\*OFFER REF: Q26000992193**  
must be stated on your PO (excl. eShop).  
Offer valid to 31.08.26 / while stock lasts.

# Aspirator.

Includes 4L PP Bottle & Hand Controller.  
**Art. No. N2400-9000**  
**OFFER: £1450\***

# Aspiration Pipettes.

Volume	Qty	Art. No.	OFFER*
2.0 ml	500	<b>E4861-0002</b>	<b>£53.00</b>
5.0 ml	200	<b>E4861-0005</b>	<b>£25.00</b>
10 ml	200	<b>E4861-0010</b>	<b>£33.00</b>

Certified RNase, DNase, DNA, Latex and Pyrogen free.  
Sterile (SAL 10<sup>-9</sup>). Individually wrapped.



# Serological Pipettes.



Vol.	Qty	Art. No.	OFFER*
1.0 ml	500	<b>E4860-0001</b>	<b>£38.50*</b>
2.0 ml	500	<b>E4860-0002</b>	<b>£46.00*</b>
5. ml	200	<b>E4860-0005</b>	<b>£25.50*</b>
10 ml	200	<b>E4860-0010</b>	<b>£33.00*</b>
25 ml	150	<b>E4860-0025</b>	<b>£41.00*</b>
50 ml	100	<b>E4860-0050</b>	<b>£50.00*</b>
100 ml	50	<b>E4860-1100</b>	<b>£60.00*</b>

Certified RNase, DNase, DNA, Latex and Pyrogen free.  
Sterile (SAL 10<sup>-9</sup>). Individually wrapped.



Replaceable rechargeable battery!

**PIPETTE CONTROLLER.**  
Pressure-responsive speed controls make the ErgoOne® FAST ideal for cell culture work such as resuspension of bacteria or cell pellets, nucleic acid precipitates and aspiration of a cell layer.

**ErgoOne® FAST**  
**Art. No. S7166-0010**  
**OFFER: £260\***

# Cryogenic Vials.

Supplied in convenient, self-standing resealable pouches.  
Certified. Sterile (SAL 10<sup>-9</sup>).



	Vol:	1.0 ml	1.8 ml	1.8 ml	3.6 ml	3.6 ml	4.5 ml	4.5 ml
	Base:	Skirted	Conical	Skirted	Conical	Skirted	Conical	Skirted
<b>EXTERNAL Threaded Tube Self-Sealing Cap, Sterile</b>	<b>E3090-6212</b>	<b>E3090-6321</b>	<b>E3090-6222</b>	n/a	n/a	n/a	<b>E3090-6242</b>	
	10 x 50	10 x 50	10 x 50				6 x 50	
	<b>£99.00*</b>	<b>£108.00*</b>	<b>£94.00*</b>				<b>£80.00*</b>	
<b>INTERNAL Threaded Tube, Silicone Seal in Cap , Sterile</b>	<b>E3110-6112</b>	<b>E3110-6121</b>	<b>E3110-6122</b>	<b>E3110-6131</b>	<b>E3110-6132</b>	<b>E3110-6141</b>	<b>E3110-6142</b>	
	10 x 50	10 x 50	10 x 50	8 x 50	8 x 50	6 x 50	6 x 50	
	<b>£119.00*</b>	<b>£131.00*</b>	<b>£109.00*</b>	<b>£99.00*</b>	<b>£120.00*</b>	<b>£103.00*</b>	<b>£103.00*</b>	

# ErgoOne® Starter Packs.



3 x ErgoOne® Single-Channel Pipettes:  
20 µl, 200 µl, 1000 µl  
3 x racks of sterile TipOne® Tips  
3 x Shelf-Adhesive Pipette Holders  
1 x Laboratory Timer

**STARTERPACK**  
**£405\***

4 x ErgoOne® Single-Channel Pipettes:  
2.5 µl, 20 µl, 200 µl, 1000 µl,  
3 x racks of sterile TipOne® Tips  
4 x Shelf-Adhesive Pipette Holders  
1 x Laboratory Timer

**STARTERPACK42.5**  
**£540\***

ErgoOne® pipettes paired with Starlab's trusted TipOne® pipette tips create the perfect pipetting partnership!

4 x ErgoOne® Single-Channel Pipettes:  
10 µl, 20 µl, 200 µl, 1000 µl,  
3 x racks of sterile TipOne® Tips  
4 x Shelf-Adhesive Pipette Holders  
1 x Laboratory Timer

**STARTERPACK410**  
**£540\***

# Cell Scrapers.

Gently raise/lift the cells with Starlab cell scrapers. The silicone blades of 220 mm and 250 mm scrapers prevent cell lysing. The flexible ABS handles on all scrapers provides easy handling inside flasks. Individually wrapped. Certified. Sterile.

220 mm straight handle, 11 mm fixed blade  
CC7600-0220  
Pack Size: 100  
SAL 10<sup>6</sup> **£136.00\***


250 mm angled handle, 11 mm pivoting blade  
CC7600-0250  
Pack Size: 100  
SAL 10<sup>5</sup> **£130.00\***

300 mm straight handle, 20 mm pivoting blade  
CC7600-0300  
Pack Size: 150  
SAL 10<sup>6</sup> **£115.00\***

# Win pizza for your team!

Enter our online competition for the chance to win a pizza lunch for your lab team. Scan the QR code or visit the link to take part.

Competition closes: 31.08.2026.



[starlab.click/win-pizza](http://starlab.click/win-pizza)

**\*OFFER REF: Q26000992193**




must be stated on your PO (excl. eShop). Offer valid to 31.08.26 / while stock lasts.

# Bottle Top Filtration.

**Flexible system!** Choose complete full assembly units for ready-to-use filtration, GL45 compatible bottle top funnels to pair with your own bottles or media bottles for versatile storage and preparation.

Superior filtration performance!



Bottle Top Filtration System					
	150 ml	250 ml	500 ml	1000 ml	
<b>Full Assembly Unit</b>					
	<b>0.22 µm PES membrane</b> Pack size: 12	<b>CC8221-3226</b> <b>£55.00*</b>	<b>CC8221-4226</b> <b>£63.00*</b>	<b>CC8221-5226</b> <b>£79.00*</b>	<b>CC8221-6226</b> <b>£123.00*</b>
<b>Funnel Tops only (GL45 compatible)</b>					
	<b>0.22 µm PES membrane</b> Pack size: 24	<b>CC8222-3226</b> <b>£65.00*</b>	<b>CC8222-4226</b> <b>£73.00*</b>	<b>CC8222-5226</b> <b>£92.00*</b>	<b>CC8222-6226</b> <b>£150.00*</b>
<b>Media Bottles only</b>					
	<b>Bottle and Cap</b> Pack size: 12	<b>CC8221-3000</b> <b>£50.00*</b>	<b>CC8221-4000</b> <b>£60.00*</b>	<b>CC8221-5000</b> <b>£60.00*</b>	<b>CC8221-6000</b> <b>£95.00*</b>

All products are certified RNase, DNase and Pyrogen free, are sterile (SAL 10<sup>-6</sup>), and supplied individually wrapped.

## CHEMGENE MEDLAB

Multi-surface disinfectant against bacteria, yeast, fungi, viruses and bacteriophages<sup>1</sup>. Denatures RNA/DNA.

**Go online for pricing!**



**Concentrate.** Economical dilution rates. Diluted solution effective for 3 months.



**1:10.** Convenient, ready to use spray and wipes.

\* not applicable to wipes.

# Syringe Filters.

Designed for efficiency and safety. Low hold-up volume means you can recover more of your valuable solution. 33 mm diameter filters supplied in easy-open blister packs showing sterile validity date (SAL 10<sup>-6</sup>).

**100 Filters: Double the usual pack size!**

0.22 µm PES Sterile.  
E4780-1226  
Pack Size: 100  
**£70.00\***

0.22 µm PVDF Sterile.  
E4780-1221  
Pack Size: 100  
**£70.00\***

0.22 µm CA Sterile.  
E4780-1223  
Pack Size: 100  
**£70.00\***

# Cell Strainers.

Cell Strainers provide a simple and reliable way to achieve uniform cell suspensions. Designed to fit perfectly into 50 ml centrifuge tubes. Gamma sterilised (SAL 10<sup>-6</sup>). Individually wrapped.

40 µm mesh. Sterile.  
CC8111-0042  
Pack Size: 50  
**£65.00\***

70 µm mesh. Sterile.  
CC8111-0072  
Pack Size: 50  
**£65.00\***

100 µm mesh. Sterile.  
CC8111-0102  
Pack Size: 50  
**£65.00\***