

Rigid Frame PCR Plates.

- Polycarbonate frame for precise automated handling
- Ultra-thin PP wells for optimal heat transfer, high reaction efficiency
- Ridged wells prevent evaporation when sealed
- White printed matrix
- RNase, DNase, DNA & PCR inhibitor free



Frame Colour	Well Colour	Cut Corner	Approx. Max. Well Capacity	Pack Size	Article No.	OFFER
96-WELL SEMI-SKIRTED						
Natural	Natural	A12	300 µl	10	E1404-0200	£30
Blue	Natural	A12	300 µl	10	E1404-0201	£30
96-WELL SKIRTED, LOW-PROFILE						
Natural	Natural	H1	200 µl	10	E1404-0100	£30
Blue	Natural	H1	200 µl	10	E1404-0101	£30
384-WELL SKIRTED						
Natural	Natural	A24	40 µl	10	E1404-0300	£36
Blue	Natural	A24	40 µl	10	E1404-0301	£36

Plate sealing.



Description	Art. No.	OFFER*
Self-Adhesive Plate Seal, Polypropylene, Clear	E2796-0793	£37
Self-Adhesive Sealing Strips, Polypropylene, Clear	E2796-2850	£29
Aluminium sealing film	E2796-0792	£220
Clear Polypropylene Seal (PCR)	E2796-0793	£150
Peelable Heat Sealing Film	E2796-3020	£125
Easi-Peel Heat Sealing Film PP	E2796-3060	£180
Easy Pierce Film (PCR)	E2796-3090	£80
Sealing Film for qPCR, self sealing	E2796-9895	£180
Sealing Foil, Aluminium	E2796-1100	£80

Black Friday week is coming. Be first to know!

A full week of savings across equipment and more is on the way.

Scan the code and we'll send you the discounts before they go live.



starlab.click/blackfriday

Storage offers.



Cardboard storage boxes



Polpropylene storage boxes



Cryogenic storage boxes

Explore the full offer range online.

Scan for full details!



starlab.click/offers

Freezer racking

SUSTAINABILITY JUST GOT COOLER

- Freezer racking helps optimise capacity by keeping freezers tidy and maximising storage space for efficient energy use
- High quality racking prevents the need for regular repair or exchange
- Manufactured in a factory powered by 100% sustainable energy
- Made from 85% recycled stainless steel
- Made in Germany. Short journey from factory to customer reduces carbon footprint



! *OFFER REF: Q26001022623

must be stated on your PO (excl. eShop).

Offer valid to 16.12.26 / while stock lasts.

Save on equipment,
spend more on research.

You are passionate
for science.

SPECIAL OFFER HIGHLIGHTS

Valid until 16.12.2026




ORDER UP TO
4PM FOR NEXT
DAY DELIVERY


IN THIS EDITION:
LAB EQUIPMENT
PIPETTES
SAMPLE STORAGE
PCR PLATES AND SEALS
FREEZER RACKING



Our lowest ever priced equipment for your lab.



Description	Art. No.	OFFER*
MIXING, ROLLING, STIRRING AND SHAKING.		
Magnetic Stirrer	N2400-2010	£145
Magnetic Stirrer with Heating	N2400-3010	£275
Rotator Mixer Multi-1, Programmable (incl. Platform N2400-5104)	N2400-5000	£585
Tube Roller	N2400-7010	£435
Rocking Shaker	N2400-8010	£860
3D Rocking Shaker	N2400-8020	£860
Orbital Shaker (a platform must be purchased separately)	N2400-8030	£755
Microplate Mixer	N2400-8040	£300
Plate Centrifuge	N2631-0100	£695
Vortex Mixer, 100-240 V	S8010-0000	£290
Vortex Mixer, Standard attachment	S8010-0010	£6
Vortex Mixer, Adapter for 1.5/2.0 ml tubes	S8010-0011	£18
Vortex Mixer, Adapter for microplates	S8010-0012	£12
Minicentrifuge	S8030-5000	£195
Minicentrifuge, tube rotor, 8 place	S8030-5010	£18.50
Minicentrifuge, rotor for PCR tubes	S8030-5011	£18.50
GUARDONE® WORKSTATION.		
GuardOne laminar flow 32in, Starlab	S8040-0000	£3,725
GuardOne laminar flow 48in, Starlab	S8040-0001	£4,900
GuardOne dead air 32in, Starlab	S8040-1000	£3,050
GuardOne dead air 48in, Starlab	S8040-1001	£3,900
LABEL PRINTER.		
Brady® M210-LAB Printer	N3000-2000	£120



- Ergonomic and lightweight
- Minimal operating forces
- Volume lock
- Universal tip fit (excl. 5.0 ml)
- 3 year warranty

Manual pipettes.

Single-Channel Manual Pipettes


S7100-0125 0.1 – 2.5 µl £135	S7100-0510 0.5 – 10 µl £135	S7100-0221 2 – 20 µl Micro £135	S7100-0220 2 – 20 µl Standard £135
S7100-1100 10 – 100 µl £135	S7100-2200 20 – 200 µl £135	S7110-1000 100 – 1,000 µl £135	S7150-5000 500 – 5,000 µl £135

8-Channel Pipettes

S7108-0510 0.5 – 10 µl £435	S7108-1100 10 – 100 µl £435	S7108-3300 30 – 300 µl £435
-----------------------------------	-----------------------------------	-----------------------------------

12-Channel Pipettes

S7112-0510 0.5 – 10 µl £525	S7112-1100 10 – 100 µl £525	S7112-3300 30 – 300 µl £525
-----------------------------------	-----------------------------------	-----------------------------------



RepeatOne S®
Dispenses up to 48 precise volumes from a single sample.
1-year warranty.

Art. No: G4026-1500
£242



ErgoOne® carousel
Holds up to 8 pipettes.
Art. No: S7200-0000
£95 / each

ErgoOne Wall Mount
Art. No: S7200-0001
£5 / each

Acrylic Stand
Holds up to 5 ErgoOne Pipettes
Art. No: E3099-0000
£30 / each



ErgoOne® Fast pipette controller.
The ErgoOne® FAST is one of the lightest pipette controllers on the market.

Art. No: S7166-0010
£254

Electronic pipettes.



- Easy to use, intuitive menu
- Lightweight, ergonomic design
- Universal tip fit
- 12-month warranty

Single-Channel ErgoOne® E Pipettes

G9001-0010 0.5 – 10 µl £370	G9001-0020 2.0 – 20 µl £370	G9001-0200 10 – 200 µl £370	G9001-1000 50 – 1,000 µl £370	G9001-5000 250 – 5,000 µl £370
-----------------------------------	-----------------------------------	-----------------------------------	-------------------------------------	--------------------------------------

8-Channel ErgoOne® E Pipettes

G9008-0010 0.5 – 10 µl £680	G9008-0020 1.0 – 20 µl £680	G9008-0100 5.0 – 100 µl £680	G9008-0300 15 – 300 µl £680	G9008-1250 50 – 1250 µl £680
-----------------------------------	-----------------------------------	------------------------------------	-----------------------------------	------------------------------------

12-Channel ErgoOne® E Pipettes

G9012-0010 0.5 – 10 µl £805	G9012-0020 1.0 – 20 µl £805	G9012-0100 5.0 – 100 µl £805	G9012-0300 15 – 300 µl £805	G9012-1250 50 – 1250 µl £805
-----------------------------------	-----------------------------------	------------------------------------	-----------------------------------	------------------------------------