

Cell culture challenges, but
you continue to grow.
**You are passionate
for science.**

SPECIAL OFFERS

Cell Culture Consumables & Equipment

Valid until 27.08.2025



WATCH VIDEO

Passionate for Science.
Starlab inspires and
supports the passion
driving your research.

Passionate for science.



CytoOne® cell culture ware.

YOUR CELLS DESERVE THE BEST

Designed for optimal performance and manufactured to the highest quality standards. Outfit your cells for true performance with CytoOne® cell culture products.

CytoOne® products are Gamma sterilised (SAL 10⁻⁵).

CytoOne® Dishes.

CytoOne® dishes feature stacking rings with vented design for thermal equalisation. Vented lids improve gas exchange.

FROM £39.53*



Size	Pack Size	Art. No.	List	OFFER*
TC Treated Dishes				
35 x 10 mm	30 x 10	CC7682-3340	£125.31	£41.35
60 x 15 mm	30 x 10	CC7682-3359	£127.50	£39.53
100 x 20 mm	30 x 10	CC7682-3394	£220.19	£85.87
150 x 20 mm	12 x 5	CC7682-3614	£134.09	£50.95

starlab.click/dishes

Convenient, small pack sizes and easy to open packaging!

CERTIFIED



CytoOne® Flasks.

All flasks are 100% pressure tested for leak-free assurance. Features stacking rims with vented design for thermal equalisation.



FROM £50.78*

Flasks with Vented Filter Caps

Size	Pack Size	Art. No.	List	OFFER*
TC Treated Flasks, Vented Filter Caps				
T-25	30 x 10	CC7682-4825	£508.40	£213.53
T-75	20 x 5	CC7682-4875	£301.01	£126.42
T-150	8 x 5	CC7682-4815	£226.43	£83.78
T-225	5 x 5	CC7682-4822	£192.35	£80.79



CytoOne® vented filter caps have a 0.2 µm hydrophobic membrane for gas exchange and to exclude bacteria and fungal contaminants

Flasks with 2-Position Plug Seal Caps

Size	Pack Size	Art. No.	List	OFFER*
TC Treated Flasks, 2-Position Plug Seal Caps				
T-150	8 x 5	CC7682-4415	£191.26	£66.94
T-225	5 x 5	CC7682-4225	£158.36	£61.76
Non-Treated Flasks, 2-Position Plug Seal Caps				
T-75	20 x 5	CC7672-4175	£193.29	£81.18
T-225	5 x 5	CC7672-4225	£153.88	£50.78



CytoOne® two-position plug seal caps allow you to choose secure manual venting or an airtight seal

starlab.click/flasks



CERTIFIED



CytoOne® Plates.

FROM £38.45*

Wells	Pack Size	Art. No.	List	OFFER*
TC Treated Plates				
6-Well	50	CC7682-7506	£98.59	£38.45
12-Well	50	CC7682-7512	£130.32	£50.82
24-Well	50	CC7682-7524	£139.37	£44.60
48-Well	50	CC7682-7548	£141.28	£55.10
96-Well	50	CC7682-7596	£142.78	£52.83

starlab.click/multiwellplates




The complete 360° chimney well design on all CytoOne® plates prevents cross-contamination and ensures uniform temperature transfer across all wells, resulting in a more uniform cell layer. Watch our video to see how CytoOne® plates are designed to prevent the edge effect so you can use all 96-wells of a plate.



WATCH VIDEO

PLATE SEALING



Porous, 50 µm self-adhesive film that allows gas exchange for cellular and bacteria cultivation. Sterile. RNase, DNase & Pyrogen free. Temp use: 4° to 40° C.

Breathable Film

Pack size: 100

Art. No. E2796-3005

List: £109.49

£75.55*

Breathable Film (Sterile)

Pack size: 50

Art. No. E2796-3015

List: £83.48

£57.60*



Clear, 50 µm self-adhesive polyester seal ideal for long-term storage & incubation. Temp use: -4.0 °C to 120 °C.

Polyester Seal

Pack size: 100

Art. No. E2796-0101

List: £58.38

£43.79*

Polyester Seal (Sterile)

Pack size: 100

Art. No. E279-0102

List: £79.62

£59.72*

starlab.click/plateseals



Cell Scrapers.

All cell scrapers feature flexible ABS handles, and are sterile and individually-wrapped.



220 mm straight handle, 11 mm fixed blade

Ideal for stem cells due to reduced shearing from gentle silicone blade. Use to gently lift cells from surface rather than scrape. Use for plates and dishes.

Art. No. CC7600-0220

Pack Size: 100

List: ~~£208.84~~ **£137.83***



250 mm angled handle, 10 mm pivoting blade

Ideal for stem cells due to reduced shearing from gentle silicone blade. Flexibility allows you to reach all corners of your small flask.

Art. No. CC7600-0250

Pack Size: 100

List: ~~£197.85~~ **£130.58***



300 mm straight handle, 20 mm pivoting blade

Large pivoting polyethylene blade for easy access and use in large flasks.

Art. No. CC7600-0300

Pack Size: 150

List: ~~£188.02~~ **£116.57***

starlab.click/cellscrapers



***OFFER REF: Q25000782740**

must be stated on your PO (excl. eShop).
Offer pricing valid to 27.08.2025.

Aspirator.

ONLY: £1,360.54*



See how
to add
a second
hand
controller.

- The powerful pump can operate two hand controllers without any loss of performance
- Chemically-resistant, high-performance pump with a long service life
- Auto ON-OFF pump switch with low noise emission for a quieter working environment
- Ergonomic, fully autoclavable hand controller
- Intuitive, reproducible adjustment of suction power via touch panel



Aspirator

Incl. everything you need to get started! 1-year warranty.

N2400-9000

£2,230.40

£1,360.54



Hand Controller for Aspirator

Incl. 1x Adapter for Pipettes, 1x Adapter for Pipette Tips, 1x Wall Mount, and 2.5 m Tubing.

N2400-9001

£202.95

£127.86



Connection Kit for Second Hand Controller

N2400-9006

£50.21

£31.63



starlab.click/aspirator

Aspiration pipettes.

- Accuracy: $\pm 2\%$
- Colour coded, paper-plastic peel-back wraps
- Gamma sterilised to SAL 10^{-6} (ISO 11137)
- Certified RNase, DNase, Pyrogen and latex free
- Non-cytotoxic and non-haemolytic USP Class VI Polystyrene



Volume	Pack Size	Art. No.	List	OFFER*
2.0 ml	500	E4861-0002	£110.72	£54.25
5.0 ml	200	E4861-0005	£51.67	£25.32
10 ml	200	E4861-0010	£66.44	£32.56



starlab.click/aspirationpipettes



Aspiration pipettes are supplied in convenient dispenser-style boxes for easy access when stacked



***OFFER REF: Q25000782740**

must be stated on your PO (excl. eShop). Offer valid to 27.08.2025.

ErgoOne® FAST pipette controller.

Use with all
glass / plastic
pipettes from
1.0 to 100 ml.



FAST!
Fills a 25 ml
pipette in
less than
4 seconds.

Pressure-responsive speed
controls make the ErgoOne®
FAST ideal for cell culture
work such as resuspension
of bacteria or cell pellets,
nucleic acid precipitates
and aspiration of a cell layer.



WATCH VIDEO

ErgoOne® FAST

Art. No. S7166-0010

List: £396.03

OFFER: £241.58*

starlab.click/ergoonefast

See our range of serological pipettes at: starlab.click/serologicals

Completely charged,
the battery can
perform 2000 x 25 ml
dosing cycles. The
battery is also very
easy to replace!



StarGuard® Protect.



**ONLY
£8.05
per 100 gloves**

CE 2777
EU Regulation
2016/425

EN ISO 374-1:2016+A1:2018
Type B
KPT

EN ISO 374-5:2016
VIRUS

AQL
0.65

Food Safe



WATCH VIDEO

TRUSTED, RELIABLE PROTECTION

- Combines strength and durability for everyday protection
- Comfort, tactile sensitivity and consistent grip
- Excellent chemical splash protection
- Longer, 250 mm glove to cover the wrist

Glove size	Pack Size	Art. No.	List	OFFER*
X-Small	10 x 100	SG-P-XS	£161.09	£80.55
Small	10 x 100	SG-P-S	£161.09	£80.55
Medium	10 x 100	SG-P-M	£161.09	£80.55
Large	10 x 100	SG-P-L	£161.09	£80.55
X-Large	10 x 100	SG-P-XL	£161.09	£80.55

starlab.click/sg-protect



***OFFER REF: Q25000782740**

must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

ErgoOne® Starter Packs

Reliability is key! Your pipette must be trusted to provide accurate and precise pipetting, day in, day out. ErgoOne® manual pipettes provide you with both the reliability you need and an ergonomic design that prevents fatigue over long periods to help keep your pipetting consistent. ErgoOne® pipettes have tighter specifications than ISO 8655, and feature a volume lock to ensure the correct volumes in your reaction.



CONTENTS	STARTERPACK	STARTERPACK42.5	STARTERPACK410
0.1 – 2.5 µl ErgoOne® Single-Ch Pipette S7100-0125		●	
0.5 – 10 µl ErgoOne® Single-Ch Pipette S7100-0510			●
2.0 – 20 µl ErgoOne® Single-Ch Pipette, Standard Cone S7100-0220	●	●	●
20 – 200 µl ErgoOne® Single-Ch Pipette S7100-2200	●	●	●
100 – 1000 µl ErgoOne® Single-Ch Pipette S7110-1000	●	●	●
TipOne® Pipette Tips (three racks of sterile tips per pack)	●	●	●
Self-Adhesive Pipette Holder (one holder per pipette supplied)	●	●	●
Laboratory Timer (one timer per pack)	●	●	●
List:	£698.63	£905.63	£905.63
OFFER:	£405.21*	£534.32*	£534.32*

ErgoOne® pipettes paired with Starlab's trusted TipOne® pipette tips create the perfect pipetting partnership!

 starlab.click/starterpacks

Cryogenic Vials.

FROM: £76.77*

- 100% polypropylene tubes
- For storage down to -196 °C; **should only be used in the gas phase of liquid nitrogen**
- Internal threaded vials have a silicone seal in the cap for the highest possible impermeability
- White writing area
- Highly visible black graduations at 100 µl increments
- Sterile (SAL 10⁻⁶)
- Certified RNase, DNase, DNA and Endotoxin free



Supplied in convenient, resealable, self-standing pouches.




EXTERNAL Threaded Tube Self-Sealing Cap, Sterile



INTERNAL Threaded Tube, Silicone Seal in Cap Sterile

CRYOGENIC VIALS (STERILE)

Images in ordering table show internal threaded tubes.



Volume:	1.0 ml	1.8 ml	1.8 ml	3.6 ml	3.6 ml	4.5 ml	4.5 ml
Base:	Skirted	Conical	Skirted	Conical	Skirted	Conical	Skirted
EXTERNAL Threaded Tube Self-Sealing Cap, Sterile	E3090-6212 10 x 50 £156.84 £101.95*	E3090-6321 10 x 50 £163.67 £109.66*	E3090-6222 10 x 50 £156.39 £95.40*	n/a	n/a	n/a	E3090-6242 6 x 50 £119.96 £76.77*
INTERNAL Threaded Tube, Silicone Seal in Cap Sterile	E3110-6112 10 x 50 £183.86 £121.35*	E3110-6121 10 x 50 £202.36 £137.60*	E3110-6122 10 x 50 £186.10 £115.38*	E3110-6131 8 x 50 £158.08 £99.59*	E3110-6132 8 x 50 £165.32 £114.07*	E3110-6141 6 x 50 £159.08 £104.99*	E3110-6142 6 x 50 £159.08 £104.99*

starlab.click/cryovials

Cryogenic Storage Boxes.



81-place & 100-place
Polycarbonate boxes with fixed dividers. Transparent lid printed with reference grid and has a one-way fit so your inventory is safe.
130 x 130 x 53 mm.

I2381-2305	£46.49
81-place box (x5)	£21.39*
I2510-5005	£52.65
100-place box (x5)	£21.59*



81-place tall
Polycarbonate boxes with fixed dividers. Matrix printed on dividers. Lift-off transparent lid.
133 x 133 x 95 mm.

E2382-6002	£87.20
Green (x5)	£43.50
E2382-6003	£87.20
Orange (x5)	£43.50
E2382-6006	£87.20
Yellow (x5)	£43.50

Available only whilst stock lasts!



25-place
Polycarbonate boxes with fixed dividers. Matrix printed on dividers. Lift-off transparent lid.
133 x 133 x 54 mm.

E2381-5001	£43.54
Blue (x5)	£21.75
E2381-5002	£43.54
Green (x5)	£21.75
E2381-5003	£43.54
Orange (x5)	£21.75
E2381-5004	£43.54
Red (x5)	£21.75

Available only whilst stock lasts!



***OFFER REF: Q25000782740**

must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

Erlenmeyer Flasks.

Autoclavable
and reusable
a minimum of
10 times!

 starlab.click/erlenmeyers

DESIGNED TO BE REUSED

Ideal for the cultivation of bacteria, yeast and fungi, as well as plant and animal cells in suspension.

PRODUCT FEATURES:

- Optically clear polycarbonate flasks
- Moulded graduations
- Leakproof and shatterproof
- Dual-function, vented polypropylene cap provides sterile air exchange or a tight seal
- Suitable for both aerobic and anaerobic cultures
- 0.22 µm PTFE cap membrane
- Individually-wrapped
- Sterility level (SAL 10⁻⁶)
- Certified RNase, DNase, DNA & Pyrogen free
- Autoclavable
- Reusable a minimum of 10 times

	Volume	Pack Size	Art. No.	List	OFFER*
FLAT BASE 	125 ml	24	E5000-1012	£207.20	£145.04
	250 ml	12	E5000-1025	£136.69	£95.68
	500 ml	12	E5000-1050	£179.40	£125.58
	1000 ml	6	E5000-1100	£125.37	£87.76
	2000 ml	6	E5000-1200	£158.10	£110.67
BAFFLED BASE 	3000 ml	3	E5000-1300	£189.27	£132.49
	125 ml	24	E5001-1012	£207.20	£145.04
	250 ml	12	E5001-1025	£136.69	£95.68
	500 ml	12	E5001-1050	£179.40	£125.58
	1000 ml	6	E5001-1100	£125.37	£87.76
	2000 ml	6	E5001-1200	£158.10	£110.67
	3000 ml	3	E5001-1300	£189.27	£132.49

All products are sterile (SAL 10⁻⁶) and individually wrapped.

DUAL-FUNCTION CAPS

Polypropylene caps with 0.22 µm PTFE membrane, ideal for both aerobic and anaerobic cultures. Remove the top section to allow sterile exchange, and then replace when ready. Avoid contamination risk and the hassle of switching caps and stock only one flask type for all cultures.



TOTALLY TUBULAR

Four, 29-place cube racks that hold a variety of tube sizes, in a support tray. Overall WxDxH: 93 x 93 x 48 mm



12344-1000
£50.96
£33.63*

COMBI MINI RACKS

80-place polypropylene racks for use with a variety of tube sizes. Ideal for use in the hood. WxDxH: 171 x 91 x 54 mm



E2345-1000	£8.87
Natural x1	£5.32*
E2345-1041	£8.97
Neon Blue x1	£4.40*
E2345-1042	£8.98
Neon Green x1	£4.31*
E2345-1043	£8.97
Neon Orange x1	£4.31*
E2345-1044	£8.98
Neon Pink x1	£4.40*
E2345-1046	£8.97
Neon Yellow x1	£4.31*
E2345-1099	£42.40
Mixed Neon x5	£27.98*

Note: Individual neon colour racks only available whilst stocks last.



***OFFER REF: Q25000782740**

must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

Bottle top filtration.

Certified
RNase, DNase,
Pyrogen,
non-cytotoxic
& heavy-metal
free.



PES Filters are ideal for pre-filtering before PFAS analysis

BOTTLE TOP FILTRATION SYSTEM

	150 ml	250 ml	500 ml	1000 ml
--	--------	--------	--------	---------

Full Assembly Unit

0.22 µm

PES membrane

Pack size: 12

CC8221-3226

List: £87.39

£56.80*

CC8221-4226

List: £98.50

£63.04*

CC8221-5226

List: £123.40

£80.21*

CC8221-6226

List: £197.01

£128.06*

0.45 µm

PES membrane

Pack size: 12

CC8221-3456

List: £87.02

£56.56*

CC8221-4456

List: £97.37

£63.29*

CC8221-5456

List: £124.53

£79.70*

CC8221-6456

List: £197.01

£128.06*

Funnel Tops only

0.22 µm

PES membrane

Pack size: 24

CC8222-3226

List: £81.08

£64.86*

CC8222-4226

List: £90.43

£73.25*

CC8222-5226

List: £114.14

£93.59*

CC8222-6226

List: £182.41

£153.22*

0.45 µm

PES membrane

Pack size: 24

CC8222-3456

List: £81.08

£64.86*

CC8222-4456

List: £90.43

£73.25*

CC8222-5456

List: £114.14

£93.59*

CC8222-6456

List: £182.41

£153.22*

Media Bottles only

Bottle & Cap

Pack size: 24

CC8221-3000

List: £65.21

£46.95*

CC8221-4000

List: £73.53

£52.94*

CC8221-5000

List: £81.56

£58.72*

CC8221-6000

List: £127.87

£92.07*

All products are sterile (SAL 10⁻⁶) and individually wrapped.

starlab.click/filtration



***OFFER REF: Q25000782740**

must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

15 ml & 50 ml centrifuge tubes.

CERTIFIED



DESIGNED FOR THE CENTRIFUGATION OF CELLS

USP Class VI polymers and certified non-toxicity ensure the viability and integrity of your cells. The sterility assurance level (SAL) of 10^{-6} fulfills the highest demands on sterility, meeting EN 556-1 standard (Sterilisation of Medical Devices).

- White writing surface and clear, black print graduations
- Longer length screw caps with sealing ring prevents leakage



15 ML CENTRIFUGE TUBES

FROM £53.78*



50 ML CENTRIFUGE TUBES

FROM £69.28*

Description	Pack Size	Art. No.	List	OFFER*
15 ml Centrifuge Tubes				
15 ml Centrifuge Tube, Conical, Loose (Sterile)	20 x 25	E1415-0100	£141.52	£53.78
15 ml Centrifuge Tube, Conical, Racked (Sterile)	20 x 25	E1415-0700	£147.16	£61.81
50 ml Centrifuge Tubes				
50 ml Centrifuge Tube, Conical, Loose (Sterile)	20 x 25	E1450-0100	£150.60	£69.28
50 ml Centrifuge Tube, Conical, Racked (Sterile)	20 x 25	E1450-0700	£158.51	£77.67



starlab.click/centrifuge tubes

Centrifuge Tube Storage.



15-place rack for 15 ml tubes

Chemically resistant epoxy-coated wire rack, blue.
127 x 76 x 176 mm

E2396-1515
(1 Rack)

List:
£38.36

OFFER:
£25.70*



8-place rack for 50 ml tubes

Chemically resistant epoxy-coated wire rack, blue.
147 x 80 x 89 mm

E2396-1550
(1 Rack)

List:
£37.19

OFFER:
£24.55*



60-place rack for 15 ml tubes

Acetal-polymer rack. White/Mixed colours*.
293 x 115 x 65 mm

E2396-9769
(5 Racks)

List:
£96.23

OFFER:
£66.40*



18-place rack for 50 ml tubes

Acetal-polymer rack. White/Mixed colours*.
293 x 115 x 65 mm

E2396-9779
(5 Racks)

List:
£96.33

OFFER:
£66.47*

*Contains one each of blue, green, violet, yellow and orange.



***OFFER REF: Q25000782740**

must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

Robotic Tips.

Biomek® i-Series • Biomek® FX • Biomek® NX •
Biomek® 3000 • Biomek® 4000 • Multimek®



RACKED TIPS (10 x 96)

20 µl Tip, Racked	E1076-4400	£67.29	£37.68*
20 µl Tip, Racked (Sterile)	E1076-4410	£79.34	£42.05*
20 µl Filter Tip, Racked (Sterile)	E1076-4418	£186.04	£98.60*
50 µl Tip, Racked	E1076-2400	£77.35	£37.90*
50 µl Tip, Racked (Sterile)	E1076-2410	£97.72	£42.02*
50 µl Filter Tip, Racked (Sterile)	E1076-2418	£186.80	£99.00*
250 µl Tip, Racked	E1076-0400	£83.48	£37.57*
250 µl Tip, Racked (Sterile)	E1076-0410	£80.02	£42.41*
175 µl Filter Tip, Racked (Sterile)	E1076-0418	£186.79	£99.00*

REFILLABLE RACKS (10 x 96)

20 µl Tip, Refillable Rack	E1076-4800	£58.91	£38.88*
250 µl Tip, Refillable Rack	E1076-0800	£73.03	£39.44*


REFILLS (10 x 96)

20 µl Tip, Refill	E1076-4700	£63.47	£37.45*
250 µl Tip, Refill	E1076-0700	£66.71	£37.36*

Certified RNase, DNase and DNA & Pyrogen free. Tested Pyrogen free. Sterile product is certified pyrogen free. Non refillable racks also compatible with Gilson Cyberlab.
Available whilst stocks last.

starlab.click/biomek

Packard MultiProbe II, MultiProbe IIHT/EX.
Tecan (with LiHa) Freedom EVO. Genesis.
MiniPrep.



Supplied in individual
blister packs (24 x 96)

TipOne® ROBOTIC TIPS (24 x 96)

50 µl Tip	R1181-5700	£518.57	£165.94*
50 µl Filter Tip	R1128-5810	£535.16	£176.60*
200 µl Tip	R1181-0700	£237.59	£133.05*
200 µl Filter Tip	R1128-0810	£259.11	£147.69*
1000 µl Tip	R1181-2700	£248.26	£148.96*
1000 µl Filter Tip	R1128-2810	£282.51	£161.03*

TipOne® CONDUCTIVE ROBOTIC TIPS (24 x 96)

50 µl Conductive Tip	R1181-5709	£518.57	£212.61*
50 µl Conductive Filter Tip	R1128-5819	£530.69	£222.89*
200 µl Conductive Tip	R1181-0709	£234.42	£131.28*
200 µl Conductive Filter Tip	R1128-0819	£290.75	£165.73*
1000 µl Conductive Tip	R1181-2709	£257.50	£144.20*
1000 µl Conductive Filter Tip	R1128-2819	£341.67	£194.75*

All tips are non-sterile and certified RNase, DNase, DNA & Pyrogen free. Also compatible with Qiagen® QIAgility® (CAS-1200) with single-channel head and Aurora Biomed® Systems with 8-channel head. Available whilst stocks last.



Seated in white hanging racks for use with robots
with Varispan 4 or 8 head probes

starlab.click/tipone-tecan



***OFFER REF: Q25000782740**
must be stated on your PO (excl. eShop).
Offer valid to 27.08.2025.

TipOne® 'closed-loop products' getting closer!

We are thrilled to announce that TipOne® racks made with recycled plastic from the TipOne® closed-loop recycling service will soon be available. The base of all TipOne® racks made in Germany¹⁾ will contain **75% TipOne® recyclate** plus 25 % recyclate from another source.

Providing a closed-loop recycling service is part of our commitment to sustainability and reducing environmental impact. With TipOne®, you not only receive the high-quality pipette tips you have come to expect, you know you are using a truly sustainable and plastic-saving pipette tip system.

How will I know when my TipOne® rack contains TipOne® recycled plastic?

Look out for the new design shelf boxes for racked tips made in Germany. Just like the blue shelf box for products made in USA²⁾, the latest shelf boxes are also made with 100 % recycled cardboard, but now also boast a design that uses a lot less ink, and have labels made from recycled paper. More small but important steps on the TipOne® sustainability journey.

starlab.click/tipone

¹⁾ TipOne® manufacturing is moving from USA to Germany. German production started January 2025; there will be a lengthy transition phase.

²⁾ During the transition phase you will receive TipOne® products from either manufacturing site.



TipOne® racks made in Germany that contain TipOne® recyclate will be supplied in a new-design shelf box!



TipOne® shelf boxes of racks made in USA.



T: +44 (0)1908 283800
info@starlab.co.uk
www.starlab.co.uk

Passionate for science.

